
POCKET PRO SERIES

PROFESSIONAL PORTABLE POCKET SCALE

Instruction Manual

Manual de Instrucciones

Bedienungshandbuch

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ACCULAB

sartorius group

- PP-62
- PP-201
- PP-401

ENGLISH

- Instruction Manual 1

ESPAÑOL

- Manual de Instrucciones 4

DEUTSCH

- Bedienungshandbuch 7

FRANÇAIS

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ITALIANO

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Thank you for purchasing Acculab Pocket Scale. Please read all of the operating instructions to ensure best performance. The unit is pre-calibrated at the factory and comes ready to use with one AA alkaline battery (not engaged).

Operating Instructions:

Note: Before using, engage AA battery (supplied with scale/battery compartment located on bottom of scale).

Slide pan (as marked on scale) completely open towards the rear of the unit (should “click” in place). Turn scale on by pressing “ON/LIGHT” key. A segment display will appear (“----” / then “0.0g” or “0.00g”) and scale is ready for operation. Turn scale “OFF” and slide pan forward when not in use to close and lock during transportation.

Key Functions / Descriptions:

ON/LIGHT Key:

Press to turn balance “ON”. Within the operation mode, press this key to “turn-on” backlight.

CAL/UNIT Key:

Press and “release” key will enable toggling between grams (g) and ounces (oz). (Cycle g, oz, gn and ct for PP-62)

Press and “hold” (for approximately 3 seconds) puts balance in calibration mode.

ZERO/OFF Key:

Press and “release” key will return display to zero.

Press and “hold” key (for approximately 3 seconds) turns balance off.

Features:

Power-Up Test: When the unit is initially turned ON, “----” will appear in the display (for approximately 2 seconds), then “0.0g” or “0.00g” depending upon model will be displayed. This is a self diagnostic test.

Low Battery: when battery is in low battery state, ie, 1.2 -0.9V, the display will show “LO” to indicate low battery. It is recommended to change battery at this time.

Auto Off: This unit comes with an auto-off feature to conserve battery power. The scale will automatically turn off after approximately 4 minutes of inactivity (note: the backlight will turn off automatically after 1 minute to conserve battery life, simply press the “ON/LIGHT” key to turn backlight “ON”).

Overload Indicator: When a load that is applied to the weighing platform exceeds the scale capacity, an “H” will appear in the display. Please remove load immediately (in order to prevent damage) and the unit should return to normal operation.

Caution: Overloading the capacity of the scale should be prevented as this may cause damage to the weight cell and void warranty. Additionally please avoid vibrations, strong electrical or magnetic fields, and rapid changes in humidity or temperature as they can adversely affect accuracy.

Calibration:

Calibration (adjustment) ensures that the weighing precision of your scale is maintained according to specifications. Calibration weight is included with package.

Calibration Procedure:

- 1 Press and "hold" CAL/UNIT key (for approximately 3 seconds) or until the display shows flashing weight value (example: "100.0" for model: PP-201). The flashing weight value indicates the correct weight that must be placed on the center of the platform for calibration (PP-201= 100gram / PP-401= 200gram / PP-62= 50gram).
- 2 Place the noted weight onto the center of the platform (note: if no weight is placed on the pan within 12 seconds an "E" will be displayed and the unit will return to weigh mode).
- 3 If weight is accepted (typically within 10 seconds or less) the display will stop flashing and briefly (4 seconds) display weight value with "g" symbol.
- 4 Four dashes "----" will then be shown indicating calibration is complete.
- 5 Remove calibration weight and display will return to zero reading.
- 6 If a calibration error has occurred an "E" will be displayed for 5 seconds and the unit will return to the weigh mode. (please repeat steps 1-5)

Important:

Please allow the scale to stabilize and wait at least 30 seconds after turning scale on, before calibration. During prolong periods of use repeat this procedure at regular intervals for optimum accuracy. Before placing the weight onto the pan, make sure the pan is thoroughly opened.

Note: Please do not open or attempt to internally calibrate the scale. This will void the warranty and may cause damage.

Care and Maintenance:

- Do not clean with liquids or aggressive cleaning agents.
- Do not use in liquids or water.
- Avoid lengthy exposure to extreme heat or cold.
- Allow scale a minimum of thirty seconds to warm up before calibration to allow internal components to stabilize.
- Keep product clean. Avoid: dust; dirt, & liquids from getting inside of scale, which can adversely affect and damage the reliability and accuracy of the scale.
- Vibrations, air currents and proximity to other electronic equipment may affect accuracy of this product. For optimum performance please ensure weighing is done on a flat, solid surface in a clean, stable environment.
- Replace low or dead batteries.

Product Specifications

MODEL NO.	PP-62	PP-201	PP-401
Capacity	65 grams	200 grams	400 grams
Readability	0.01 grams 0.001 oz 0.2 gn 0.1 ct	0.1 grams 0.005 oz	0.1 grams 0.01 oz
Weighing Modes	g, oz, gn, ct	g, oz	g, oz
User Calibration	at 50g	at 100g	at 200g
Cal. Weight (Included)	50g	100g	200g
Auto Off	No Activity for 4 Minutes		
Low Battery Indication	Between 0.9 ~1.2 V, Display "LO"		
Over Load Indication	Display Shows "H" with Beep		
Under Load Indication	Four Dashes "----"		
Power	One AA Battery		
Scale Dimension (WxLxH)	75 mm x 90 mm x 23 mm / 2.94" x 3.5" x 0.88" (Closed)		
Pan Dimension (WxL)	73 mm x 66 mm / 2.88" x 2.5"		
Estimate Battery Life	12 hrs with Backlite, 20 hrs Normal Operation		
Error or Display Codes			
"LO"	Low Battery		
"H"	Overload Indicator		
"E"	Calibration Error		
"----"	Under Load Indicator		

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