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## Index

<b>1 Safety- and user instructions</b>	<b>3</b>
1.1 General safety instructions	3
1.2 Use	4
1.3 General user instructions	4
<b>2 Commissioning</b>	<b>5</b>
<b>3 Functional characteristics</b>	<b>6</b>
3.1 Switch on	6
3.2 Air-purging	6
3.3 Titration	7
3.4 Dosing	7
3.5 Draining	8
3.6 Charging	8
<b>4 Titrexsoft</b>	<b>9</b>
<b>5 Calibration mode</b>	<b>10</b>
5.1 Calibration without software and PC	10
<b>6. Cleaning and maintenance</b>	<b>11</b>
6.1 Disassembly/Assembly	11
6.2 Cleaning	11
<b>7 Accessories</b>	<b>12</b>

# 1 Safety- and user instructions

## 1.1 General safety instructions

Highest safety precautions should be used when dispensing corrosive, poisonous, radioactive or hazardous chemicals.

- General operation and safety instructions must be followed. Eye protection, protective clothing and glove must be worn.
- Every user must study these instructions carefully before using the instrument.
- The instrument must only be used for intended purpose and utilized the component tolerances of material resistance.
- Always work with the instrument so as not to endanger the operator or other persons.
- Always check the dispenser for leaks before using. Make sure that the suction and ejection cannula are firmly in position.
- Do not apply force to the instrument.
- Do not drop the instrument.
- Do not immerse the instrument in liquids.
- Only use original spares and accessories.
- Repairs must only be carried out by the manufacturer's trained personnel. Do not open the instrument. Infringements will void your warranty!
- Do not use the instrument in explosive environments.
- Do not expose the instrument to temperatures exceeding 50°C. Do not autoclave!
- Keep the instrument away from fire or flames.
- We only consider ourselves responsible in regard to safety, reliability and performance of the device if:
  - The dispenser has been used in accordance to the instruction manual
  - Only the parts supplied by the manufacturer have been used

## 1.2 Use

The TITREX Digital Burette is intended to be used for titration and dosing within following limits:

15 – 40 °C of reagent and instrument

When using according to the instructions, only these parts of the Digital Burette will come into direct contact with following chemically resistant material:

- FEP (Suction hose, discharge tube, return dosing tube)
- Tefzel® (Valve block)
- Borosilicate glass (Valves, KPG Cylinder)
- Ceramic, Platinum-Iridium (balls and valves)

Never use the instrument with

- Hydrofluoric acid
- liquids which corrode PTFE, FEP, Tefzel®, Platinum-Iridium or where catalytic changes to these material can occur
- crystallizing, particulating or precipitating solutions (e.g. suspensions)

The including adapters are made of polypropylene and shall only be used with non PP-aggressive media!

## 1.3 General instructions


In case the burette remained for longer period screwed to the flask and was not in regular use, we recommend to purge the instrument prior to again titration or dosing.

**Caution !!** The instrument is not equipped with an „off switch“!

The TITREX will automatically switch off when you do not use it for 3 minutes.

The display of the instrument indicates „OFF“ before switching off. Simply press button „D“ or button „T“ to reactivate it. In the display now appears the last titration or dosing value.

## 2 Commissioning

Put supplied Discharge Tube into the opening and clip-on the protective sleeve (1). Fix the Discharge Tube by pressing it slightly into the protective sleeve (2). Pay attention that the bulge is placed between valve body and fixing bridge of the protective sleeve (3).

### Installation of suction tube

Cut the tube to desired length and out it into the smaller socket at the bottom part of the instrument. Check for seating

### Installation of reverse dosing tube

Cut the reverse dosing tube to the desired length, too (we recommend to cut on the same length as the suction tube to avoid the building of bubbles and foam). Afterwards, put the tube on the bigger socket of the instrument.

## 3 Functional characteristics

### 3.1 Switch on

Switch on the instrument by pressing button „D“ or „T“. During pressing „dos“ or „tit“ appears on the display.

### 3.2 Air-purging

Turn the selector switch to position Purging (90°-turn of ejection cannula) and move handwheels several times clockwise and anticlockwise until no air bubbles appear in the reverse dosing tube. Observe that the piston during purging is moved again to the lower stop. Suck few amounts of liquids (turn handwheels upwards) and put selector switch to position discharge Turn handwheels upwards (discharge) until the discharge tube is completely filled with liquid. Then again turn selector switch to DOS/TIT for dosing or titrating and re-dose remaining liquid in the draining position.

### 3.3 Titration


The titration mode measures the emitted volume in 0,1 ml steps. Turn the Discharge Tube in discharge position.

Now you can titrate 99,99 ml without further knob-pushing.

Of you reached the 99,99 ml point, you have to set the display "00,00" before you can measure titration beyond. Push button "T" and shortly after button "Res". Now you can continue with the titration.

You can aspirate new liquid during titration as oft as necessary. If you want to aspirate certain plenty of liquid you can measure by pressing button „D“ (For more information, see "Dosing").

Press „T“ again before next titration.

### 3.4 Dosing

The dosing mode measures the aspirated volume in 0,1 ml steps.

During dosing mode a dot is blinking between 1st and 2nd digit.

Firstly, turn the plunger completely down and set display "00,00". Push therefore button "D" and shortly after button "Res".

Now aspirate your desired volume and then turn the plunger completely down to eject the measured volume.

## 3.5 Draining

To drain the discharge tube without losses turn the selector switch to position draining (180°-turn of ejection cannula). The remaining liquid in the discharge tube flows back into the flask.

This position avoids accidental titrating or dosing.

## 3.6 Charging

If necessary, charge the built-in NiMH batteries with the supplied charging unit. First put the charging unit into the outlet and then the plug into the bush at the back of the burette.

Please have the batteries changed after 3 years.


# 4 Titrexsoft

If you want to have best working results directly on your computer please use our Titrexsoft 2 software which gives you full access and control to the memory of the digital burette Titrex and more features:

- supports all manual titration applications
- automatic graph creation
- protocol export function.

Please use the manual for Titrexsoft 2 from your scope of supply!

**NEW DEVELOPED**

The screenshot displays the Titrexsoft 2.0 software interface. At the top, the DOS value is 46.44 and the TIT value is 40.77. The interface includes a menu bar with options like File, Reports, and Device. Below the menu, there are input fields for Name, Date, and Y-Axis. The main area is divided into a 'Manual Mode' section with a table of data and a 'Selection graphic' section with a line graph.

Delete	ID	X	Y	Time	Comment
<input type="checkbox"/>	0	2.00	0.5000		
<input type="checkbox"/>	1	3.00	0.6000		
<input type="checkbox"/>	2	4.00	0.7000		
<input type="checkbox"/>	3	5.00	0.7500		
<input type="checkbox"/>	4	6.00	0.7600		
<input type="checkbox"/>	5	7.00	0.7700		
<input type="checkbox"/>	6	8.00	0.7890		
<input type="checkbox"/>	7	9.00	0.8120		
<input type="checkbox"/>	8	10.00	0.8300		
<input type="checkbox"/>	9				

The graph shows a line plot of Y values against X values, with the Y-axis ranging from 0 to 1.3 and the X-axis from 2 to 10. The data points from the table are plotted, showing a curve that rises and then levels off.

## 5 Calibration

The instrument is factory-calibrated to 20°C aqueous solutions. For solutions of varying density and viscosity, the TITREX can be adjusted by the user in case of changed conditions, provided that an analytical balance with a resolution of  $\leq 1$  mg is available.

### 5.1 Calibration without software and PC


- Weight out 5 titrations of nominal volume (50.00 ml) and calculate the volume with temperature correction factor (order our justification procedure for Titrex).
- Average the 5 resultant values.
- Depress and hold button “RES” until in the display appears „CAL“, release button and in the display the dots flashes between the first and second digit as well as between the third and fourth digit.
- Adjust with „D“ and „T“ the calculated average .
- To activate the calibration value press button “MEM” .
- Press button „RES to return to operation mode.
- During the work with this calibration value, in the display flashes the dot between 2nd and 3rd digit.
- In order to return to the factory settings, depress button „RES” until in the display „CAL“ appears, then deactivate the calibration value by pressing the button “MEM” and return to operation mode via pressing „RES”.

# 6 Cleaning and maintenance

## 6.1 Disassembly/Assembly

- (8) Suction valve
- (9) Valve star
- (10) Valve ball
- (11) Ejection valve

Place the valve ball (10) and the valve star (9) into the suction valve (8). Check that the spikes of the valve star (9) point upwards.

## 6.2 Cleaning

The instrument should be cleaned regularly but always when

- changing reagent,
- storing for longer terms,
- using crystallizing media.

### Warning !!

**When cleaning observe all safety instructions. Components of the instrument are contaminated with titrating liquid! Adequate protective clothing, safety glasses and gloves must be worn!**

- Drain instrument and place it along with the bottle inside a sink.
- Unscrew instrument from the bottle and pull out suction tube out of liquid.
- Demount tube, let it drain and mount it again or tap carefully at the inner flask for draining.
- Turn hand wheels like at air purging short up and down until there comes only a modicum of liquid.
- Finally rinse TITREX with distilled water or ethanol.

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**witeg Labortechnik GmbH**