HT-AX

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Contents

1 Introduction	8
2 Safety Instructions	8
2.1 Description of Symbols	8
2.2 Instructions	8
3 Product description	9
3.1 Key features.....	9
3.2 General survey	10
4 Packing contents.....	10
5 Positioning the device and commissioning.....	10
6 Job description	11
7 Maintenance and cleaning.....	11
8 Technical data	11
9 Certificate.....	12

1 Introduction

Thank you for purchasing our **HT-AX Overhead Stirrer**. You did a good choice buying a **Wisd** product.

Before use, please read this manual and use it properly.

To cover the user's health and avoid the product's destruction you have to follow the safety instructions of the second chapter.

Please keep this operating manual in a save place for future reference and hand it out if you consign the device to a third person.

Witeg reserves to modify the capsule or make amendments for the improvement of quality and efficiency without prior notice.

2 Safety Instructions

2.1 Description of Symbols

**WARNING**

Information for avoiding injuries or fatal accidents.

**ATTENTION**

Safety instruction to avoid a damage of the product.

**ATTENTION**

Information about technical parameters.

**Environment**

Information about the operating environment.

**WARNING of explosions**

Safety instruction to avoid explosions.

2.2 Instructions

For a safety and faultless using of the device it is indispensable that you read the following instructions purposely and follow them by the use.

- Connect the device only to the suitable power supply. Please have a look to the identification plate.
- Use the device only in electric circuits with ground fault circuit interrupter.
- Check the power supply cord before you stick it into the electrical outlet.
- Do not use power bars to extend the power supply cord.
- Only use extension cables with sufficient cable cross section.

- Do not touch the device with wet hands during the usage because this might cause an electrical shock.
- If unexpected sound, smell or smog is generated by the device, pull out the main plug and contact the manufacturer or your supplier immediately.
- Disconnect the power supply plug before cleaning or maintenance.
- If you use radiation or contaminated samples the responsibility is only up to you.

- Use this device only if you read the whole safety instructions.
- This device should only be used by persons who are familiarized with safety instructions made for laboratories.
- Pay attention that no liquid reaches the controller or inward of the capsule.
- You have to keep the device clean if you want to get a longer life of it.
- Wear corresponding protective work clothing during the usage.
- Pull out the power cable if you didn't want to use the device for a longer time.

- Do not use this device in areas with potentially explosive atmosphere.

- Make sure that the device can not slip away.
- The device could get damaged if the vibrations are too strong.
- To avoid heat congestions or fires make sure that there is enough space for air circulations.
- Avoid direct solar radiation.
- Do not use the device near highly combustible materials.
- Do not set the device outdoors.
- It is superiorly to set the device in a good illuminated place to avoid accidents engendered by handling errors.
- To set the device in wet or dusty places could cause overheating, short circuits or fires.

- All circuits and electrotechnical components are **witeg** owned technology. In fact of this it is necessary that only **witeg** authorized personal repair the device. Contraventions exclude the warranty.

3 Product description


3.1 Key features

The **HT-AX Overhead Stirrer** is designed for chemical, biological, medical and pharmaceutical applications. The device is specially suited for low viscous liquids getting high viscous during the stirring process.

Key features are as follows:

- Long life by simply and solid construction.
- Constant speed in spite of viscosity fluctuation.
- Motor provided with overheat protection.
- Stepless variable adjustability of every stirrer concerning to the continuous bored chuck.
- Smooth increase of speed for higher safety.

3.2 General survey

4 Packing contents

quantity	DH.WOS01251	DH.WOS01271	DH.WOS01252	DH.WOS01272
1	HT-50AX	HT-120AX	HT-50AX	HT-120AX
1	Allen key			
1	body holder			
1	chuck key			
1	operating manual			
1			KA1177 stand	
1			CL200 clamp	
1			PL020 stirring tool	

5 Positioning the device and commissioning

1. Follow especially the safety instructions marked with following symbol:
2. Remove the package and perhaps soiling carefully.
3. If there are any unexpected transport damages, please inform your forwarder and your costumer service.
4. Attach the clamp to the tripod rod.
5. Stick the device at the clamp and tighten the screw.
6. Check that each connection fits close.
7. Fix the stirrer in the chuck with included chuck key.
8. Before you connect the device with the electrical outlet, please let some time laps until the device is acclimatised.

9. Make sure that you follow all safety instructions.
10. Connect the device with the power outlet.

6 Job description

Attention!

When knob for speed adjustment is not set to "0", the impeller starts turning after switching the device on.

1. Switch the device on.
2. Adjust speed.

7 Maintenance and cleaning

- Please read the whole safety instructions before you start to maintain or clean the device.
- Pull the power plug out of the electrical outlet.
- Only use care cleaning agents and soft cleaning rags.
- Organic lotions, strong chemicals and rough cleaning rags can damage the device.
- If you not want to use the device for a long time, pull of the power plug and keep it dry in the package.
- Only use original spare parts.

8 Technical data

		HT-50AX	HT-120AX
Stirring capacity (H₂O)		Max. 40 lit.	Max.60 lit.
Max. viscosity (mPas)		100.000	150.000
Motor rating (Input / Output)		83 / 75 W	130 / 120 W
Maximal torque (Ncm)		300	500
Max. Shaft diameter (mm)		10	
Speed range (U/min)		0 - 1.000	
Operating environment	Ambient temp. (°C)	0 ~ 40	
	Max. perm. humidity	80 %	
Safety device		Over heat protection	
Dimension (BxLxH mm)		90x180x130	90x210x130
Weight (kg)		3,5	4
Other		- smooth start and change of stirring speed - stable running in spite of viscosity fluctuation	
Power supply		AC 230 V, 1N~, 50/60 Hz	

9 Certificate

Certificate

for

The quality and all features were checked by the manufacturer before the shipment.

We grant from date of purchase

two years guaranty.

This certificate excludes damages by natural disasters or incorrect usages by the costumer.

Please look on your account and complete following table:

Article	HT-AX Overhead Stirrer
Typ	
Serialno.	
Date	

witeg Labortechnik GmbH

Am Bildacker 16
D-97877 Wertheim
TEL: +49-9342-9301-0
FAX: +49-9342-9301-77
Email: info@witeg.de
www.witeg.de