LABmax airless

Flaschenaufsatz-Dispenser
Bottle-top dispenser

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GL 32/25 : Art. Nr. 5.377.025

Index

1 Safety- and user instructions	17
1.1 General safety instructions	17
1.2 Warranty	17
1.3 General product description	18
1.4 To be considered before initial operation	18
2 Commissioning	19
3 Functional characteristics	20
3.1 Air-pruiging	20
3.2 Volume adjustment	20
3.3 Dispensing	20
3.4 Dispensing lock	20
4 Cleaning	21
4.1 Description of the components	21
4.2 Cleaning	22
4.3 Sterilization	23
4.4 Disassembly	24
4.5 Assembly	25
4.6 Maintenance	26
5 Accessories	26

1 Safety- and user instructions

1.1 General safety instructions

Highest safety precautions should be used when dispensing corrosive, poisonous, radioactive or hazardous chemicals. Study these instructions carefully!

- General operation and safety instructions must be followed. Eye protection, protective clothing and glove must be worn.
- The instrument must only be used for intended purpose and utilized the component tolerances of material resistance. If unsure about dispensing a specific solution, please contact the manufacturer.
- Always work with the instrument so as not to endanger the operator or other persons.
- Never use force. Breakage of any part will lead to hazardous exposure for the user as well as others.
- Always check the dispenser for leaks before using. Make sure that the suction and ejection cannula are firmly in position.
- Do not apply force to the instrument.
- Do not drop the instrument and do not immerse the instrument in liquids.
- If the valve is faulty the cannula might drip. Only use original spares and accessories.
- Repairs must only be carried out by the manufacturer's trained personnel. Do not open the instrument. Infringements will void your warranty!
- Do not use the instrument in explosive environments.
- Do not expose the instrument to temperatures exceeding 50°C. Do not autoclave!
- Keep the instrument away from fire or flames.
- Check all screw fittings approximately one hour after each assembly/disassembly if they are fitted tightly, because leakages might occur due to material expansion conditional upon temperature.
- Never use damaged or deformed cannulas.
- The ejection cannula (6) should never face the user and a collection vessel should be placed underneath.

1.2 Warranty

We only consider ourselves responsible in regard to safety, reliability and performance of the device if the dispenser has been used in accordance to the instruction manual and only the parts supplied by the manufacturer have been used. Witeg warrants the Labmx airless for **12 months** from the date of purchase.

1.2 General product description

Our new dispenser-generation "Labmax airless" has an automatic air-purging system which warrants a completely stress-free handling of the media.

Other than with conventional devices there is no return flow of the media in to the bottle during air-purging, therefore, no air bubbles get into the media at reduced filling level.

Through the newly designed canal system the downward moving plunger always displaces the air first before dispensing the liquid. By that, the amount of substance used for air-purging is drastically reduced.

The air-purging of the "Labmax airless" takes place automatically at normal dispense mode so that no manual air-purging, such as opening and closing of the valve cap is necessary. Another prominent feature is the patented cannula reflux dispensing which prevents unnecessary loss of media and serves as an additional safety lock. With that, the use of a closing cap is no longer needed.

With the adapter, the device fits on a variety of bottles.


1. To be considered before initial operation

- Check that the device has not been damaged in transit.
- Push the suction tube (20) on tightly.
- Do not carry the device at the outer housing (2).
- Push ejection cannula (6) on tightly and slide the cannula holder (7) into the formed grooves.
- When screwing on the bottle, do not grip the device at its outer housing (2) Hold the housing at the bottom screw cap (as well during unscrewing).
- Do not pump before the device has been assembled completely and a collecting vessel has been placed underneath.

2 Commissioning

First push the ejection cannula into the valve head and then slide the cannula holder from above into the tongue.

Make sure it fits tightly!

Push the reflux tube (21) into the big socket and the suction tube (20) into the small socket. Control the length of the suction tube. Now you can screw the "Labmax airless" onto the bottle.


Make sure it fits tightly!

3 Functional characteristics

3.1 Air-purging

Before starting to pump, it is essential to place a collecting vessel under the cannula. The valve axle has to be in the dispensing position (arrow up).

Now perform short pumping movements until there are no more air bubbles. Perform constant and smooth movements and make sure that the outer housing (2) is always pushed all the way down.

3.2 Volume adjustment

The quick-volume-adjustment is performed with the adjusting switch. This assures a maximum of exact reproduction of liquid quantities. To set the required dispensing quantity, turn the calibration screw counter clockwise until the adjusting switch can be moved easily. Tighten the calibration screw clockwise after the indicator has been set to the required value.

3.3 Dispensing


Avoid personal injury from chemicals. To protect yourself wear eye protection and use appropriate safety equipment and clothing.

Please follow all safety instructions and observe operating procedures.

Raise the plunger until it stops at the point set by the adjusting switch. By pressing down the plunger to the lowest point the set quantity will be dispensed into the collecting vessel. The plunger movements should be smooth and constant to achieve an exact ejection volume.


3.4 Dispensing lock


The ejection cannula of the "Labmax airless" can be drained safely and easily. Turn the valve axle 180° to the left until it stops (arrow downward). This way it is possible to lock the dispenser which prevents any further drop out. The closing caps are no longer needed. If slightly tilted backwards the ejection cannula drains into the bottle.

4 Cleaning

4.1 Description of the components

4.2 Cleaning

1. Follow the safety instructions!
2. Push the plunger downwards until it stops and then turn the valve axle (12) to the draining position.
3. Let the remaining reagent in the ejection tube (6) flow back into the reagent bottle.
4. Place the device with the bottle in a suitable sink.
5. Unscrew the device from the bottle and slightly tap the inside of the bottle with the suction tube (20) in order to drain it.

Disassemble the "Labmax airless" until you can securely clean it from all reagent build up.

The cleaning can be performed with distilled water and alcohol.

4.3 Sterilization

After removal of reflux tube (21) and suction tube (20) the dispenser can be steam-sterilized (121 °C, 2 bar, 15 minutes) according to DIN EN 285. Place the device on a cloth and avoid any contacts with hot metal surfaces. The discharge tube (6) has to be attached to the protective sleeve (7).

In order to prevent loss of adjustment due to heat expansion of the different materials, the quick lock knob (3) has to be set to minimum 2/10 of its maximum volume. Before you use the dispenser, let it reach room temperature (about 2 hours cooling time). After autoclaving, check all screw fittings for tightness.

After autoclaving please check all parts for deformations or other changes. Do not use any deformed or leaking parts. In case of deformation, please return the device including the autoclave protocols. After autoclaving the dispenser must be checked and recalibrated if necessary.


4.4 Disassembly

1. Follow the safety instructions!
2. Push the glass plunger (17) downwards until it stops and then turn the valve axle (12) to the draining position.
3. Let the remaining reagent in the ejection cannula (6) flow back into the reagent bottle.
4. Place the device with the bottle into a suitable sink.
5. Unscrew the device from the bottle and carefully tap with the suction tube (20) from the inside against the bottle, so that it is drained as well.
6. Rinse the "Labmax airless" with distilled water.
7. Remove the suction tube (20) and the reflux tube (21).
8. Remove the cannula holder (7) and then the cannula (6).
9. Loosen the fix screw (1).
10. Pull the glass plunger (17) out of the PTFE cylinder (18). Before this pump once, then the glass plunger can be handled easier.
11. Unscrew the adaptor ring (4) from the outer housing (2).
12. Remove the calibration screw (3.1) and then the remaining parts of the volume rocker switch (3.2 and 3.3).
13. Remove the outer housing (2).
14. Unscrew the seal nut (15) and place the PTFE seal ring (16) aside.
15. Pull off the adaptor ring (4) and loosen the cap nut (5).
16. Pull the PTFE cylinder (18) out of the valve head (19), keep in mind that the valve star (9) and the ceramic ball (10) might drop out.
17. Take the valve star (9) and the ceramic ball (10) out of the suction valve.
18. Loosen the closing cap (13) and remove the O-ring (14 a).
19. Pull the valve axle (12) out of the valve head (19).

4.5 Assembly

For a problem-free functionality please observe the tight fit between valve head (19) and cap nut (5). If necessary please retighten.

1. Push the valve axle (12) into the valve head (19).
2. Assemble the O-ring (14 a) and the closing cap (13).
3. Place the ceramic ball (10) and the valve star (9) in the suction valve. Check that the spikes of the valve star (9) point upwards.
4. Assemble the PTFE cylinder (18) on the valve head (19). Check that the notches of the PTFE cylinder (18) are placed exactly over the domes of the valve head (19).
5. Tighten the PTFE cylinder (18) with the cap nut (5). Check that all parts fit tightly.
6. Position the adaptor ring (4) with the thread upwards over the cap nut (5).
7. Position the PTFE seal ring (16) with the groove upwards on the PTFE cylinder (18) and screw the seal nut (15) tightly on the thread.
8. Place the device horizontally and perform the following steps consecutively. Prior to that slide the outer housing (2) onto the PTFE cylinder (18).

9. Place the device upright and push the glass plunger (17) into the PTFE cylinder (18) until it stops.
10. Screw the fix screw (1) into the outer housing (2).
11. Push the ejection cannula (6) into the valve axle (12) and fix it with the cannula holder (7).
12. Push the reflux tube (21) and the suction tube (20) into the valve head (19).
13. Screw the device onto the reagent bottle.

4.6 Maintenance

When the "Labmax airless" is screwed open, make sure that the valves are surrounded by liquid.

You can keep the dispenser smooth-running by rinsing it with distilled water or alcohol.

In order to ensure that the device meets its accuracy, from time to time check-ups can be done with (2x)-distilled water on semi-micro-scales. Mind temperature dependence ($1 \text{ ml water at } 20 \text{ }^{\circ}\text{C} \approx 0,998\text{g}$).

5 Accessories

Accessories	Specifications	Art.-No.
Fix screw(1)	2,5 - 10 ml	5 370 251
	25 - 100 ml	5 370 252
Outer Housing (2)	2,5 ml	5 370 83001
	5 ml	5 370 83002
	10 ml	5 370 83003
	25 ml	5 370 83004
	50 ml	5 370 83005
	100 ml	5 370 83006
Adjusting switch (3)	2,5 - 10 ml	5 370 83300
	25 - 100 ml	5 370 83300
Adaptor ring (4)	2,5 & 5 ml	5 370 83231
	10 ml	5 370 83233
	25 ml	5 370 83234
	50 ml	5 370 83235
	100 ml	5 370 83236

Accessories	Specifications	Art.-No.
Cap nut (5)	2,5 ml	5 370 83241
	5 ml	5 370 83242
	10 ml	5 370 83243
	25 ml	5 370 83244
	50 ml	5 370 83245
	100 ml	5 370 83246
Ejections cannula (6) + cannula holder-set (7)	2,5 - 10 ml	5 375 003
	25 - 100 ml	5 375 004
Valve star (9)	2,5 - 10 ml	5 370 377
	25 - 100 ml	5 370 378
Ceramic ball (10)	2,5 - 10 ml	5 370 373
	25 - 100 ml	5 370 374
Valve axle (12)	2,5 - 10 ml	5 370 200
	25 - 100 ml	5 370 201
Closing cap (13)	2,5 - 100 ml	5 370 380
O-Ring a (14 a)	2,5 - 100 ml	5 370 381
O-Ring b (14 b)	2,5 - 10ml	5 370 387
Seal nut (15)	2,5 ml	5 370 83211
	5 ml	5 370 83212
	10 ml	5 370 83213
	25 ml	5 370 83214
	50 ml	5 370 83215
	100 ml	5 370 83216
PTFE seal ring (16)	2,5 ml	5 370 83221
	5 ml	5 370 83222
	10 ml	5 370 83223
	25 ml	5 370 83224
	50 ml	5 370 83225
	100 ml	5 370 83226

Accessories	Specifications	Art.-No.
Glass plunger (17)	2,5 ml	5 370 83101
	5 ml	5 370 83102
	10 ml	5 370 83103
	25 ml	5 370 83104
	50 ml	5 370 83105
	100 ml	5 370 83106
PTFE cylinder (18)	2,5 ml	5 370 83251
	5 ml	5 370 83252
	10 ml	5 370 83253
	25 ml	5 370 83254
	50 ml	5 370 83255
	100 ml	5 370 83256
Valve head (19) -Set (incl. 9,10,12,13,14 a+b)	2,5 - 10 ml	5 370 360
	25 - 100 ml	5 370 361
Suction tube (20) 310 mm	2,5 - 10 ml	5 375 000 N 310
	25 - 100 ml	5 375 001 310
Suction tube (20) 390 mm	2,5 - 10 ml	5 375 000 N 390
	25 - 100 ml	5 375 001 390
Reflux tube (21)	2,5 - 10 ml	5 375 000 N 40
	25 - 100 ml	5 375 001 40

6 Optional accessories

1. Borosilicate glass one-way valve with plastic coating, avoids emitting of gases and strong smells

Article-No. 5.377.300

2. Borosilicate glass tube with plastic coating for calcium chloride, ensures dry atmosphere in storage vessel (delivery without Calcium chloride)

Article-No. 5.377.310

3. Filter against dust and bacteria in storage vessel, pore size 2 μm / 4,5 μm

Article-No. 8.214.002 / 8.214.001

4. Spare adapters

GL 45/32 : Art. No. 5.497.600

GL 32/45 : Art. No. 5.377.045

GL 32/40 : Art. No. 5.377.040

GL 32/38 : Art. No. 5.377.038

GL 32/28 : Art. No. 5.377.028

GL 32/25 : Art. No. 5.377.025

Liquid Handling

LABmax airless

witeg Labortechnik GmbH
Am Bildacker 16
97877 Wertheim
Germany

Liquid Handling

DE-M

www.witeg.de

Phone.: +49 9342 / 9301-0

Fax: +49 9342 / 9301-77

www.witeg.de

info@witeg.de

5 371 XXX_2016-06-16