
SMHS

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Contents

1 Introduction	12
2 Safety Instructions	12
2.1 Description of Symbols	12
2.2 Instructions	12
3 Product description	14
3.1 Key features.....	14
3.2 General survey	14
3.3 Operating surface.....	15
3.4 Description of components	15
4 Packing contents.....	15
4.1 Standard scope of delivery.....	15
4.2 Additional Accessories	15
5 Positioning the device and commissioning.....	16
6 Job description	16
6.1 Menu structure	16
6.3 Set temperature.....	17
6.4 Justification	17
6.5 Set delay time	17
6.6 Key lock	18
6.7 Storage function.....	18
7 Maintenance and cleaning.....	18
8 Troubleshooting Guide	19
9 Technical data.....	20
10 Certificate.....	21

1 Introduction

Thank you for purchasing our **SMHS Multi-Point Hotplate Stirrer**. You did a good choice buying a product from witeg Labortechnik.

Before use, please read this manual and use it properly.

To cover the user's health and avoid the product's destruction you have to follow the safety instructions of the second chapter.

Please keep this operating manual in a save place for future reference and hand it out if you consign the device to a third person.

Witeg reserves to modify the capsule or make amendments for the improvement of quality and efficiency without prior notice.

2 Safety Instructions

2.1 Description of Symbols

**WARNING**

Information for avoiding injuries or fatal accidents.

**ATTENTION**

Safety instruction to avoid a damage of the product.

**ATTENTION**

Information about technical parameters.

**Environment**

Information about the operating environment.

**WARNING of explosions**

Safety instruction to avoid explosions.

**ATTENTION to high temperatures**

Instructions to avoid injuries caused by high temperatures.

2.2 Instructions

For a safety and faultless using of the device it is indispensable that you read the following instructions purposely and follow them by the use.

- Connect the device only to the suitable power supply. Please have a look to the identification plate.
- Use the device only in electric circuits with ground fault circuit interrupter.

- Check the power supply cord before you stick it into the electrical outlet.
- Do not use power bars to extend the power supply cord.
- Only use extension cables with sufficient cable cross section.
- Do not touch the device with wet hands during the usage because this might cause an electrical shock.
- If unexpected sound, smell or smog is generated by the device, pull out the main plug and contact the manufacturer or your supplier immediately.
- Disconnect the power supply plug before cleaning or maintenance.
- If you use radiation or contaminated samples the responsibility is only up to you.

- The hotplate could be scorching during and after the usage.
- Do not put heat-sensitive items at the hot plate.
- Warn all attendees if you leave the device unattended.
- Only use adequate sized vessels and avoid an overboil of the liquids.
- Do not clean the scorching hotplate.

- Use this device only if you read the whole safety instructions.
- This device should only be used by persons who are familiarized with safety instructions made for laboratories.
- Pay attention that no liquid reaches the controller or inward of the capsule.
- You have to keep the device clean if you want to get a longer life of it.
- Wear corresponding protective work clothing during the usage.
- Pull out the power cable if you didn't want to use the device for a longer time.

- Make sure that the device can not slip away.
- The device could get damaged if the vibrations are too strong.
- To avoid heat congestions or fires make sure that there is enough space for air circulations.
- Avoid direct solar radiation.
- Do not use the device near highly combustible materials.
- Do not set the device outdoors.
- It is superior to set the device in a good illuminated place to avoid accidents engendered by handling errors.
- To set the device in wet or dusty places could cause overheating, short circuits or fires.

- Reparations and amendments should only be made by **witeg** authorized people. Contraventions exclude the warranty.

- Do not use this device in areas with potentially explosive atmosphere.

3 Product description

3.1 Key features


The **SMHS Multi-Point Hotplate Stirrer** is designed for chemical, biological, medical and pharmaceutical high temperature stirring-applications. The device was developed for simplest usage and a very long life.

Key features are as follows:


- Long life by simply and solid construction.
- Ceramic coated heating plate for high chemical resistance.
- Suitable for most laboratory stirring and heating applications up to 350 °C.
- Motor provided with overload and overheat protection.
- Power bar graph, speed, temperature and timer showed on illuminated LCD.
- Storage function for set time and speed values.
- Turn-off delay from one minute up to over 99 hours in one minute steps.
- Signal sound after end of turn-off delay.
- Smooth increase of speed.
- Key lock against inadvertent changing of set speed and timer values.
- Connection for external temperature sensor available to measure the liquid's temperature directly.
- Integrated socket for sensor stand.

3.2 General survey

(The **SMHS-3** is showed in the next pictures)

3.3 Operating surface

3.4 Description of components

SV: set value (upper line, in key lock menu no SV)

PV: present value, measured by the device (lower line, in key lock menu also upper line)

LCD: liquid crystal display

MODE button

- For keyboard shortcut and change between the menu items.

Jog-Shuttle knob

- Change SVs by turning.
- Push the knob for start and stop or probate and save set values.
- Also used for keyboard shortcut.

Main switch

- On / off switch.

Stirring lamp

- Is blinking in menu items for speed modulation.

Heating lamp

- Shines if the heater is activated.

Keyboard shortcut

- For reaching the sub menu items.

4 Packing contents

4.1 Standard scope of delivery

quantity	DH.WMH03503	DH.WMH03506
1	SMHS-3	SMHS-6
1	operating manual	

4.2 Additional Accessories

order number	denomination
DH.WHP503060	SS500 PT100 temperature sensor
DH.WHP503061	SS510 PTFE-coated PT100 temperature sensor
DH.WHP503011	RD100 stainless steel stand rod, Ø 12mm, length 450mm
DH.WHP503012	CL200 clamp + sensor holder


5 Positioning the device and commissioning

1. Follow especially the safety instructions marked with following symbol:
2. Remove the package carefully.
3. If there are any unexpected transport damages, please inform your forwarder and your costumer service.
4. If you want to use the device without additional accessories jump to point 10 otherwise follow the next points.
5. Attach the stand rod to the socked.
6. Fix the clamp at the stand rod.
7. Attach the sensor holder to the clamp. Have a look at the next picture if you are unsure. Ensure that the bigger drill hole is upside.
8. Stick the sensor into the holder and connect the cable with the device's sensor socked backwards.
9. Check that each connection fits close.
10. Before you connect the device with the electrical outlet, please let some time laps until the device is acclimatised.
11. Make sure that you follow all safety instructions.
12. Connect the device with the power outlet.
13. Switch the device on.
14. Set your values and start the process (look at job description).

6 Job description

6.1 Menu structure

6.2 Set speed

- First menu item after switching the device on is the menu “set speed”. The stirring lamp is blinking. You can recognize the menu at “RPM” behind SV and PV.
- Speed is from 80 rpm till 1500 rpm adjustable in 5 rpm-steps.
- If there is showed “O RPM” in the SV-line, there is no value stored and you can start the stirring process by turning the Jog-Shuttle knob clockwise.

- If "STIRRING" is blinking on the LCD there is a speed value stored. The saved value is showed in the SV you can start the stirring process by pushing the Jog-Shuttle knob.
- In both cases is the speed during the stirring process variable. To change the SV you only have to turn the Jog-Shuttle knob clockwise to increase the speed and counter-clockwise to reduce it.
- To stop the stirring process manually, you have to push the Jog-Shuttle knob again or decrease the speed till "0" by turning the Jog-Shuttle knob counter-clockwise.
- The set speed is increasing gradually, so it could takes some time till the SV is climaxed.

6.3 Set temperature

- When you push the MODE button in the menu item "set speed" you reach the menu item "set temperature". You can recognize this menu at "°C" behind the SV and PV.
- The temperature is adjustable from 25 °C till 350 °C in 0,5 °C-steps.
- Push the MODE button once again if you want to work without the heating function. Then you reach the next menu (see point 7.5) without collimated temperature value.
- Choose your desired temperature by turning the Jog-Shuttle and confirm the set value by pushing the Jog-Shuttle.
- For changing the confirmed value again, you first have to push the Jog-Shuttle knob again. "HEATING" starts to blink at the LCD and the temperature is adjustable again.
- When the heating lamp is shining, the heater is activated.

6.4 Justification

If the temperature measured by you differs from the temperature the device is showing (PV) you have the possibility to compensate the divergence by the offset. So you can balance differences of about ± 20 °C.

The Offset is calculated as follows:

Temperature measured by you – showed Temperature (PV) = Offset

- To reach this menu you have use the keyboard shortcut in the menu item "set temperature".
- Set the calculated value by turning the Jog-Shuttle. The Offset will be added to the PV immediately.
- Save the offset by holding the Jog-Shuttle pushed for two seconds until the values in PV and SV line starts blinking.
- Push the MODE button to leave the justify menu.

6.5 Set delay time

- The delay timer serves as automatically stop function. It turns off the stirring and heating process after elapsing of the timeframe you set.
- To reach the timer menu you have to push the MODE button in the menu item "set temperature". You can recognize this menu by the blinking timer lamp.
- By renewed using of the MODE button you reach the key lock without setting delay time.
- The delay time is adjustable from one minute till over 99 hours (00:01 - 99:59) in one minute steps.
- Turn the Jog-Shuttle knob to set your delay time. Push the Jog-Shuttle if you reached the desired value to confirm it. With the confirmation the timeframe is elapsing.
- If you want to change the value after confirmation, you have to push the Jog-Shuttle once again. The PV initiates blinking and you can modify the value again.
- To confirm the value and start the delay time push the Jog-Shuttle yet again.

6.6 Key lock

- To avoid accidental changing of the set values, push the MODE button in the menu item "set delay time".
- In this menu item can you see the actually temperature in the SV-line and the speed in the PV-line. The function of the Jog-Shuttle knob is deactivated.
- Press the MODE button once more to leaf the key lock. Than you are in the menu item "set speed".

6.7 Storage function

The **SMHS Multi-Point Hotplate Stirrer** includes a storage function for speed and delay values. That is practical if you want to work with the same parameters often. Then you need not to set them in every usage.

- For storing values you only need to hold the Jog-Shuttle knob pushed for two seconds. Successful storing is displayed by blinking of the SV for three times.
- If you saved a speed value, "STIRRING" is blinking and if you saved a temperature value "HEATING" is blinking after the device's restart.
- If you saved the delay time it is not showed at the LCD.
- Every stored parameter is erasable. Turn the value with the Jug-Shuttle till "0" and then hold the knob pressed for two seconds.
- If you want to store a new value it is not necessary to erase the old one.

7 Maintenance and cleaning

- Please read the whole safety instructions before you start to maintain or clean the device and pull the power plug out of the electrical outlet.
- Only use care cleaning agents and soft cleaning rags.
- Organic lotions, strong chemicals and rough cleaning rags can damage the device.
- If you not want to use the device for a long time, pull of the power plug and keep it dry in the package.
- Only use original spare parts.

8 Troubleshooting Guide

Situation	Confirmations and Solutions
Power Failure	✓ Check the supplied voltage is present. ✓ Check the power connection and the power cord for proper connection and integrity. ✓ Check if the MAIN Power Switch is "ON". ✓ Fuse is blown.
Unit Not Heating	✓ Check if the MAIN Power Switch is "ON". ✓ Check if the set temperature (SV) confirmed with the Jog-Dial knob (by turn the Jog-Dial knob cannot change the set temperature – see point 6.3).
Heating, but unit not reaching temperature or temperature is fluctuated	✓ Problem at the sensor or main controller.
Unit Not Stirring	✓ Check if the MAIN Power Switch is "ON". ✓ Check the set RPM.
ERROR 1	✓ Hotplate defect.
ERROR 2	✓ Motor defect.

If other problems arise, or one of the above not solves problems, please contact your official agent or the manufacturer.

Reparations and amendments should only made by [witeg](#) authorized people. Contraventions exclude the warranty.

9 Technical data

		SMHS-3	SMHS-6
Number of plates		3 × 1	3 × 2
Stirring capacity (H ₂ O)		max. 20 l / place	
Motor rating (per plate)		input 9 W / output 4 W	
Speed range (rpm)		80 - 1.500	
Speed control resolution		5 rpm	
Temp.	Range (°C)	max. 350	
	Accuracy	±0.3 °C at set temperature	
Temp. control resolution		0.5 °C	
Temp. display resolution		0.1 °C	
Heating power		400W / Plate (total 1.2 kW)	400W / plate (total 2.4 kW)
Controller		digital fuzzy controller with Jog-Shuttle knob (Turn + Push)	
Display		illuminated LCD	
Timer		99h 59min	
Plate Size (mm)		140×140	
Sockets for external temp. Sensor		3	3
Sockets for Stand rod		2	4
Material	Plate	ceramic coated aluminium	
	Body	powder coated steel	
Dimension (w×d×h mm)		510×240×100	510×410×180
Weight (kg)		10	18
Others		- storage function (temperature and timer), - key lock (Jog-Shuttle input disabled)	
Power supply		AC 230 V, 1N~, 50/60 Hz	

10 Certificate

Certificate

for

The quality and all features were checked by the manufacturer before the shipment.

We grant from date of purchase

two years guaranty.

This certificate excludes damages by natural disasters or incorrect usages by the costumer.

Please look on your account and complete following table:

Article	SMHS Multi-Point Hotplate Stirrer
Typ	
Serialno.	
Date	

witeg Labortechnik GmbH

Am Bildacker 16
D-97877 Wertheim
TEL: +49-9342-9301-0
FAX: +49-9342-9301-77
Email: info@witeg.de
www.witeg.de