SWON

Smart Ofen Bedienungsanleitung

Smart Oven Operating Manual

gilt für / written for
SWON -32/ -50/ -105/ -155/ -W105/ -W155

witeg Labortechnik GmbH

(English manual at page 26)

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Contents

1 Introduction	27
2 Safety Instructions	27
2.1 Description of symbols	27
2.2 Instructions	27
3 Packing contents	29
4 Product description	29
4.1 Key features	29
4.2 General survey	30
4.3 Operating surface	30
4.4 Term description	31
5 Positioning the device and commissioning	31
6 Installation of WiRe™ App & Web Service	31
7 Product Usage	36
7.1 Temperature Setting	36
7.2 Settings	36
7.3 CHART	43
7.4 OFFSET	44
7.5 Error Message.....	45
8 Maintenance and cleaning	46
9 Troubleshooting Guide	47
10 Technical data	48
11 Certificate	49

1 Introduction

Thank you for purchasing our **SWON Smart Oven**. You did a good choice buying a **Wisdom** product.

Before use, please read this manual and use it properly.

To cover the user's health and avoid the product's destruction you have to follow the safety instructions of the second chapter.

Please keep this operating manual in a save place for future reference and hand it out if you consign the device to a third person.

Witeg reserves to modify the capsule or make amendments for the improvement of quality and efficiency without prior notice.

2 Safety Instructions

2.1 Description of symbols

**WARNING**

Information for avoiding injuries or fatal accidents.

**ATTENTION**

Safety instruction to avoid a damage of the product.

**ATTENTION**

Information about technical parameters.

**Environment**

Information about the operating environment.

**WARNING of explosions**

Safety instruction to avoid explosions.

**ATTENTION to high temperatures**

Instructions to avoid injuries caused by high temperatures.

2.2 Instructions

For a safety and faultless using of the device it is indispensable that you read the following instructions purposely and follow them by the use.

- Connect the device only to the suitable power supply. Please have a look to the identification plate.

- Use the device only in electric circuits with ground fault circuit interrupter.
- Check the power supply cord before you stick it into the electrical outlet.
- Do not use power bars to extend the power supply cord.
- Only use extension cables with sufficient cable cross section.
- Do not touch the device with wet hands during the usage because this might cause an electrical shock.
- If unexpected sound, smell or smog is generated by the device, pull out the main plug and contact the manufacturer or your supplier immediately.
- Disconnect the power supply plug before cleaning or maintenance.
- If you use radiation or contaminated samples the responsibility is only up to you.

- Pay attention if you open the door, hot air is coming out.
- Do not put heat-sensitive items at the device.
- Warn all attendees if you leave the device unattended.
- Do not clean the scorching device.

- Repairs and amendments should only be made by **witeg** authorized people. Contraventions exclude the warranty.

- Use this device only if you read the whole safety instructions.
- This device should only be used by persons who are familiarized with safety instructions made for laboratories.
- Pay attention that no liquid reaches the controller or inward of the capsule.
- You have to keep the device clean if you want to get a longer life of it.
- Wear corresponding protective work clothing during the usage.
- All circuits and electrotechnical components are **witeg** owned technology. In fact of this it is necessary that only **witeg** authorized personnel repair the device.
- Pull out the power cable if you don't want to use the device for a longer time.

- It is important that the device stands up on a horizontal and firm subsurface to avoid vibrations.
- If you set the device in a site with vibrations, make sure that the device can not slip away.
- The device could get damaged if the vibrations are too strong.
- To avoid heat congestion or fires make sure that there is enough space for air circulations.
- Avoid direct solar radiation.
- Do not use the device near highly combustible materials.
- Do not set the device outdoors.

- It is superiorly to set the device in a good illuminated place to avoid accidents engendered by handling errors.
- To set the device in wet or dusty places could cause overheating, short circuits or fires.

- Do not use this device in area with potentially explosive atmosphere.

3 Packing contents

quantum	denomination
1	SWON Smart Oven
2	stainless Steel Shelves
1	operating manual
1	connection cable

4 Product description

4.1 Key features


The **Wisdom SWON Smart Oven** are outfitted with high performance heating mechanism that is optimized respectively for capacity of chamber, power of heating element and air circulation type (Gravity Convection or Forced Convection). The temperature displayed on the controller is accurately consistent with the actual temperature in the very center of the chamber or the average temperature in the entire area of the chamber, which maintains precisely within fluctuation specification. In addition, the temperature in the entire area of the chamber shows a very high stability keeping within a very slight variation specification in terms of temperature uniformity. Simultaneously, through the optimization of Fuzzy-PID controller, the 'heat-up time' to reach the Set Temperature after operating and the 'recovery time' to get back to the set temperature after opening and closing of the door are very short. The newly developed high performance heating mechanism has achieved high degree of temperature accuracy and uniformity in the chamber. This is achieved by the heating elements installed on three sides - that is, the left, right and bottom parts of a chamber - which directly heat the stainless-steel chamber. This either enhances its material thermal conductivity to maximum (in case of the Incubator) or the heated air from the pre-heating zone in air-jacket type that is built in the three sides outside a chamber circulates first.

Key features are as follows:

- Suitable for Drying, Baking, Conditioning, Curing, Pre-Heating and Aging.
- The Best Temp. Uniformity & Accuracy by High Performance Heating Mechanism : 3 Side-Heating.
- Optimized Air Flow by Gravity Convection.
- Ergonomic Design 4" Full Touch Screen TFT LCD (Smart-Lab™ Controller).
- WiRe™ App & Web Service.
- Automatic Data Recording.
- Self-Diagnostic Function.
- Compact New Body Design.
- Variable Program Setting for Multiple Patterns with Multiple Segments.
- Internet Connectivity with WiFi.
- Controller SW can be upgraded through USB port.

- The Recorded Data could be transferred to PC by USB memory.
- Corrosion Resistant 304 Stainless Steel Chamber.
- Over Temperature & Over Current Protection, Sensor Error Detection.

4.2 General survey

Main power

- On / off switch.

Safety switch

- The safety switch turns the heating off if your set maximal temperature is overshoot.
- Set the maximal temperature rough 30 – 40°C over the desired maximal temperature, otherwise the device maybe not reaches your set value.

Fuse

- Protect unit from High Current

RS232C Connector

- Monitoring and Controlling with PC by connecting RS232C cable.

4.3 Operating surface

4.4 Term description

NO	Names	Description
①	SET TEMP	Adjust and check the setting temperature.
②	CURR TEMP	Display the current temperature.
③	START/ STOP	Controls operation starts and stops
④	Program Mode control	Sets programs and controls operations
⑤	CHART	Monitors operational status in graph
⑥	SETTINGS	Sets various functions of device
⑦	Error	Notifies device's error
⑧	WiFi connectivity	Indicates internet connectivity
⑨	Date	Indicates current date, time, operation hours
⑩	Operation Information	Indicates operation information
⑪	Operation ON/OFF ICON	Indicates ON/OFF in operation with icons
⑫	OFF SET	Set compensating temperature values


	Heater ON		Heater OFF		In operation		Stopping
---	-----------	---	------------	---	--------------	---	----------

5 Positioning the device and commissioning

1. Follow especially the safety instructions marked with following symbol: 
2. Remove the package carefully.
3. If there are any unexpected transport damages, please inform your forwarder and your customer service.
4. Before you connect the device with the electrical outlet, please let some time laps until the device is acclimatised.
5. Make sure that you follow all safety instructions.
6. Connect the device with the power outlet.
7. Switch the device on.
8. Set your values and start the process (look at job description).

6 Installation of WiRe™ App & Web Service

- After installing this product in a static position, the user who uses WiRe™ App & Web Service, please use this product after installing WiRe Service before using the product installed.
- Proceed with the installation from the Touch Screen in the front.
- Because the Wireless LAN for using the internet is not installed in the product when shipping it at the factory, please connect by finding the wireless LAN among the components of the product.

- This description is for the installation method for users who want to use WiRe™ App & Web Service. Only matters that should be installed when use WiRe Service are written,

and the detailed description on the product is listed in detail on "Instruction Manual".

- The sequence of the setting of WiRe™ Service during installation is as described in the following instructions. Please follow the steps in order as described.
- Find the "SETTING" located on top-left of the main display to make settings required for the connection of WiRe™ Service.
- "SETTING" menu is disabled during operation. Stop operating it by touching "STOP" button on the top right and check if the SETTING menu is activated.

- The "SETTING" has submenus of "System", "WiRe", and "WiRe Users" prepared for the connection to WiRe™ Service. Please follow following numbers of each instruction in order.


- The initial page of the description screen is the same as follows, and to press "BACK" button of the left top when the set-up is completed is switched back to the main screen.
- This description is listed in order to be set when install WiRe™ App & Web Service. Please progress in order of this description.
- First, the item that you must set is Current Time. Please fill in the correct Current Time, before that the product operates through Time Set Menu.

① **System Settings:**

1) Current Time Set


- The present time is set to input correct test time of user. Touch the box displaying time and date on the right side of the display.
- If you randomly change the current time since supplied from the factory, please be sure to use after re-applying power to the equipment. If you have changed the time of product behavior it may also cause problems with some functions. Therefore if the current time changed be sure to use the device again after turning it off.
- Please select the “Automatic” with the Check Box unchecked.
- The “Select time zone” will synchronize the product with the time of selected zone (countries) for users using the ‘Mobile App’ to control the WiRe™ Service. Please be sure to select the time zone wherein the product is being used.

- The Check Box on the right prepared for the “Use 24-hour format” must be checked for the correct operation of the product, and users should never change it.
- After completing the setting of the product, push the return button of blue icon to go back to previous page.

2) WiFi Setup Set

- The ‘WiRe™ Service’ will only be available by the (wireless) internet connection. Users should make the product be connected to internet via wireless network.


- Click the ‘Check Box’ of Wi-Fi window to make the Wi-Fi connection (by typing set password if requested). In the case unable to detect Wi-Fi signals, please check the mounting of USB dongle prepared for the Wi-Fi connection (Check the blinking of the blue LED).
- After completing the Wi-Fi connection to the product, push the blue return button located on top-right of the display to go back to previous page.

② WiRe Settings

- Slide the “WiRe Mode” button to the right to connect to the “WiRe™ Service”. This will only be enabled by the stabilized Wi-Fi connection.
- In cases of the unstable or disconnected Wi-Fi connection, the “WiRe Mode” will not be switched normally. Be sure to secure stable Wi-Fi connection to the product.
- After completing the internet connection, users can find the “Serial Number”, “Model”, and “Organization Name” as shown in the following display.

➤ WiRe Users

1) Registration of Additional Users

- The users registered through existing ‘Smart-Lab’ device can use existing IDs to be registered in this product. When the user is logged in with the registered ID from the ‘Mobile App’ then the list of devices already registered will be appearing.


2) Deletion of Users

- Registered users may be deleted from the list by pushing the selected button corresponding to users to be deleted therefrom.



3) Registration of New User

- Users who are using the 'Smart-Lab' device and 'WiRe Service' for the first time shall register users' information. Input the 'ID', 'Password', and 'Name' of the user for the registration.
- The user ID can be composed of minimum 6 to maximum 20 lowercase alphanumeric characters. In cases of the input of duplicated ID or other characters beyond the lowercase alphanumeric characters, the message showing the disability of ID setting will be popped up.
- The Password can comprise minimum 6 to maximum 30 uppercase or lowercase alphanumeric or special characters. The input with less than 6 or beyond 30 characters thereof or the input of characters other than predefined characters will encounter the pop up message showing the disability of Password setting.
- Once the 'ID', 'Password', and 'Name' are set in the list then they cannot be modified thereafter.

4) Setting of Authority of each User

- Users registered in the list can be set as an administrator or as a general user. The 'Admin' (the administrator) is the only one who can entirely administrate the 'Remote Control System'. Up to 100 users can be registered in the list as a 'User' allowed the monitoring of the status of device. (The remote control is not allowed for the 'User'.)
- The 'ID' among the list of registered users is selected to assign an authority set either as 'Admin' (the administrator) or as 'User' (general users).
- Since the authority of 'Admin' was prepared for the only one 'ID' of 'Admin', the existing 'ID' of 'Admin' will be turned into the 'ID' of 'User' with the appointment of new 'ID' of 'Admin'.

- The appearance of the message showing the disability to invoke the information of registered users may come from the unstable Wi-Fi connection. Check the status of Wi-Fi connection and then reenter into the menu of 'WiRe Users'.


- On completion of above procedure, the product will be connected to the 'Mobile App' enabling the service of 'WiRe™ Service'.

7 Product Usage

- Set the SAFETY knob to 30~40°C higher than operation temperature.
- Connect the power cord to the proper power socket and Turn on the Main Switch.
- All the operations should be operated through the touch screen.
- For secure operation, the major function of controller is protected by password. Default password is "0000" and we recommend change it in "Set Up".

7.1 Temperature Setting

- You can set the temperature by entering the number directly on the Temperature Setup screen by pressing "SET" button on the right where the set temperature is shown.

- You can set the temperature between 25°C to 230°C. Press "ENTER" after entering the set temperature.


7.2 Settings

- Functions embodied in this product can be set in the "SETTING" menu. Push the return button located on top-right of all display to go back to previous display.
- The 'SETTING' menu will be deactivated during the operation of the product. Push the "STOP" button located on top-right of the display to stop the operation and to activate the 'SETTING' menu.

- Following are the sequentially arranged instructions prepared for the operation of the product. Please follow the instructions described in order.
- To use the 'WiRe™ Service' of the product, please refer to the "WiRe™ Service" included in the installation manual.

① Operation Settings

(1) Running Mode After Power Failure

- It sets a mode for power is abruptly cut off and turned back up due to power supply issues such as power outage.
- **Stop mode:** It stands ready at the value prior to the power outage and maintains the operation on hold.
- **Restart mode:** It operates depending on (4) Operation Mode.

A) Fixed Mode: When reauthorizing the power, it automatically operates at a temperature set that it previously operated on and the operation time on main screen starts from 0 again.

B) Program Mode: It restarts from the first segment of the previous pattern it operated on.


- **Continue mode:** It operates depending on (4) Operation Mode.

A) Fixed Mode: When reauthorizing the power, it automatically operates at a temperature set that it previously operated on and the operation time on main screen continues its counting.

B) Program Mode: It restarts the count from the remaining time of first segment that stopped and it restarts.



(2) Chart Display Cycle

- You can set the temperature cycle shown on Chart. User can select from 30 seconds, 1 minute and 2 minutes.

(3) Reserve Mode

- Set date and time you want to operate the device and it will start at the reserved time. Minimum time required to use this mode is 5 minutes.
- It automatically turns off after being applied one time.

(4) Operation Mode

- There are two choices: the Fixed Mode that operates within the fixed temperature range and the Program Mode that operates in various temperature in sequence by patterns.

*** Fixed Mode**


- If Fixed is selected, Fixed Operation time Mode window is activated.

- Set Fixed Operation Time Mode to ON and press Fixed Operation Time to set the value of timer. It operates for a time it set when you start the product from main screen.
- If Fixed Operation Time Mode is set to OFF, the set temperature keeps operating. Limit of Operation Time is 3 minutes ~ 23 hours and 59 minutes.

*** Program Mode (in order)**

- If the program is selected, Program Mode window shown below is activated.
- Program Mode can select one of two different methods to reach a target temperature.
- The following describes how to set the product operation. Please refer to (5) Program for how to set Program Mode in detail.

A) Ramp: It can set time to reach a target temperature.


Per one target temperature, it needs a segment time setting to maintain a target temperature after setting the segment time to reach the target temperature. If the maintaining time for target temperature is not needed, do not make second segment for the target temperature.

B) Normal: It operates at the maximum performance to reach a target temperature.

② System Settings

- (1) Buzzer
 - You can select ON/OFF for Error alarm.
- (2) Touch Sound
 - You can select ON/OFF for Touch Button sound.
- (3) Data Initialization
 - It initializes Operation Settings stored and Chart Data.


③ About


- You can check information of the device as Serial Number, Model, Temp Range and a version of Touch App.

④ Security Settings


- (1) Key Lock
 - It requests password to access some key menu settings in menu. Key Lock is activated by setting button to the right.
- (2) Password Change
 - It allows password to change. Enter 4 digits password stored in "Current Password" (factory setting "0000") and type in "New Password" and "Retype Password" to finalize the process.


⑤ **Program (In order mode)**


- You can select and store operation program as pattern and segment.
- You can select and use one of several programs stored when using the product.
- Program is a collection of Pattern and pattern is a collection of Segment.



- Pattern can store up to 120 and a maximum of 100 segments can be stored in 1 Pattern.
- Pattern is created by using + button, and in case you want to delete, press the number of Pattern you want to delete and verify the Pattern number selected in white color and delete using – button. If Pattern is deleted, the Segment stored within the Pattern also is deleted.
- If you want to set Segment within Pattern, select the Pattern number you want to select as Segment, create by using “+” button next to Segment. If you want to delete, press the corresponding Segment number and verify the selection in white color, delete using “–” button.
- Pattern and Segment created can be immediately verified on screen in graph, and can Enter or change values by clicking the graph.
- Enter the value of Segment after creating Pattern and Segment. You can verify the temperature and time set for the Segment selected when you select Segment from the Segment window and touch.

- You can reset temperature by clicking (1) Temperature window of Segment and can reset the timer by clicking (2) Operation time window.

- If OK is pressed after entering the value set, you can immediately verify the changes from graph.
- Segment graph is displayed depending on the operation method selected from ramp or normal from Program Mode of Operation Mode.

- When setting with Program Mode, it is activated when (1) Program menu is created and operation is stopped. (2) STEP and (3) buttons are created in the operation of Program Mode.


- (1) Program button: it enters into the Program setting screen rapidly.
- (2) Step button: It skips to next Segment while in operation.
- (3) Hold button: It temporarily stops while in operation.

⑥ Repeat

You can set to repeat the operation within the selected Pattern by entering Repeat after selecting Pattern from Program Settings screen.

- If Pattern is created from Program Settings, you can set and store Repeat Mode regardless of Operation Mode. However, when operating with Program Mode, you must set Repeat Mode to ON.

(1) You set the number of repetition after verifying the Pattern number of Select Pattern. You can set to repeat up to 100 times, and can be set to “∞” infinity by pressing “-” button at 1.

(2) You can edit the order of operation after setting the number of repeat. Select Start Segment and End Segment by pressing “+” button from Sequence list, and set the number of repetition between the Sequences selected.

(A) Start Segment: Select first Segment for repeat operation

(B) End Segment: Select final Segment for repeat operation


(C) Segment Count: Set the number of repeats from Start Segment till End Segment

- You can enter maximum Sequence of 10


(3) Next Pattern enters Pattern that operates after the repeat operation of selected repeat Pattern is completed.

- Next Pattern can select from stored Pattern other than selected Pattern with Repeat operation, it does not allow to operation in continuation of other Pattern after the Repeat operation if set to none.


⑦ Wait

- In operation of Program Mode, if Segment of selected Pattern has not reached the temperature set within the time set, Wait Time allows to wait until the temperature reaches Temperature Zone (set error of the target temperature ± value) without moving into the next Segment.
- If not reached the temperature set within Remaining Time in Ramp Mode, Wait Time starts.


- Wait time can be set to 23 hours and 59 minutes and Temperature Zone can be set to ±10°C.
- If the target temperature does not reach Temperature Zone within Wait Time, it moves on to the next Segment.


7.3 CHART

- Temperature data is stored at internal memories at every 10 seconds and maximum 6 months of data are stored. Please use Data Backup once every 6 months.
- Temperature scope shown on graph can be changed using the set button on the low left. Maximum and minimum values can be changed within limits shown on a pop up window.

- After inserting a USB into the port on top of controller, it stores in the USB if the store icon on top right is pressed.
- Error message appears if USB is not recognized. Press “OK” button and remove USB and reinsert and press the store icon.

- You can see the detailed information such as date, time, temperature, Error, etc. if you press a point on graph. If you drag right and left or press other points while pressing the point, detailed information of the location emerges on a detail display.
- You can check the information on graph by moving at a small increment of the point using arrows on the detail display. The area showing detail information is displayed in a straight line on the overall graph in yellow vertical lines.


7.4 OFFSET

In case a user must test by matching temperature gauge used as a standard for test with temperature gauge shown on the controller, you can set the value of offset and increase or decrease the temperature displayed on the screen.

- Limit of offset values are -20.0°C ~ $+20.0^{\circ}\text{C}$, and current temperature after setting offset value displays the value that includes the offset.
- Once you set the offset value, the value is stored and the offset value is reflected in the current temperature even you press ON/OFF button.

7.5 Error Message

Wait screen shows a warning Error! in case of failure and it shows which Error it is if you press the icon.

ERROR Message	Cause	Measures
Temp Sensor Error	Informs when temperature sensor fails	Contact A/S Center
Door Open	Informs if door is left open over 10 sec.	Check the door
Control Response Error	Informs communication errors of device and Touch Board occur	Check the connection of Touch Board


* **Touch Screen is rebooted when the reset button is pressed. Press Reset Button by inserting appropriate tools (paper clips) into Reset hole on the right side of Touch Screen as shown below in case Touch Screen fails.**

8 Maintenance and cleaning

- Please read the whole safety instructions before you start to maintain or clean the device.
- If you not want to use the device for a long time, pull of the power plug and keep it dry in the package.
- The Air Duct is controlled by user's needs, and if it is unused, it should be closed the entrance and then stored in order to prevent the external contaminants entering from the entrance.
- Before cleaning the unit, disconnect the power cord from the wall socket. Otherwise, it may cause an electric shock or fire.
- To clean the unit, a neutral detergent and soft cloth is recommended. Do not use coarse cloth, strong chemicals or organic solutions.
- Do not use water, Benzene, Thinner or any alcohol for cleaning the product. It may cause discoloration, damage, an electric shock or fire.
- If you expect damage by a chemical during cleaning, call technical support team.
- Do not pour water directly into the unit. It may cause an electric shock or fire.
- During cleaning of the chamber, be careful cleaning the corners to avoid injury.

Sample loading method

- Do not place samples on the floor of internal area of chamber. It may cause fire and overheat.
- Do not fully pack with sample for an easier air circulation inside the product and maintain a certain distance.
- Do not touch the inside of product or floor surface. It may cause injury if not sufficiently cooled.

9 Troubleshooting Guide

Situation	Confirmations and Solutions
Power Failure	<ul style="list-style-type: none"> ✓ Check the supplied voltage is present. ✓ Check the power connection and the power cord for proper connection and integrity. ✓ Check if the MAIN Power Switch is "ON". ✓ Fuse is blown.
Unit Not Heating	<ul style="list-style-type: none"> ✓ Check if the MAIN Power Switch is "ON". ✓ Check the Safety switch is set higher (min. 30°C-40°C) than the set temperature (SV).
Heating, but unit not reaching temperature or temperature is fluctuated	<ul style="list-style-type: none"> ✓ Check the Set Temperature is higher than ambient temperature (recommended min set temp: 10°C higher than ambient temp). ✓ Check the Safety switch is set higher (min. 30°C-40°C) than the set temperature (SV). ✓ Check the Ventilation openings are closed. ✓ Problem in the sensor or main controller.

If other problems arise, or one of the above not solves problems, please contact your official agent or the manufacturer.

Reparations and amendments should only made by [witeg](#) authorized people. Contraventions exclude the warranty.

10 Technical data

		SWON-32	SWON-50	SWON-105 SWON-W105	SWON-155 SWON-W155
Dimension (wxdxh)	Internal (mm)	310x290x360	370x350x420	485x409x535	550x474x600
	External (mm)	458x557x664	518x647x724	653x715x882	718x777x947
Capacity		32 l	50 l	105 l	155l
Heater		500 W	650 W	1,4 kW	1,6 kW
Power Consumption		502 W	674 W	1390 W	1566 W
Energy Consumption, at 100°C		90 Wh	117 Wh	187 Wh	204 Wh
at 150°C		181 Wh	230 Wh	337 Wh	398 Wh
Temp.	Range	ambient +5 °C up to 230°C			
	Accuracy	±0.5°C at 100°C, ±0.6°C at 150°C			±0.5°C at 100°C ±0.8°C at 150°C
	Uniformity	±1,2°C at 100°C, ±2,2°C at 150°C			
Control resolution		±0.1°C			
Heat-up Time		30min. to 100°C 43min. to 150°C	30min. to 100°C 45min. to 150°C	35min. to 100°C 50min. to 150°C	
Sensor		PT100			
Recovery Time (Door open 30sec)		8min. to 100°C 10min. to 150°C	12min. to 100°C 10min. to 150°C		
Controller		Smart-Lab™ Controller			
Display		4"Full Touch Screen TFT LCD			
Timer		99hr 59min, (Delay & Continuous Function Included)			
Material	Internal	Stainless Steel (#304)			
	External	Powder Coated Steel			
	Insulation	Glass Wool			
Shelves		2 Stainless Steel Shelves			
Ventilation		Stainless Steel Cap			
		Dia. 40mmx1ea		Dia. 40mmx2ea	
Circulation		Gravity Convection Type			
Safety Device		Over Temp. and Over Current Protector, Sensor Error Detector			
Others		Storage Function (Temp., Timer), Alarm (Error Status and Timer-end)			
Net Weight (kg)		36	44	69	78
Packing Size (wxdxh) (mm)		562x664x747	622x754x807	743x819x1.068	808x884x1.133
Shipping Weight (kg)		43	59	79	92
Power Supply		AC 230 V, 1N~, 50/60 Hz			

11 Certificate

Certificate

for

The quality and all features were checked by the manufacturer before the shipment.

We grant from date of purchase

two years guaranty.

This certificate excludes damages by natural disasters or incorrect usages by the costumer.

Please look on your account and complete following table:

Article	SWON Smart Oven
Typ	
Serialno.	
Date	

witeg Labortechnik GmbH

Am Bildacker 16
D-97877 Wertheim
TEL: +49-9342-9301-0
FAX: +49-9342-9301-77
Email: info@witeg.de
www.witeg.de