VM-10

Vortex Mixer

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Contents

1 Introduction	7
2 Safety Instructions	7
2.1 Description of symbols	7
2.2 Instructions	7
3 Product description	8
3.1 Key features	8
3.2 General survey	9
4 Package content.....	9
5 Positioning the device and commissioning.....	9
6 Job description.....	10
7 Maintenance and cleaning.....	10
8 Technical data.....	10
9 Certificate.....	11

1 Introduction

Thank you for purchasing our **VM-10 Vortex Mixer**. You did a good choice buying a **Wisdom** product.

Before use, please read this manual and use it properly.

To cover the user's health and avoid the product's destruction you have to follow the safety instructions of the second chapter.

Please keep this operating manual in a safe place for future reference and hand it out if you consign the device to a third person.

Witeg reserves to modify the capsule or make amendments for the improvement of quality and efficiency without prior notice.

2 Safety Instructions

2.1 Description of symbols

**WARNING**

Information for avoiding injuries or fatal accidents.

**ATTENTION**

Safety instruction to avoid a damage of the product.

**ATTENTION**

Information about technical parameters.

**Environment**

Information about the operating environment.

**WARNING of explosions**

Safety instruction to avoid explosions.

2.2 Instructions

For a safety and faultless using of the device it is indispensable that you read the following instructions purposely and follow them by the use.

- Connect the device only to the suitable power supply. Please have a look to the identification plate.
- Use the device only in electric circuits with ground fault circuit interrupter.
- Check the power supply cord before you stick it into the electrical outlet.
- Do not use power bars to extend the power supply cord.

- Only use extension cables with sufficient cable cross section.
- Do not touch the device with wet hands during the usage because this might cause an electrical shock.
- If unexpected sound, smell or smog is generated by the device, pull out the main plug and contact the manufacturer or your supplier immediately.
- Disconnect the power supply plug before cleaning or maintenance.
- If you use radiation or contaminated samples the responsibility is only up to you.

- Do not use this device in areas with potentially explosive atmosphere.

- Use this device only if you read the whole safety instructions.
- This device should only be used by persons who are familiarized with safety instructions made for laboratories.
- Pay attention that no liquid reaches the controller or inward of the capsule.
- You have to keep the device clean if you want to get a longer life of it.
- Wear corresponding protective work clothing during the usage.
- Pull out the power cable if you didn't want to use the device for a longer time.

- Make sure that the device can not slip away.
- To avoid heat congestions or fires make sure that there is enough space for air circulations.
- Avoid direct solar radiation.
- Do not use the device near highly combustible materials.
- Do not set the device outdoors.
- It is superiorly to set the device in a good illuminated place to avoid accidents engendered by handling errors.
- To set the device in wet or dusty places could cause overheating, short circuits or fires.

- Reparations and amendments should only be made by **witeg** authorized people.
Contraventions exclude the warranty.

3 Product description

3.1 Key features

The **Wisdom VM-10 Vortex Mixer** is designed for chemical, biological, medical and pharmaceutical mixing-applications. The device was developed for simplest usage and a very long life.

Key features are as follows:

- Simplest operating.
- Fast commissioning.
- Long life by simply and solid construction.
- Continuous run and tip operation mode.
- Over heat and over current protection.
- Powerful and quiet operation.

3.2 General survey

4 Package content

quantum	denomination
1	VM-10 Vortex Mixer
1	Power cord
1	Mixing platform
1	Mixing head
1	Operating manual

5 Positioning the device and commissioning

1. Follow especially the safety instructions marked with following symbol:
2. Remove the package carefully.
3. If there are any unexpected transport damages, please inform your forwarder and your customer service.
4. Remove perhaps soiling and package material.
5. Before you connect the device with the electrical outlet, please let some time laps until the device is acclimatised.
6. Make sure that you follow all safety instructions.
7. Connect the device with the power outlet.
8. Switch the device on and start the process (look at job description).

6 Job description

Continuous run: Set the switch “ON” for continuous mixing.

Tip operation mode: Press on the shaking unit to mix.

Mixing intensity: Turn the knob for speed to adjust mixing intensity.

7 Maintenance and cleaning

- Please read the whole safety instructions before you start to maintain or clean the device.
- Pull the power plug out of the electrical outlet.
- Only use care cleaning agents and soft cleaning rags.
- Organic lotions, strong chemicals and rough cleaning rags can damage the device.
- If you not want to use the device for a long time, pull of the power plug and keep it dry in the package.
- Only use original spare parts and accessories.

8 Technical data

	VM-10
Dimension (mm) (wxdxh)	160 x 128 x 145
Shaking motion	Orbital
Speed (rpm)	3300
Weight (kg)	3
Power supply	AC 230 V, 1N~, 50/60 Hz

9 Certificate

Certificate

for

The quality and all features were checked by the manufacturer before the shipment.

We grant from date of purchase

two years guaranty.

This certificate excludes damages by natural disasters or incorrect usages by the costumer.

Please look on your account and complete following table:

Article	VM-10 Vortex Mixer
Serialno.	
Date	

witeg Labortechnik GmbH
Am Bildacker 16
D-97877 Wertheim
TEL: +49-9342-9301-0
FAX: +49-9342-9301-77
Email: info@witeg.de
www.witeg.de