WCR / WCL

Kälte-Umwälzthermostat
Bedienungsanleitung
Refrigerated Bath Circulator
Operating Manual

gilt für / written for
WCR -P8/ -P12/ -P22/ -P30
WCL -P8/ -P12/ -P22/ -P30

witeg Labortechnik GmbH

(English manual at page 13)

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Contents

1 Introduction	14
2 Safety Instructions	14
2.1 Description of symbols.....	14
2.2 Instructions	14
3 Product description	16
3.1 Key features.....	16
3.2 General survey	17
3.3 Operating surface.....	17
3.4 Term description	17
4 Packing contents.....	18
5 Positioning the device and commissioning.....	18
6 Job description	19
6.1 Menu structure	19
6.2 Set turn-off delay.....	19
6.3 Set turn-on delay.....	19
6.4 Set temperature.....	20
6.5 Set Offset	20
6.6 Set Pump Power	20
6.7 Key lock	20
6.8 Storage function.....	20
7 Maintenance and cleaning.....	21
8 Troubleshooting Guide.....	21
9 Technical data.....	22
10 Certificate.....	25

1 Introduction

Thank you for purchasing our **Refrigerated Bath Circulator**. You did a good choice buying a product from witeg Labortechnik.

Before use, please read this manual and use it properly.

To cover the user's health and avoid the product's destruction you have to follow the safety instructions of the second chapter.

Please keep this operating manual in a save place for future reference and hand it out if you consign the device to a third person.

Witeg reserves to modify the capsule or make amendments for the improvement of quality and efficiency without prior notice.

2 Safety Instructions

2.1 Description of symbols

**WARNING**

Information for avoiding injuries or fatal accidents.

**ATTENTION**

Safety instruction to avoid a damage of the product.

**ATTENTION**

Information about technical parameters.

**Environment**

Information about the operating environment.

**WARNING of explosions**

Safety instruction to avoid explosions.

**ATTENTION to high temperatures**

Instructions to avoid injuries caused by high temperatures.

**ATTENTION to low temperatures**

Instructions to avoid injuries caused by low temperatures.

2.2 Instructions

For a safety and faultless using of the device it is indispensable that you read the following instructions purposely and follow them by the use.

- Connect the device only to the suitable power supply. Please have a look to the identification plate.

- Use the device only in electric circuits with ground fault circuit interrupter.
- Check the power supply cord before you stick it into the electrical outlet.
- Do not use power bars to extend the power supply cord.
- Only use extension cables with sufficient cable cross section.
- Do not touch the device with wet hands during the usage because this might cause an electrical shock.
- If unexpected sound, smell or smog is generated by the device, pull out the main plug and contact the manufacturer or your supplier immediately.
- Disconnect the power supply plug before cleaning or maintenance.
- If you use radiation or contaminated samples the responsibility is only up to you.

- Warn all attendees if you leave the device unattended.
- Pay attention when you lift up the lid, hot condensate may drop down or hot steam could come out.
- Do not clean the scorching or frozen device.
- Pay attention and only work with gloves, otherwise you may freeze on the device's metal parts or get burns.
- Do not put thermal-sensitive items at the device.

- Use this device only if you read the whole safety instructions.
- This device should only be used by persons who are familiarized with safety instructions made for laboratories.
- Pay attention that no liquid reaches the controller or inward of the capsule.
- You have to keep the device clean if you want to get a longer life of it.
- Wear corresponding protective work clothing during the usage.
- Pull out the power cable if you didn't want to use the device for a longer time.

- It is important that the device stands up on a horizontal and firm subsurface to avoid vibrations.
- If you set the device in a site with vibrations, make sure that the device can not slip away.
- The device could get damaged if the vibrations are too strong.
- To avoid heat congestion or fires make sure that there is enough space for air circulations.
- Avoid direct solar radiation.
- Do not use the device near highly combustible materials.
- Do not set the device outdoors.
- It is superior to set the device in a good illuminated place to avoid accidents engendered by handling errors.
- To set the device in wet or dusty places could cause overheating, short circuits or fires.
- Make sure the device stands upright within the last 24 hours.
- Do not tilt the device.

- Do not use this device in area with potentially explosive atmosphere.

- Repairs and amendments should only be made by **witeg** authorized people. Contraventions exclude the warranty.
- Make sure, that the safety switch is adjusted up to 20-30°C over max. temperature. Otherwise the heating may not work.

3 Product description

3.1 Key features

The **Refrigerated Bath Circulator** is designed for chemical, biological, medical and pharmaceutical applications. The device was developed for simplest usage and a very long life.

Key features are as follows:

- Simplest operating by MAIN button and Jog-Dial knob.
- Long life by simply and solid construction.
- Highest temperature accuracy by digital fuzzy controller.
- Temperature and delay time showed on illuminated LCD.
- Storage functions for set time and temperature values.
- Turn-off and on delay from one minute up to over 99 hours in one minute steps adjustable.
- Signal sound after end of turn-off delay.
- Easy to justify.
- Key lock against inadvertent changing of set speed and timer values.
- Over temperature protection by safety switch.
- A lot of laboratory devices are mountable very quickly.
- Low filling level sensor.

3.2 General survey

3.3 Operating surface

3.4 Term description

SV: set value (upper line, in key lock menu items no SV)

PV: present value, measured by the device (lower line, in key lock menu item also upper line)

Main switch

- On / off switch.

Timer lamp

- Is blinking in menu items for timer modulation
- Shines continuous in the key lock menu item.

Heating lamp

- Shines if the heater is activated.

Refrigerator lamp

- Shines continuous in the key lock menu item when refrigerator runs.

Jog-Dial knob

- Change SV by turning.
- Push the knob for start and stop or probate set values.

Safety switch

- The safety switch turns the heating off if your set maximal temperature is overshoot.
- Set the maximal temperature rough 20 – 30 °C over the desired maximal temperature, otherwise the device maybe not reaches your set value.

MAIN button

- For changing between the different menu items.

SUB button

- For reaching the sub menu Items.

4 Packing contents

quantity	denomination
1	Refrigerated Bath Circulator
1	power supply cord
1	stainless steel lid
1	silicon tube
1	operating manual


5 Positioning the device and commissioning

1. Follow especially the safety instructions marked with following symbol:
2. Remove the package and perhaps soiling carefully.
3. Do not tilt the device. The cooling agent must not disperse in the cooling unit, otherwise the device could get damaged during the operation. Make sure that the device stands upright for the last 12 hours.
4. If there are any unexpected transport damages, please inform your forwarder and your costumer service.
5. Set the protection switch 30 °C over desired bath temperature otherwise the heating coil does not work.
6. Before you connect the device with the electrical outlet, please let some time laps until the device is acclimatised.
7. Make sure that you follow all safety instructions.
8. Connect the device with the power outlet.
9. Fill the operating liquid (also distilled water) in.
10. Switch the device on and start the process (look at job description).

6 Job description

6.1 Menu structure

6.2 Set turn-off delay

The turn-off delay timer serves as automatically stop function. It turns off the heating or refrigerator after the timeframe you set is elapsed.

- First menu item after switching the device on is the menu "set turn-off delay". The timer lamp is blinking.
- The delay time is adjustable from one minute till over 99 hours (00:01 - 99:59) in one minute steps.
- "00:00" stands for endurance run.
- Turn the Jog-Dial to choose your desired value and probate it by pushing the Jog-Dial. The delay time starts to run by probating the value.
- If you want to change the confirmed value, you have to push the Jog-Dial once again.

6.3 Set turn-on delay

The turn-on delay timer serves as automatically start function. It turns on the heating or refrigerator after the timeframe you set is elapsed.

- To reach this sub menu item you have to push the SUB button in the menu "set turn-off delay"
- The delay time is adjustable from one minute till over 99 hours (00:01 - 99:59) in one minute steps.
- Turn the Jog-Dial to choose your desired value and probate it by pushing the Jog-Dial. The delay time starts to run by probating the value.
- If you want to change the confirmed value, you have to push the Jog-Dial once again.
- Turn the Dial knob clockwise to increase the time and counter-clockwise to reduce it.
- The turn-off delay time will not elapse before the turn-on delay time is executed.

6.4 Set temperature

- To reach this menu item you have to push the MAIN button in the menu “set turn-off delay”
- The temperature is adjustable in 0,1°C-steps.
- Turn the Jog-Dial to choose your desired value and probate it by pushing the Jog-Dial. The delay time starts to run by probating the value.
- If you want to change the confirmed value, you have to push the Jog-Dial once again.
- Turn the Dial knob clockwise to increase the temperature and counter-clockwise to reduce it.
- The heating lamp shines, if the heater is activated.

The refrigerator only starts cooling if the liquid's temperature is fallen off to less than 39 °C. That is to protect the cooling unit from an overheat destruction.

6.5 Set Offset

If the temperature measured by you differs from the temperature the device is showing (PV) you have the possibility to compensate the divergence by the Offset. So you can balance differences of about $\pm 10^{\circ}\text{C}$.

The Offset is calculated as follows:

Temperature measured by you – showed Temperature (PV) = Offset

- To reach this menu you have to push the SUB button in the menu item “set temperature”.
- You can recognize this menu at the PV showed in the SV-line.
- If you turn the Jog-Dial knob.
- Set the calculated value by turning the Jog-Dial. The Offset will be added to the PV immediately.
- Save and confirm the Offset by pushing the MAIN button, you will also reach the previous menu.

6.6 Set Pump Power

- To reach this menu item you have to push the MAIN button in the menu “set temperature”.
- To set pump power turn the Jog Dial knob. The performance is adjustable in steps of 1-5. To confirm you have to press the Jog Dial knob.
- If you want to change the confirmed value, you have to push the Jog-Dial once again.

6.7 Key lock

- To avoid accidental changing of the set values, push the MAIN button in the menu item “set temperature”.
- In this menu item can you see the actually temperature in the SV-line and the remaining time of turn-off delay in the PV-line. The function of the Jog-Dial knob is deactivated.
- Press the MAIN button once more to leave the key lock. Than you are in the menu item “set turn-off delay”.

6.8 Storage function

The **Refrigerated Bath Circulator** includes a storage function for temperature and time values. That is practical if you want to work with the same parameters often. Then you didn't need to set them in every usage.

- For storing values you only need to hold the Jog-Dial knob pushed for two seconds. Successful storing is displayed by blinking of the SV for three times.
- Any stored parameters are erasable. Turn the value with the Jug-Dial till “0” and then hold the knob pressed for two seconds.
- If you want to overwrite a value it is not necessary to erase the old one. Set the desired

value and use the Jog-Dial knob for storing.

7 Maintenance and cleaning

- Please read the whole safety instructions before you start to maintain or clean the device.
- Pull the power plug out of the electrical outlet.
- Only use care cleaning agents and soft cleaning rags.
- Organic lotions, strong chemicals and rough cleaning rags can damage the device.
- If you not want to use the device for a long time, pull of the power plug and keep it dry in the package.
- Only use original spare parts.

8 Troubleshooting Guide

Situation	Confirmations and Solutions
Power Failure	✓ Check the supplied voltage is present. ✓ Check the power connection and the power cord for proper connection and integrity. ✓ Check if the MAIN Power Switch is "ON". ✓ Fuse is blown.
Unit Not Heating	✓ Check if the MAIN Power Switch is "ON". ✓ Check if the set temperature (SV) confirmed with the Jog-Dial knob (by turn the Jog-Dial knob cannot change the set temperature – see point 6.4) ✓ Check the Safety switch is set higher (min. 10°C-20°C) than the set temperature (SV).

Heating, but unit not reaching temperature	✓ Check the Safety switch is set higher (min. 10°C-20°C) than the set temperature (SV). ✓ Problem in the sensor or main controller.
Unit Not Refrigerating	✓ Check if the MAIN Power Switch is "ON". ✓ Check if the set temperature (SV) confirmed with the Jog-Dial knob (by turn the Jog-Dial knob cannot change the set temperature – see point 6.4) ✓ Check that the set temperature (SV) is lower than 35 °C.
Refrigerating but it doesn't fall to the SV temperature	✓ Problem with the the cooling circuit.
ERROR 1	✓ Temperature Sensor or main controller Error.
ERROR 3	✓ Low level error. Please fill water in.

If other problems arise, or one of the above not solves problems, please contact your official agent or the manufacturer.

Reparations and amendments should only made by **witeg** authorized people. Contraventions exclude the warranty.

9 Technical data

		WCR-P8	WCR-P12	WCR-P22	WCR-P30
Dimension (wxdxh)	Bath (mm) Effective space	175x153x150	230x181x150	290x220x200	330x290x200
	Bath (mm)	175x320x150	230x350x150	290x390x200	330x460x200
	Extern (mm)	330x474x752	400x523x807	483x563x850	494x690x906
Capacity		8 Lit.	12 Lit.	22 Lit.	30 Lit.
Temp.	Range	-25 °C~+150 °C			
	Accuracy	± 0,1 °C			
	Uniformity	±0.2°C at -10°C			
	Sensor	PT100			
Refrigerator		1/3 HP	1/3 HP	1/2 HP	3/4 HP
Cooling	at +20°C	300 W	380 W	680 W	850 W

Capacity	at 0°C	200 W	250 W	500 W	620 W
	at -20°C	110 W	130 W	300 W	400 W
Refrigerant		CFC-Free(R-404A) Refrigeration System			
Heating		1 kW	1.5 kW	2 kW	2.2 kW
Total power Consumption		1.6 kW	2 kW	2.7 kW	3 kW
Controller		Digital fuzzy controller with Jog-Dial knob (turn + push)			
Display		Digital LCD with Back-Light Function			
Timer & Alarm		99hr 59min (Delay & continuous run), Error status & Timer-end			
Material	Intern	Stainless steel			
	Extern	Powder coated steel			
	Insulation	40mm Urethane Foaming			
Pump		Max Flow Rate : 25Lit/min, Max Pressure : 7.5psi Pump Max. Height: 4 m			
Compressor Control		Intelligent Control of Refrigerator for preventing Compressor Overload			
Safety Circuit		Over-Temp./-Current protection, Sensor Error Detection, Leakage Breaker			
Others		Storage function (temp., timer); Key lock (Jog-Dial input disabled); Easy justification			
Packing Size (mm)		560x420x922	630x530x972	670x615x1022	740x625x1072
Gross Weight		45 kg	54 kg	63 kg	70 kg
Stromversorgung		AC 230 V, 1N~, 50/60 Hz			

		WCL-P8	WCL -P12	WCL -P22	WCL -P30
Dimension (wxdxh)	Bath (mm)	175x153x150	230x181x150	290x220x200	330x290x200
	Effective space				
	Bath (mm)	175x320x150	230x350x150	290x390x200	330x460x200
		330x474x752	400x523x807	483x563x850	494x690x906
Capacity		8 Lit.	12 Lit.	22 Lit.	30 Lit.
Temp.	Range	-35 °C~+150 °C			
	Accuracy	± 0,1 °C			
	Uniformity	±0.2°C at -20°C			
	Sensor	PT100			
Refrigerator		1/2 HP	1/2 HP	7/8 HP	1 HP
Cooling Capacity	500 W	530 W	860 W	1030 W	850 W
	340 W	370 W	620 W	760 W	620 W

	200 W	210 W	400 W	500 W	400 W
Refrigerant		CFC-Free(R-404A) Refrigeration System			
Heating		1 kW	1.5 kW	2 kW	2.2 kW
Total power Consumption		1.7 kW	2.2 kW	3 kW	3.3 kW
Controller		Digital fuzzy controller with Jog-Dial knob (turn + push)			
Display		Digital LCD with Back-Light Function			
Timer & Alarm		99hr 59min (Delay & continuous run), Error status & Timer-end			
Material	Intern	Stainless steel			
	Extern	Powder coated steel			
	Insulation	40mm Urethane Foaming			
Pump		Max Flow Rate : 25Lit/min, Max Pressure : 7.5psi Pump Max. Height: 4 m			
Compressor Control		Intelligent Control of Refrigerator for preventing Compressor Overload			
Safety Circuit		Over-Temp./-Current protection, Sensor Error Detection, Leakage Breaker			
Others		Storage function (temp., timer); Key lock (Jog-Dial input disabled); Easy justification			
Packing Size (mm)		560x420x922	630x530x972	670x615x1022	740x625x1072
Gross Weight		50 kg	55 kg	65 kg	78 kg
Stromversorgung		AC 230 V, 1N~, 50/60 Hz			

10 Certificate

Certificate

for

The quality and all features were checked by the manufacturer before the shipment.

We grant from date of purchase

two years guaranty.

This certificate excludes damages by natural disasters or incorrect usages by the costumer.

Please look on your account and complete following table:

Article	Refrigerated Bath Circulator
Typ	
Serialno.	
Date	

witeg Labortechnik GmbH

Am Bildacker 16
D-97877 Wertheim
TEL: +49-9342-9301-0
FAX: +49-9342-9301-77
Email: info@witeg.de
www.witeg.de