WVB

**Präzises digitales
Durchsichtthermostat
Bedienungsanleitung
Digital Precise
Viscosity Bath
Operating Manual**

gilt für / written for
WVB -30

witeg Labortechnik GmbH

(English manual at page 11)

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Contents

1 Introduction	12
2 Safety Instructions	12
2.1 Description of symbols	12
2.2 Instructions	12
3 Positioning the device and commissioning	14
4 Product description	14
4.1 Key features	14
4.2 Operating surface	15
4.3 Term description	15
5 Packing contents	16
6 Job description	16
6.1 Menu structure	16
6.2 Set turn-off delay	16
6.3 Set turn-on delay	17
6.4 Set temperature	17
6.5 Temperature Offset (Calibration)	17
6.6 Key lock	17
6.7 Storage function	18
7 Maintenance and cleaning	18
8 Troubleshooting Guide	18
9 Technical data	19
10 Certificate	20

1 Introduction

Thank you for purchasing our **WVB Digital Precise Viscosity Bath**. You did a good choice buying a **Wisd** product.

Before use, please read this manual and use it properly.

To cover the user's health and avoid the product's destruction you have to follow the safety instructions of the second chapter.

Please keep this operating manual in a save place for future reference and hand it out if you consign the device to a third person.

Witeg reserves to modify the capsule or make amendments for the improvement of quality and efficiency without prior notice.

2 Safety Instructions

2.1 Description of symbols

**WARNING**

Information for avoiding injuries or fatal accidents.

**ATTENTION**

Safety instruction to avoid a damage of the product.

**ATTENTION**

Information about technical parameters.

**Environment**

Information about the operating environment.

**WARNING of explosions**

Safety instruction to avoid explosions.

**ATTENTION to high temperatures**

Instructions to avoid injuries caused by high temperatures.

2.2 Instructions

For a safety and faultless using of the device it is indispensable that you read the following instructions purposely and follow them by the use.

- Connect the device only to the suitable power supply. Please have a look to the identification plate.

- Use the device only in electric circuits with ground fault circuit interrupter.
- Check the power supply cord before you stick it into the electrical outlet.
- Do not use power bars to extend the power supply cord.
- Only use extension cables with sufficient cable cross section.
- Do not touch the device with wet hands during the usage because this might cause an electrical shock.
- If unexpected sound, smell or smog is generated by the device, pull out the main plug and contact the manufacturer or your supplier immediately.
- Disconnect the power supply plug before cleaning or maintenance.
- If you use radiation or contaminated samples the responsibility is only up to you.

- Do not put heat-sensitive items on the device.
- Warn all attendees if you leave the device unattended.
- Pay attention when you lift up the lid, hot condensation water may drop down or hot steam could come out.
- Do not clean the scorching device.

- Reparations and amendments should only made by **witeg** authorized people. Contraventions exclude the warranty.
- Make sure, that the safety switch is adjusted up to 20-30°C over max. temperature. Otherwise the heating maybe doesn't work.

- Use this device only if you read the whole safety instructions.
- This device should only used by persons who are familiarized with safety instructions made for laboratories.
- Pay attention that no liquid reaches the controller or inward of the capsule.
- You have to keep the device clean if you want to get a longer life of it.
- Wear corresponding protective work clothing during the usage.
- Pull out the power cable if you didn't want to use the device for a longer time.

- Make sure that the device can not slip away.
- The device could get damaged if the vibrations are to strong.
- To avoid heat congestions or fires make sure that there is enough space for air circulations.
- Avoid direct solar radiation.
- Do not use the device near highly combustible materials.
- Do not set the device outdoors.
- It is superiorly to set the device in a good illuminated place to avoid accidents engendered by handling errors.

- To set the device in wet or dusty places could cause overheating, short circuits or fires.

- Do not use this device in areas with potentially explosive atmosphere.

3 Positioning the device and commissioning

1. Follow especially the safety instructions marked with following symbol:
2. Remove the package carefully.
3. If there are any unexpected transport damages, please inform your forwarder and your customer service.
4. Remove perhaps soiling and package material.
5. Set the protection switch 20 - 30 °C over desired water temperature otherwise the heating coil does not work.
6. Fill enough distilled water in.
7. Before you connect the device with the electrical outlet, please let some time laps until the device is acclimatised.
8. Make sure that you follow all safety instructions.
9. Connect the device with the power outlet.
10. Switch the device on and start the process (look at job description).

4 Product description


4.1 Key features

The **Wisd WVB Digital Precise Viscosity Bath** is designed for chemical, biological, medical and pharmaceutical applications. Used for Precise Measurement with Capillary Viscometers. Because of its long life and simplest usage it is more and more integrating in modern laboratories.

Key features are as follows:

- Simplest operating by patented combination of MODE button and Jog-Shuttle knob.
- Fast commissioning.
- Long life by simply and solid construction.
- Highest temperature accuracy by digital fuzzy controller.
- Temperature and delay time showed on illuminated LCD.
- Storage function for set time and temperature values.
- Turn-off and turn-on delay from one minute up to over 99 hours in one minute steps adjustable.
- Signal sound after end of turn-off delay.
- Easy to justify.
- Key lock against inadvertent changing of set speed and timer values.
- Over temperature protection by safety switch.
- Powerful Circulation Pump ensures Highest Uniformity of Temperature.
- Stainless Steel Lid Contains Holes for Insertion of Viscometer Holders.

4.2 Operating surface

4.3 Term description

SV: set value (upper line, in key lock menu item no SV)

PV: present value, measured by the device (lower line, in key lock menu item also upper line)

Main switch

- On / off switch.

Timer lamp

- Is blinking in menu items for timer modulation
- Shines continuous in the key lock menu item if turn-off delay time runs.

Heating lamp

- Shines if the heating is activated.

Jog-Shuttle knob

- Change SVs by turning.
- Push the knob for start and stop or probate set values.
- Hold it pushed for storing SVs.
- Used for key combination.

Safety switch

- The safety switch turns the heating off if your set maximal temperature is overshoot.
- Set the maximal temperature rough 20 – 30 °C over the desired maximal temperature, otherwise the device maybe not reaches your set value.

MODE button

- For changing between the different menu items.

Key combination

- Push and hold the MODE button first and in addition the Jog-Shuttle knob to reach the sub menu items, then release both buttons.


5 Packing contents

quantity	denomination
1	WVB Digital Precise Viscosity Bath
1	Flat Lid for Viscometer Holder
1	shelf
1	operating manual

6 Job description

6.1 Menu structure

The **Wisd.2 WVB Digital Precise Viscosity Bath** proffers the genial combination of MODE button and Jog-Shuttle knob. That allows a simple and comfortable usage. The menu structure is well-arranged and pictured as follows.

6.2 Set turn-off delay

The turn-off delay timer serves as automatically stop function. It turns off the heating when timeframe is elapsed.

- First menu item after switching the device on is the menu “set turn-off delay”. The timer lamp is blinking.
- The delay time is adjustable from one minute till over 99 hours (00:01 - 99:59) in one minute steps.
- “00:00” stands for endurance run.

- Turn the Jog-Shuttle to choose your desired value and probate it by pushing the Jog-Shuttle. The delay time starts to run by probating the value.
- If you want to change the confirmed value, you have to push the Jog-Shuttle once more. The PV starts blinking and is adjustable again.
- To stop the delay timer you can set "00:00" or push the MODE button when the value in PV line is blinking.

6.3 Set turn-on delay

The turn-on delay timer serves as automatically start function. It turns on the heating or turn-off delay timer when timeframe is elapsed. The turn-on delay time is only adjustable if the turn-off timer is not activated.

- To reach this sub menu item you have use the key combination (look at 3.4) in the menu "set turn-off delay". **LCD shows "d-tr"**.
- The delay time is adjustable from one minute till over 99 hours (00:01 - 99:59) in one minute steps.
- Turn the Jog-Shuttle to choose your desired value and probate it by pushing the Jog-Shuttle. The delay time starts to run by probating the value.
- If you want to change the confirmed value, you have to push the Jog-Shuttle once again.
- Turn the Shuttle knob clockwise to increase the time and counter-clockwise to reduce it.
- The turn-off delay time will not elapse before the turn-on delay time is executed.

6.4 Set temperature

- To reach this menu item you have to push the MODE button in the menu "set turn-off delay"
- The temperature is adjustable from 0 °C till 100 °C in 0,1 °C-steps.
- Turn the Jog-Shuttle to choose your desired value and probate it by pushing the Jog-Shuttle. The heating process starts shortly after probating the value when there is no turn-on delay time running.
- If you want to change the confirmed value, you have to push the Jog-Shuttle once more. The PV starts blinking and is adjustable again.
- To stop the heating process you have to push the MODE button when the value in PV line is blinking.
- The heating lamp shines, if the heater is activated.

6.5 Temperature Offset (Calibration)

If the temperature measured by you differs from the temperature the device is showing (PV) you have the possibility to compensate the divergence by the Offset. So you can balance differences of about ± 10 °C.

The Offset is calculated as follows:

Temperature measured by you – showed Temperature (PV) = Offset

- To reach this menu you have to use the key combination in the menu item "set temperature".
- Set the calculated value by turning the Jog-Shuttle. The Offset will be added to the PV immediately.
- Confirm the Offset by pushing the MODE button, you will also reach the previous menu.

6.6 Key lock

- To avoid accidental changing of the set values, push the MODE button in the menu item

”set temperature”.

- In this menu item can you see the actually temperature in the SV-line and the remaining time of turn-off delay time in the PV-line. The function of the Jog-Shuttle knob is deactivated.
- Press the MODE button once more to leave the key lock. Than you are in the menu item “set turn-off delay”.

6.7 Storage function

The **WVB Digital Precise Viscosity Bath** includes a storage function for temperature and timer values. That is practical if you want to work with the same parameters often. Then you didn't need to set them in every usage.

- For storing values you only need to hold the Jog-Shuttle knob pushed for two seconds. Successful storing is displayed by blinking of the SV for three times.
- Any stored parameters are erasable. Turn the value with the Jug-Shuttle till “0” and then hold the knob pressed for two seconds.
- If you want to overwrite a value it is not necessary to erase the old one. Set the desired value and use the Jog-Shuttle knob for over-storing.

7 Maintenance and cleaning

- Please read the whole safety instructions before you start to maintain or clean the device.
- Pull the power plug out of the electrical outlet.
- Only use care cleaning agents and soft cleaning rags.
- Organic lotions, strong chemicals and rough cleaning rags can damage the device.
- If you not want to use the device for a long time, pull of the power plug and keep it dry in the package.
- Only use original spare parts.

8 Troubleshooting Guide

Situation	Confirmations and Solutions
Power Failure	✓ Check the supplied voltage is present. ✓ Check the power connection and the power cord for proper connection and integrity. ✓ Check if the MAIN Power Switch is “ON”. ✓ Fuse is blown.
Unit Not Heating	✓ Check if the MAIN Power Switch is “ON”. ✓ Check if the set temperature (SV) confirmed with the Jog-Shuttle knob (by turn the Jog-Shuttle knob cannot change the set temperature – see point 6.4) ✓ Check the Safety switch is set higher (min. 10°C-20°C) than the set temperature (SV).
Heating, but unit not reaching temperature	✓ Check the Safety switch is set higher (min. 10°C-20°C) than the set temperature (SV). ✓ Problem in the sensor or main controller.
ERROR 1	✓ Temperature Sensor Error.

If other problems arise, or one of the above not solves problems, please contact your official agent or the manufacturer.

Reparations and amendments should only made by **witeg** authorized people. Contraventions exclude the warranty.

9 Technical data

		WVB-30
Dimension (w × d × h)	Effective Space	330 x 300 x 260 mm
	Bath	350 x 300 x 300 mm
	External	585 x 375 x 420 mm
Heater		2 kW
Capacity		30 Lit.
Temp.	Range	Ambient temperature + 5 °C up to 100°C
	Accuracy	±0.1 °C
	Sensor	PT 100
Controller		Digital Fuzzy controller with Jog-Shuttle knob (turn + push)
Display		Illuminated LCD
Timer		99 h 59 min for turn-off and turn-on delay
Safety Device		Over temp. and over current protector, sensor control
Pump		Magnetic Drive, Max. Flow 30 Lit./min
Lid		stainless steel
Viewing Window Size		230 x 180 mm
Others		- Storage function - Key lock (Jog-Shuttle input disabled) - Easy to justify - Signal sound (sensor error and timer-end)
Power Supply		AC 230 V, 1N~, 50/60 Hz

10 Certificate

Certificate

for

The quality and all features were checked by the manufacturer before the shipment.

We grant from date of purchase

two years guaranty.

This certificate excludes damages by natural disasters or incorrect usages by the costumer.

Please look on your account and complete following table:

Article	WVB Digital Precise Viscosity Bath
Typ	
Serialno.	
Date	

witeg Labortechnik GmbH

Am Bildacker 16
D-97877 Wertheim
TEL: +49-9342-9301-0
FAX: +49-9342-9301-77
Email: info@witeg.de
www.witeg.de