WITOPET

FRANCAIS / ENGLISH / DEUTSCH / ESPANOL

FRANCAIS

1 GENERALITES	4
2 REGLAGE DU VOLUME	6
3 UTILISATION	7
4 PRE-RINCAGE DU CONE	8
5 SOLUTIONS DENSES ET VISQUEUSES.....	8
6 CONES DE PRELEVEMENT	9
7 RECOMMANDATIONS	10
8 MAINTENANCE	11
9 EQUILMENT AUXILIARE	12
10 PIECES DETACHES	13
11 RECALIBRATION.....	13

1 GENERALITES

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1 INTRODUCTION

The pipette is a volumetric instrument designed to measure and transfer liquids precisely and safely. It can measure and transfer volumes from 0.2 μl to 5000 μl depending on the model.

The pipette has a digital counter which displays the volume. The volume is adjusted by turning the black knurled adjustment ring, Fig. 1B, and is continuously adjustable within the volume range for the pipette. The maximum volume of the pipette is shown on the push-button and corresponds to the model number. Fig. 1A.

The volume range from 0.2 μl to 5000 μl is covered by seven pipette models.

Model	Volume range (μl)		
PW2	0,2	-	2
PW10	0.5	-	10
PW20	2	-	20
PW100	10	-	100
PW200	20	-	200
PW1000	100	-	1000
PW5000	1000	-	5000

PW2,PW10	Measurement and transfer of micro-volumes, DNA sequencing and enzyme-assay applications.
PW20,PW100,PW200 et PW1000	Measurement and transfer of general aqueous solution, acids and bases
PW5000	Measurement and transfer of large volumes

The pipette uses disposable polypropylene tips. Fig. 1E. The disposable tips ensure maximum safety for the user and avoid cross contamination between samples. To protect the user from contaminated tips, the pipette is equipped with a built-in tip ejector. Fig. 1D.

When using very narrow test tubes, it maybe necessary to remove the tip ejector. It is simply removed by pulling it away from the pipette.

The pipette is a high quality instrument which offers excellent accuracy and precision. The values for accuracy and precision given in the following table were obtained using original tips. These figures are only guaranteed when original tips are used.

SPECIFICATIONS				
Model	Artivle no.	Volume	R[%]	VK[%]
PW2	5 402 002	Min 0,2	$\pm 12,0$	$\leq 6,0$
		1,0	$\pm 2,7$	$\leq 1,3$
		Max 2,0	$\pm 1,5$	$\leq 0,7$
PW10	5 402 010	Min 0,5	$\pm 4,0$	$\leq 4,0$
		5,0	$\pm 1,0$	$\leq 0,8$
		Max 10,0	$\pm 0,5$	$\leq 0,4$
PW20	5 402 020	Min 2	$\pm 4,0$	$\leq 3,0$
		10	$\pm 1,0$	$\leq 0,7$
		Max 20	$\pm 0,8$	$\leq 0,4$
PW100	5 402 100	Min 10	$\pm 1,6$	$\leq 0,8$
		50	$\pm 1,0$	$\leq 0,4$
		Max 100	$\pm 0,8$	$\leq 0,3$
PW200	5 402 200	Min 20	$\pm 1,2$	$\leq 0,6$
		100	$\pm 1,0$	$\leq 0,3$
		Max 200	$\pm 0,6$	$\leq 0,3$
PW1000	5 402 901	Min 100	$\pm 0,9$	$\leq 0,50$
		500	$\pm 0,7$	$\leq 0,50$
		Max 1000	$\pm 0,6$	$\leq 0,50$
PW5000	5 402 950	Min 1000	$\pm 0,6$	$\leq 0,3$
		2500	$\pm 0,6$	$\leq 0,3$
		Max 5000	$\pm 0,5$	$\leq 0,2$

These specifications are obtained in forward mode, using a gravimetric method with the temperature of the distilled water, tips and all other conditions stabilised between 19°C and 21°C. Number of measurements - minimum 15. The values given include all components of error due to both normal handwarming and the changing of the tip.

These specifications were obtained using tips of best quality.

Performance tests : The specifications for the volumetric performance of the pipette are derived from extensive gravimetric tests as prescribed by the Deutsches Institut für Normung e.V. To verify the performance of your pipette, request document DIN 12650.

2 SETTING THE VOLUME

The counter displays three figures to be read from top to bottom. Additional to the figures on the lower wheel are printed graduations of 2/10 which facilitate a volume setting in the range of 1/1000 of each pipette model.

Examples of the meaning of the black and red figures :

	PW10		PW20
Pipettes PW2,PW10,PW20: red figures = 1/10	0		1
	7		2
	5	red	5
	7.5 μ l		12.50 μ l

	PW100		PW200
Pipettes PW100 and PW200: black figures only = μ l	0		1
	7		2
	5		5
	75 μ l		125 μ l

	PW1000		PW5000
Pipettes PW1000 and PW5000: red figures at the top = ml	0	red	1
	7		2
	5		5
	0.75 ml		1.25 ml

The volume of the pipette is set by turning the black knurled adjustment ring, Fig. 1B. To obtain maximum accuracy and reproducibility always change the volume from the high to the lower value :

- For decreasing the volume turn the adjustment ring slowly to the required setting, making sure not to undergo the mark
- When increasing the volume turn the adjustment ring until the lower figure wheel comes 1/3 of a turn above the required setting. Now turn the ring slowly backward until the setting reaches the desired value. Be sure not to undergo this mark.

3 OPERATION

Place a tip on the shaft of the pipette. See Section 6 for the appropriate tip. Press the tip on firmly using a slight twisting motion to ensure a positive, airtight seal.

Important: Never aspirate liquids into the pipette without a tip attached.

Aspiration

Press the pushbutton to the first positive stop. Fig. 2A. Holding the pipette vertically, immerse the tip into the sample liquid. The depth to which the tip is immersed depends on the model.

PW2, PW1010	≤ 1 mm
PE20, PW100	2-3 mm
PW200, PW1000	2-4 mm
PW5000	3-6 mm

Release the pushbutton slowly and smoothly to aspirate the sample. Fig. 2B. Wait one second and then withdraw the tip from the liquid. Wipe any droplets away from the outside of the tip using a medical wipe.

Avoid touching the orifice of the tip.

Dispensing

Place the end of the tip against the inside wall of the vessel at an angle of 10 to 40 degrees. Press the pushbutton smoothly to the first stop, Fig. 2C. Wait one second. Press the pushbutton to the second stop to expel any remaining liquid. Fig. 2D. Keeping the pushbutton pressed to the end, remove the pipette by drawing the tip along the inside surface of the vessel. Release the pushbutton. Fig. 2E.

Eject the tip by pressing the tip ejector button. Fig. 2F. It is necessary to change the tip if a different liquid is being sampled or if the volumeter setting is changed.

4 PRE-RINSING

When pipetting liquids of higher viscosity or lower surface tension than water (e.g. sera or organic solvents), a film of liquid is formed on the inside wall of the pipette tip. This film can create an error. Since the film remains relatively constant in successive pipetting operations with the same tip, this error can be avoided by forming the film before transferring the first sample. This is done by aspirating a sample and dispensing it back into the same vessel. Since the film is already formed, all of the following samples will have better accuracy and repeatability.

This pre-rinsing operation should be repeated when the volume to be aspirated is changed or when a new tip is used.

5 DENSE AND VISCOUS LIQUIDS

The pipette specifications of accuracy and precision are based on pipetting distilled water. The handling of liquids with physical qualities of density, viscosity and surface tension differing extremely from water may need a gravimetrically checked compensation of the volume setting.

Normally the degree of error resulting from heavy or viscous liquids is negligible if the pipetting is done slowly and carefully. It is most important to give the liquids some time to follow the change of pressure by holding the pipette tip in its position for at least 2 sec. after the aspiration and the blow out stroke.

If in extreme cases this method of operation does not result in accurate values, a compensation could be achieved as follows :

Weigh the liquid pipetted when the pipette is set to the nominal volume.

Then calculate the set-off from the nominal volume :

$$\text{Corr. val.} = 2 \times \text{nom. val.} - \frac{\text{weight of liquid}}{\text{density of liquid}}$$

Check this operation once again and correct if necessary. Note the corrected value for further pipetting the same kind of liquid.

6 PIPETTE TIPS

High quality tips are made from high performance polypropylene and their quality guarantees the precision and accuracy associated with the pipette. Strict control is maintained throughout the manufacturing process to ensure the highest quality.

The accuracy and precision figures for the pipette are only guaranteed when quality tips are used. The use of inferior quality tips will seriously degrade the performance of the pipette.

Some tips are colour coded and correspond to the colour on the pushbutton.

White tips 10

These tips are used for volumes between 0.2 μl and 10 μl . They are used with the PW2, PW10 models.

Yellow tips 200 and white tips 200 with reference markers

These tips are used for volumes between 2 μl and 200 μl . They are used with the PW20, PW100, and the PW200 models.

Blue tips 1000

These tips are used for volumes between 100 μl and 1000 μl . They are used with the PW1000 model.

White tips 5000

These tips are used for volumes between 1000 μl and 5000 μl . They are used with the PW5000 model

7 RECOMMENDATIONS

The recommendations below will ensure maximum accuracy and precision from your pipette.

- Make sure to operate the pipette slowly and smoothly
- The depth of immersion in the sample liquid should be the minimum necessary and should remain constant during aspiration.
- The pipette should be held in a vertical position.
- Change the tip when volume setting is changed or when a different liquid is to be aspirated.
- Change the tip if a droplet remains on the end of the tip from the previous pipetting operation.
- Each new tip should be pre-rinsed with the liquid to be pipetted.
- Liquid should never enter the pipette shaft. To prevent this :
 - Press and release the pushbutton slowly and smoothly.
 - Never turn the pipette upside down.
 - Never lay the pipette on its side when there is liquid in the tip.
- Never force the volume setting beyond its recommended limits.
- When pipetting liquids with temperatures different from the ambient temperature, it is recommended to pre-rinse the tip several times before use.
- Do not pipette liquids with temperatures above 70°C.
- When pipetting acids or corrosive solutions which emit vapours, it is recommended to disassemble the shaft and to rinse the piston and seal with distilled water after finishing the pipetting operation.

8 TROUBLESHOOTING

If there is a leak or the pipette is not accurate, check for the following :

- The connecting nut is loose. Fig. 3B. **Tighten the connecting nut.**
- The shaft is cracked or scored. **Remove the tip ejector and inspect the shaft. For the PW2, PW10 and PW20 models, if the shaft is damaged, the piston may also be damaged. Replace the damaged parts (see Section 10). When reassembling the pipette, the nut should be hand tightened**
- Chemical damage to the piston and seals **Replace the piston and seals (See Section 10). Wash the inside of the shaft using distilled water.**
- Improper reassembly.
Disassemble the pipette and reassemble the pipette correctly. Fig. 3.

If there is liquid in the shaft, clean the pipette as follows :

- Remove the tip ejector. Unscrew the shaft nut and rinse the shaft, piston, seal and O-ring with distilled water. Dry these parts and reassemble the pipette.

Note : If necessary, the shaft, shaft nut and the tip ejector can be autoclaved for 20 minutes at 121°C and at 1 bar pressure. The seal and the O-ring should not be autoclaved.

The outside of the pipette can be cleaned using a cloth dampened with isopropanol.

If there is an air bubble when the sample is aspirated :

- Eject the sample into its original vessel.
- Check that the tip is properly immersed in the same liquid.
- Pipette the sample more slowly.
If the bubble appears for the second time, replace the tip.

If a problem continues after carrying out the above steps, send the pipette to your representative.

Before returning a pipette to a representative, please ensure that the instrument is completely free of any contamination (chemical, microbial or radioactive).

9 ACCESSORIES

Cat. no.	Description
5 408 010	White tips 10 µl, pack of 1000 pcs
5 409 005	White tips 10 µl, 96 pcs in a rack
5 408 200	Yellow tips 200 µl, pack of 1000 pcs
5 409 200	Yellow tips 200 µl, 96 pcs in a rack
5 410 000	Blue tips 1000µl, pack of 1000 pcs
5 409 210	Blue tips 1000 µl, 96 pcs in a rack
5 412 100	White tips 5000 µl, pack of 150 pcs
5 406 000	Stand for 5 pipettes

10 SPARE PARTS

All the spare parts indicated in Fig. 3, that is :

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can be ordered from our representatives (type of the pipette and name of the part for this pipette should be specified).

Warning : Replacing the piston will require recalibration of the pipette. We therefore recommend that this work be carried out by your representative.

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