

Bosshead


▪ suitable for all rods from Ø12mm - Ø23mm

Type	Material	Mounting	Order number
hook-type 90°	stainless steel	for rods up to Ø12mm	7 400 470
DIN12985, 90°, with screw	aluminium	for rods up to Ø16mm	7 400 500
DIN12985, 90°, with screw	cast iron	for rods up to Ø16mm	7 400 502
double-cross 90° & 180°	stainless steel	for rods up to Ø16mm	7 400 511
heavy-duty	cast aluminium	for rods up to Ø23mm	7 400 440



Clamp

Type	Material	Opening	Measurement	Order number
retort clamp 2 jaws PVC-coated	stainless steel	30mm		7 400 660
retort clamp 2 jaws PVC-coated	stainless steel	40mm		7 400 661
retort clamp 2 jaws PVC-coated	stainless steel	60mm		7 400 662
retort clamp 4 fingers silicone-coated	stainless steel	30-100 mm	190 x Ø12 mm	7 400 712
retort clamp 4 fingers silicone-coated	stainless steel	50-150 mm	190 x Ø12 mm	7 400 413
universal clamp 4 fingers silicone-coated cork-liner, hinge 180°	aluminium	0-80 mm	270 x Ø10 mm	7 400 720
condenser clamp 4 fingers PVC-coated	cast iron	50-100 mm	150 x Ø10 mm	7 400 730
flask clamp 2 jaws PVC-coated	chromed	12-45 mm	130 x Ø10 mm	7 400 830
burette clamp	chromed	4-20 mm	130 x Ø10mm	7 400 800
thermometer clamp with hook and bosshead	chromed	6-16 mm		7 400 810



Ring

Type	Material	Ring Ø	Length	Order number
closed	stainless steel	50 mm	160 mm	7 400 900
		70 mm	200 mm	7 400 901
		100 mm	220 mm	7 400 902
		120 mm	220 mm	7 400 903
		140 mm	220 mm	7 400 904
closed with bosshead	stainless steel	50 mm	70 mm	7 400 925
		70 mm	70 mm	7 400 926
		100 mm	70 mm	7 400 927
		120 mm	70 mm	7 400 928
		140 mm	70 mm	7 400 929

