

# Miscellaneous

## Plate stand with stainless steel rod

- plate with silicone rubber pad (ST110, ST120) or heavy duty safety stand (ST200, KA1177)


Type	Rod (included)	Plate	Article	Order number
Silicone rubber pad	Ø16xh800 mm with thread	20 x 35 cm	ST110	DH.WOS501510
Silicone rubber pad	Ø23xh1000 mm with thread	30 x 40 cm	ST120	DH.WOS501515
Heavy duty safety stand	Ø23xh1200 mm	39 x up to 75 cm adjustable	ST200	DH.WOS501520
Heavy duty safety stand	Ø23xh800 mm	34 x 29 cm	KA1177	DH.WOS1177



## Standard plate stand, thread M10

- with 1 adjustable foot, green painted steel, delivered without rod

Plate	Order number
210 x 130 mm	7 400 012
250 x 160 mm	7 400 014
315 x 200 mm	7 400 017



## Rods for plate stands

- with thread or without thread
- made of stainless steel

Type	Rod	Article	Order number
without thread	Ø12 x h500 mm		7 400 180
without thread	Ø12 x h600 mm		7 400 181
without thread	Ø12 x h750 mm		7 400 182
without thread	Ø12 x h1000 mm		7 400 183
without thread	Ø12 x h1250 mm		7 400 184
without thread	Ø12 x h1500 mm		7 400 185
without thread	Ø16 x h800 mm	RD501	DH.WOS502001
without thread	Ø16 x h1000 mm	RD502	DH.WOS502002
without thread	Ø16 x h1200 mm	RD503	DH.WOS502003
without thread	Ø23 x h800 mm	RD531	DH.WOS502031
without thread	Ø23 x 1200 mm	RD533	DH.WOS502033
with thread M10	Ø12 x h500 mm		7 400 160
with thread M10	Ø12 x h600 mm		7 400 161
with thread M10	Ø12 x h750 mm		7 400 162
with thread M10	Ø12 x h1000 mm		7 400 163
with thread M10	Ø12 x h1250 mm		7 400 164
with thread M10	Ø12 x h1500 mm		7 400 165
with thread M8	Ø16 x h800 mm	RD511	DH.WOS502011
with thread M8	Ø16 x h1000 mm	RD512	DH.WOS502012
with thread M8	Ø16 x h1200 mm	RD513	DH.WOS502013
with thread M8	Ø23 x h800 mm	RD521	DH.WOS502021
with thread M8	Ø23 x h1000 mm	RD522	DH.WOS502022
with thread M8	Ø23 x h1200 mm	RD523	DH.WOS502023



## Stand rods for hotplates and magnetic stirrer with hotplate

For device(s)	Design	Measurements	Article	Order number
HP-30A, HP-30D MSH-30A, MSH-30D	without thread	Ø12.7xh450 mm	RD100	DH.WHP503011
HP-20A, HP-20D, HP-500D MSH-20A, MSH-20D MSH-50A, MSH-50D MSH-500D, SMHS	with thread	Ø12xh450 mm	RD200	DH.WHP503036

