

Consumables

Reagent bottles - round - stopper - made of soda glass

- conical shoulder
- narrow or wide neck
- with glass or PE stopper
- clear glass or amber stained

Capacity	Socket	Neck	Order number			
ml	ST		Glass stopper		PE stopper	
			clear glass	amber stained	clear glass	amber glass
50	14/15	narrow	5 868 050	5 869 050	5 870 050	5 871 050
100	14/23	narrow	5 868 100	5 869 100	5 870 100	5 871 100
250	19/26	narrow	5 868 250	5 869 250	5 870 250	5 871 250
500	24/29	narrow	5 868 500	5 869 500	5 870 500	5 871 500
1000	29/32	narrow	5 868 001	5 869 001	5 870 001	5 871 001
2000	29/32	narrow	5 868 002	5 869 002	5 870 002	5 871 002
50	24/20	wide	5 873 050	5 874 050	5 875 050	5 876 050
100	29/22	wide	5 873 100	5 874 100	5 875 100	5 876 100
250	34/24	wide	5 873 250	5 874 250	5 875 250	5 876 250
500	45/27	wide	5 873 500	5 874 500	5 875 500	5 876 500
1000	60/64	wide	5 873 001	5 874 001	5 875 001	5 876 001



Bottles - square - screw thread - narrow neck - made of soda glass

- with or without PE-coating **NEW**
- clear glass or amber stained
- suitable tamper-proof caps: 5 897 032, 5 897 045

Capacity	Thread	Diameter	Height	Order number		
ml	GL	mm	mm	clear glass	clear glass PE-coating*	amber glass*
100	32	49	119	5 896 100	5 897 400	5 897 100
250	32	64	155	5 896 250	5 897 450	5 897 250
500	32	77	186	5 896 500	5 897 600	5 897 500
1000	45	97	223	5 896 001	5 897 701	5 897 001

* with pouring ring

Bottles - square - screw thread - wide neck - made of soda glass

- with or without PE-coating
- clear glass or amber stained
- suitable tamper-proof caps: 5 899 032, 5 899 045, 5 899 054, 5 899 060

Capacity	Thread	Diameter	Height	Order number		
ml	GL	mm	mm	clear glass	clear glass PE-coating	amber glass
100	32	49	111	5 898 100	5 899 110	5 899 100
250	45	64	146	5 898 250	5 899 260	5 899 250
500	54	76	173	5 898 500	5 899 510	5 899 500
1000	60	98	213	5 898 001	5 899 011	5 899 001



Tamper-evident caps for square bottles

- made of PP
- blue

Thread	Suitable for	Order number	
GL		narrow neck	wide neck
32	narrow neck: 100, 250, 500ml wide neck: 100ml	5 897 032	5 899 032
45	narrow neck: 1000ml wide neck: 250ml	5 897 045	5 899 045
54	wide neck: 500ml	-	5 899 054
60	wide neck: 1000ml	-	5 899 060

