

Consumables

Centrifuge tubes - with screw cap - heavy wall

- made of borosilicate glass 3.3

Capacity	Length	Outer Ø	Bottom	Order number
ml	mm	mm		
15	98	16	round	4 945 201
25	98	24		4 945 202
50	98	34		4 945 203
80	98	44		4 945 204
15	98	16	conical	4 945 301
25	98	24		4 945 302
50	98	34		4 945 303
80	98	44		4 945 304



Centrifuge tubes graduated - with screw cap - heavy wall

- made of borosilicate glass 3.3

Capacity	Length	Outer Ø	Bottom	Order number
ml	mm	mm		
15	98	16	round	4 945 401
25	98	24		4 945 402
50	98	34		4 945 403
80	98	44		4 945 404
15	98	16	conical	4 945 501
25	98	24		4 945 502
50	98	34		4 945 503
80	98	44		4 945 504



Centrifuge tubes - BIOGEN

- made of borosilicate glass 3.3, heavy wall
- with screw cap

Volume	Length	Outer Ø	Division	Thread	Bottom	Order number
ml	mm	mm	ml	GL		
12	137	17	1	18	conical	4 946 001
15	140	17	1	18	conical	4 946 002
40	145	28	0-10 : 0.5 10-40 : 1.0	25	conical	4 946 003
50	140	28	0-10 : 0.5 10-50 : 1.0	25	shallow cone	4 946 101
10	97	20	-	18	round	4 947 001
35	100	28	-	18	round	4 947 002
50	140	28	-	25	round	4 947 003
140	152	60	0-3 : 0.5 3-10 : 1.0 10-20 : 2.0 20-50 : 3.0 50-125 : 5.0	-	conical	4 947 101
140	152	60	-	32	conical	4 947 201
590	175	95	-	32	conical	4 947 202
760	175	95	-	45	shallow cone	4 947 301
140	150	60	0-3 : 0.5 3-10 : 1.0 10-20 : 2.0 20-50 : 3.0 50-125 : 5.0	32	conical	4 947 401
200	173	60	0-3 : 0.5 3-10 : 1.0 10-20 : 2.0 20-50 : 3.0 50-125 : 5.0	32	conical	4 947 402

