

## Test tubes

pack = 1000 pcs.

- made of clear PS or translucent PP

Volume ml	Ø x H mm	Material	Bottom	Order number pack
4	11 x 70	PS	round	5 449 001
4	12 x 55	PS	round	5 449 002
5	12 x 75	PS	round	5 449 003
5	12 x 75	PP	round	5 449 004
4,5	12 x 75	PS	conical	5 449 005
12	16 x 100	PS	round	5 449 006
12	16 x 100	PP	round	5 449 007
12	17 x 103	PS	conical	5 449 009
12	18 x 115	PS	conical	5 449 010
12	18 x 115	PS	round	5 449 011




## Centrifuge tubes - conical bottom

- with and without rim
- with and without graduation

pack = 100 pcs.



Capacity ml	Dimensions mm	Wall thickness mm	Made of	Rim	Graduation	Order number pack
~15	112 x 17	0,8 - 1,0	AR-glass	-	-	4 860 001
~15	112 x 17	0,8 - 1,0	AR-glass	yes	-	4 860 004
~15	112 x 17	0,8 - 1,0	Borosilicate	yes	-	4 860 005
~10	98 x 17	0,8 - 1,0	AR-glass	-	-	4 865 001
~10	98 x 17	0,8 - 1,0	AR-glass	yes	-	4 865 002
~15	112 x 17	0,8 - 1,0	AR-glass	-	10ml:0.1	4 870 001
~15	112 x 17	0,8 - 1,0	AR-glass	yes	mark at 15ml	4 870 004
~15	112 x 17	0,8 - 1,0	AR-glass	-	-	4 875 001
~15	112 x 17	0,8 - 1,0	AR-glass	yes	15ml:0.1	4 875 004
~15	112 x 17	0,8 - 1,0	Borosilicate	yes	-	4 875 005



## Centrifuge tubes - conical bottom 60° - heavy wall

- made of borosilicate glass 3.3

Capacity ml	Length mm	Outer Ø mm	Order number
15*	100	16	4 945 001
25	100	24	4 945 002
50	100	34	4 945 003
80	100	44	4 945 004
85	115	40	4 945 005

\* conical bottom 30°




## Centrifuge tubes graduated - heavy wall

- made of borosilicate glass 3.3



Capacity ml	Length mm	Outer Ø mm	Bottom	Graduation	Order number
15*	100	16	round	10ml:0.1	4 940 102
25*	100	24		20ml:0.1	4 940 103
50*	100	34		45ml:0.1	4 940 104
100*	100	44		85ml:0.1	4 940 105
250*	147	56		200ml:5	4 940 107
15	98	16	conical	10ml:0.1	4 945 101
25	98	24		20ml:0.1	4 945 102
50	98	34		45ml:0.1	4 945 103
100	98	44		85ml:0.1	4 945 104



\* according to DIN 58970 (Part 2)