

Test tubes

pack = 1000 pcs.

- made of clear PS or translucent PP

Volume	Ø x H	Material	Bottom	Order number
ml	mm			pack
4	11 x 70	PS	round	5 449 001
4	12 x 55	PS	round	5 449 002
5	12 x 75	PS	round	5 449 003
5	12 x 75	PP	round	5 449 004
4,5	12 x 75	PS	conical	5 449 005
12	16 x 100	PS	round	5 449 006
12	16 x 100	PP	round	5 449 007
12	17 x 103	PS	conical	5 449 009
12	18 x 115	PS	conical	5 449 010
12	18 x 115	PS	round	5 449 011




Centrifuge tubes - conical bottom

- with and without rim
- with and without graduation

pack = 100 pcs.



Capacity	Dimensions	Wall thickness	Made of	Rim	Graduation	Order number
ml	mm	mm				pack
~15	112 x 17	0,8 - 1,0	AR-glass	-	-	4 860 001
~15	112 x 17	0,8 - 1,0	AR-glass	yes		4 860 004
~15	112 x 17	0,8 - 1,0	Borosilicate	yes		4 860 005
~10	98 x 17	0,8 - 1,0	AR-glass	-		4 865 001
~10	98 x 17	0,8 - 1,0	AR-glass	yes	10ml:0.1 mark at 15ml	4 865 002
~15	112 x 17	0,8 - 1,0	AR-glass	-		4 870 001
~15	112 x 17	0,8 - 1,0	AR-glass	yes	15ml:0.1	4 870 004
~15	112 x 17	0,8 - 1,0	AR-glass	-		4 875 001
~15	112 x 17	0,8 - 1,0	AR-glass	yes		4 875 004
~15	112 x 17	0,8 - 1,0	Borosilicate	yes		4 875 005



Centrifuge tubes - conical bottom 60° - heavy wall

- made of borosilicate glass 3.3

Capacity	Length	Outer Ø	Order number
ml	mm	mm	
15*	100	16	4 945 001
25	100	24	4 945 002
50	100	34	4 945 003
80	100	44	4 945 004
85	115	40	4 945 005

* conical bottom 30°




Centrifuge tubes graduated - heavy wall

- made of borosilicate glass 3.3



Capacity	Length	Outer Ø	Bottom	Graduation	Order number
ml	mm	mm			
15*	100	16	round	10ml:0.1	4 940 102
25*	100	24		20ml:0.1	4 940 103
50*	100	34		45ml:0.1	4 940 104
100*	100	44		85ml:0.1	4 940 105
250*	147	56	conical	200ml:5	4 940 107
15	98	16		10ml:0.1	4 945 101
25	98	24		20ml:0.1	4 945 102
50	98	34		45ml:0.1	4 945 103
100	98	44		85ml:0.1	4 945 104



* according to DIN 58970 (Part 2)