


# Consumables

## Culture flasks (BIOGEN)

- Erlenmeyer shaped
- straight neck for metal caps, Ø 38mm
- made of borosilicate glass

Capacity	Flask Ø	Height	Order number
ml	mm	mm	
50	51	85	5 560 050
100	64	112	5 560 100
200	79	136	5 560 200
250	85	145	5 560 250
300	87	161	5 560 300
500	105	183	5 560 500
1000	131	232	5 560 001
2000	166	305	5 560 002
3000	187	310	5 560 003
5000	220	365	5 560 005



## Culture bottles (BIOGEN)

- straight neck for metal caps, Ø 38mm
- made from Duran® tubing



Capacity	Flask Ø	Height	Order number
ml	mm	mm	
100	50	125	5 562 100
300	71	170	5 562 300
500	83	208	5 562 500
1000	105	240	5 562 001



## Culture flasks (Fernbach)

- Fernbach type (regular or quad battled)
- straight neck for metal caps, Ø 38mm
- made of borosilicate glass

Capacity	Flask Ø	Height	Type	Order number
ml	mm	mm		
1800	200	155	regular	5 568 050
1800	200	155	quad battled	5 568 100



## Caps and stopper for culture flasks with neck Ø 38mm

- for neck Ø 38mm
- metal caps made of stainless steel or aluminium
- cap stopper made of plated anodised coloured aluminium with internal spring in nickel-chromed steel

Type	Material	Color	for neck Ø	Order number
			mm	
Metal cap	stainless steel	grey	38	5 570 001
	aluminium	grey	38	5 570 002
Cap stopper	plated anodised aluminium with internal spring	silver	37 - 39	5 570 011
		blue	37 - 39	5 570 012
		red	37 - 39	5 570 013

