


New Arrival

Digital Hotplate Stirrer MSH-R

- With round heating plate Ø 135 mm
- RS232 connector
- Digital temp. control up to 340°C
- 3-Year Warranty



Also available as Set

Warranty
3
Year

Order Now



E-Mail: info@witeg.de
Tel. : +49 9342 9301-0

witeg Labortechnik GmbH
Am Bildacker 16
97877 Wertheim Deutschland



www.witeg.de

witeg's Safety Stirring Beaker 2000 ml

Art. Nr.: 5 500 314 / 5 500 316

- With an absolutely flat bottom, witeg's Safety Stirring Beaker can **save up to 30% on energy** and be more effective/faster heating.
- The bottom rim provides a solid **no-slipping** Beaker with a **strong handle** - easy and safe to handle even when hot!
- The **plane bottom** ensures a steady and silent movement of the magnetic stirring bar
- Fits **all round magnetic stirrers with a Ø 135 mm - Ø 145 mm**

Exclusive

Efficiency and Safety, built in!

Magnetic stirring bars PTFE or Glass:

Art. Nr.: ab 7 21* *** / 7 211 0** G

- cylindrical, with ring, oval or triangular
- diameters between Ø3 mm - Ø20 mm
- lengths between 6 mm - 80 mm
- made from PTFE or **borosilicate glass (no PFAS)**

Also available:

- Magnetic stirring bar retrievers
- Magnetic stirring bar set with a retriever


Order Now

Precision heating, turn after turn.

Digital Hotplate Stirrer MSH-R

Art. Nr.: 6 242 300

- Different **heating modes - Fast, Standard, Stable**
- Digital temperature control up to **340°C**
- Digital speed control up to 1500rpm
- Max. stirring quantity of H₂O is 20L
- **“HOT” warning** flashes if plate temperature is **above 50°C - even after switching off.**
- High resolution LCD display shows actual temperature and speed
- Brushless DC motor is **maintenance free** and **explosion-proof**
- **ceramic coated plate** for high chemical resistance - Ø 135 mm
- RS232 connector - provides **PC control** and data transmission



Digital Hotplate Stirrer MSH-R Set

Art. Nr.: 6 242 301

The complete Set includes:

- MSH-R Digital Hotplate stirrer
- Stand
- External Temperature sensor with an **accuracy of $\pm 0.2^{\circ}\text{C}$**
- Magnetic stirring bar 30 x 6 mm PTFE



Order Now

Work smarter. Work Safer.

Product	Model	Order number
Digital Hotplate Stirrer	MSH-R	6 242 300
Digital Hotplate Stirrer Set	MSH-R Set	6 242 301
witeg's Safety Stirring Beaker 2000 ml for Plate Ø 135 mm	Ø 145 mm	5 500 314
witeg's Safety Stirring Beaker 2000 ml	Ø 165 mm	5 500 316
Magnetic Stirring Bars PTFE	30 x 6 mm - cilindrical	7 211 007
	30 x 6 mm - cilindrical with ring	7 212 006
	25 x 8 mm - triangular	7 213 002
	25 x 12 mm - oval	7 214 002
Magnetic Stirring Bars glass	25 x 8 mm - cilindrical	7 211 020 G
Magnetic Stirring Bar retriever	300 x 10 mm	7 218 003
Magnetic Stirring Bar set	with retriever	7 213 112

