

Spare parts & Accessories

Bottles made of amber soda-lime glass with thread and plastic coating

Volume ml	Thread GL	Width x Height mm	suitable for	Order number
100*	28	Ø51x100	MINISPENSOR	5 372 128
100	32	50x125	LABmax, TITREX	5 372 101
250	32	65x160	LABmax, TITREX	5 372 025
500	32	80x195	LABmax, TITREX	5 372 050
1000	45	95x230	LABmax, TITREX	5 372 100
2500	45	Ø140x230	LABmax, TITREX	5 372 250

*100ml & 2500ml bottle with round shape

Bottle made of borosilicate glass 3.3 with threads, plastic coating and screw cap

- straight neck GL 45 (closed screw cap)
- lateral neck GL 18 (screw cap with bore 11mm)

Volume ml	Thread GL	Ø x Height mm	suitable for	Order number
2000	45	Ø136x265	LABmax, TITREX	5 372 300

Standard joint adapters, made of borosilicate glass 3.3

Thread GL	Cone ST	suitable for	Order number
32	19/26	LABmax, TITREX, MINISPENSOR	0 550 029
32	24/29	LABmax, TITREX, MINISPENSOR	0 550 030
32	29/32	LABmax, TITREX, MINISPENSOR	0 550 031
45	29/32	LABmax, TITREX, MINISPENSOR	0 550 040
45	45/40	LABmax, TITREX, MINISPENSOR	0 550 041

Flexible discharge tube, coiled, with protective sleeve, made of FEP / PP

Length cm	Outer Ø mm	Inner Ø mm	suitable for	Order number
180	3.0	2.0	LABmax 2.5-10 ml	5 375 002
180	6	4	LABmax 25-100 ml, TITREX	5 375 012

Discharge tube with protective sleeve, made of FEP / PP

Length mm	Type	Nominal volume	suitable for	Order number
80	bent	2.5-10 ml	LABmax	5 375 003
120	bent	25-100 ml	LABmax	5 375 004
120	bent with fine tip	25-100 ml	TITREX	5 497 400
	bent with fine tip	100-2000 µl	MINISPENSOR	5 372 755

