

# Verbrauchsartikel

## Reagenzgläser FIOLAX

- hergestellt aus FIOLAX Borosilikatglas 3.3
- mit Bördelrand
- runder Boden



Pack = 100 St.

Länge	Außen Ø	Bestellnummer
mm	mm	Pack
70	8	4 902 001
100	10	4 902 002
100	12	4 902 003
100	16	4 902 004
130	14	4 902 005
160	16	4 902 006
180	18	4 902 007
180	20	4 902 008
200	25	4 902 009
200	30	4 902 010



## Einmal-Kulturröhrchen

- hergestellt aus AR-Glas
- glatter Rand
- runder Boden



Pack = 250 St.

Länge	Außen Ø	Bestellnummer
mm	mm	Pack
70	10	4 925 001
75	12	4 925 002
100	12	4 925 003
100	16	4 925 004
125	16	4 925 005
130	13	4 925 008
160	16	4 925 006
180	18	4 925 007



## Kulturröhrchen, glatter Rand

- hergestellt aus Borosilikatglas 3.3
- glatter Rand
- runder Boden



Pack = 100 St.

Länge	Außen Ø	Wandung	Bestellnummer
mm	mm	mm	Pack
75	12	0,8-1,0	4 903 001
100	12	0,8-1,0	4 903 002
130	16	1,0-1,2	4 903 003
160	16	1,0-1,2	4 903 004
180	18	1,0-1,2	4 903 005
150	20	1,0-1,2	4 903 006
150	25	1,0-1,2	4 903 008
200	25	1,0-1,2	4 903 009
200	30	1,0-1,4	4 903 010
100	10	0,8-1,0	4 903 013
130	14	0,8-1,0	4 903 014
180	18	1,0-1,2	4 903 020



## Reagenzgläser, starkwandig (~1mm)

- hergestellt aus AR-Glas
- glatter Rand
- runder Boden



Pack = 100 St.

Länge	Außen Ø	Bestellnummer
mm	mm	Pack
40	8	4 904 002
70	10	4 904 003
75	12	4 904 004
85	15	4 904 006
100	10	4 904 007
100	12	4 904 008
100	14	4 904 009
100	16	4 904 010
130	14	4 904 012
160	16	4 904 014
180	18	4 904 015



## Kulturröhrchen mit Bördelrand

- hergestellt aus Borosilikatglas 3.3
- DIN 12395
- runder Boden



Pack = 100 St.

Länge	Außen Ø	Wandung	Bestellnummer
mm	mm	mm	Pack
70*	8	0,8-1,0	4 903 100
75	12	0,8-1,0	4 903 101
100	12	0,8-1,0	4 903 102
130	16	1,0-1,2	4 903 103
160	16	1,0-1,2	4 903 104
180	18	1,0-1,2	4 903 105
150	20	1,0-1,2	4 903 106
150*	25	1,0-1,2	4 903 108
200	25	1,0-1,2	4 903 109
200	30	1,0-1,4	4 903 110
100	10	0,8-1,0	4 903 113
130	14	0,8-1,0	4 903 114
180	18	1,0-1,2	4 903 120



## Reagenzgläser mit Ausguss und Graduierung

- hergestellt aus Borosilikatglas 3.3
- mit Bördelrand und Ausguss
- runder Boden



Kapazität	Teilung	Länge	Außen Ø	Bestellnummer
ml	ml	mm	mm	
10	0,1	150	12	4 906 001
20	0,2	210	14	4 906 002
25	0,2	180	18	4 906 003
30	0,5	180	20	4 906 004



\* nicht gemäß DIN 12395