



TotalLab™ Quant


TotalLab Quant is a suite of core analysis tools for quantitative image analysis applications in the Life Sciences. Tools are included for analysis of both bands and spots from a variety of samples.

Following modules are included:


- 1D electrophoresis gel and Western blot analysis
- array / dot blot / slot blot analysis
- colony counting & basic 2D spot measurement
- general feature-based image analysis




1D DNA gel




1D Protein gel



Western blot



Protein array





Colony counting

The user then has full control of the visualisation tools and data display - outputting only those data fields that are of importance as well as the images of choice.

Features:

- compatible with 21CFR part 11 compliance module
- automatic detection of lanes and bands
- automatic background subtraction
- image manipulation tools
- lane templates
- molecular weight calibration
- quantity calibration & normalisation
- profile deconvolution
- reports
- array analysis module
- colony counting module
- toolbox for general analysis



Software

TotalLab™ Quant (included WGD-20S / WGD-30S)

Order number


DH.WGD00100

Phoretix 1D

Phoretix 1D analyses gel in quantitative and qualitative ways. It can import TIF images from different sources.

Features:

- automatic detection of lanes and bands
- automatic background subtraction
- image manipulation tools
- lane templates
- molecular weight calibration
- quantity calibration & normalisation
- profile deconvolution
- Rf calibration
- band picking
- band pattern matching - single gel
- dendrogram - single gel
- reports
- array analysis module
- colony counting module
- toolbox for general analysis



Phoretix 1D Pro

Phoretix 1D Pro includes all features of Phoretix 1D as well as a database for storing and controlling of 1D experiments, gels and lanes.

Additional features:

- band pattern matching - lanes across multiple gels
- band pattern queries
- dendrogram - lanes from multiple gels
- data archive and search facility
- classification and identification tools

Software

Phoretix 1D

Order number

DH.WGD00120

Phoretix 1D Pro

DH.WGD00121