

Consumables


Erlenmeyer flask with socket (conical flask)

- made of borosilicate glass 3.3
- with graduation
- clear glass, amber stained or economy
- according to DIN 12387

delivered without stopper!

Capacity ml	Socket ST	Outer Ø mm	Height mm	Order number		
				clear glass	amber stained	clear glass economy
5*	14/23	30	45	0 619 005	0 619 005 BR	0 619 005 E
10*		33	60	0 619 010	0 619 010 BR	0 619 010 E
25		42	70	0 619 025	0 619 025 BR	0 619 025 E
50		51	80	0 619 050	0 619 050 BR	0 619 050 E
100		64	100	0 619 100	0 619 100 BR	0 619 100 E
25	19/26	42	70	0 620 025	0 620 025 BR	0 620 025 E
50		51	80	0 620 050	0 620 050 BR	0 620 050 E
100		64	100	0 620 100	0 620 100 BR	0 620 100 E
200*		79	130	0 620 200	0 620 200 BR	-
250*		85	140	0 620 250	0 620 250 BR	0 620 250 E
50*	24/29	51	80	0 621 050	0 621 050 BR	0 621 050 E
100		64	100	0 621 100	0 621 100 BR	0 621 100 E
200*		79	130	0 621 200	0 621 200 BR	0 621 200 E
250*		85	140	0 621 250	0 621 250 BR	0 621 250 E
300*		87	150	0 621 300	0 621 300 BR	0 621 300 E
500*	29/32	105	170	0 621 500	0 621 500 BR	0 621 500 E
1000*		131	220	0 621 001	0 621 001 BR	0 621 001 E
50		51	80	0 622 050	0 622 050 BR	0 622 050 E
100		64	100	0 622 100	0 622 100 BR	0 622 100 E
200*		79	130	0 622 200	0 622 200 BR	0 622 200 E
250	45/40	85	140	0 622 250	0 622 250 BR	0 622 250 E
300*		87	150	0 622 300	0 622 300 BR	0 622 300 E
500		105	170	0 622 500	0 622 500 BR	0 622 500 E
1000		131	220	0 622 001	0 622 001 BR	0 622 001 E
2000*		166	280	0 622 002	-	-
200*	45/40	79	130	0 624 200	0 624 200 BR	-
250		85	140	0 624 250	0 624 250 BR	0 624 250 E
300*		87	150	0 624 300	0 624 300 BR	0 624 300 E
500		105	170	0 624 500	0 624 500 BR	0 624 500 E
1000		131	220	0 624 001	0 624 001 BR	0 624 001 E
2000*	166	280	0 624 002	0 624 002 BR	0 624 002 E	
3000	187	310	0 624 003	0 624 003 BR	-	
5000	220	365	0 624 005	0 624 005 BR	-	

* special size in addition to DIN 12387



Stopper - hexagonal plate

- according to DIN 12254 / DIN 12252



For socket	Order number		
ST	PE-stopper	Borosilicate 3.3 ST-hollow glass stopper	Borosilicate 3.3 ST-hollow glass stopper amber stained
14/23	0 306 014	0 302 014	0 303 014
19/26	0 306 019	0 302 019	0 303 019
24/29	0 306 024	0 302 024	0 303 024
29/32	0 306 029	0 302 029	0 303 029
45/40	0 306 045	0 302 045	0 303 045

