

Filter holder filtration set



Model	Rocker 300 & LF30	
Scope of supply	oil free vacuum pump Rocker 300 filtration set LF30 silicone tube (1m)	
Pump	110V (23 l/min, -680 mm Hg)	220V (20 l/min, -680 mm Hg)
Filtration set	300 ml PES funnel LF30 1200 ml PC waste bottle overflow protection spin-lock connections	
Ideal for	general filtration in laboratories air sampling suspended solid test	
Order number 110V	00829 167305-11	-
Order number 220V	-	00829 167305-22



Model	Rocker 300 & MF31	
Scope of supply	oil free vacuum pump Rocker 300 filtration set MF31 silicone tube (1m)	
Pump	110V (23 l/min, -680 mm Hg)	220V (20 l/min, -680 mm Hg)
Filtration set	300 ml PES funnel MF31 1000 ml PC receiving flask magnetic connections	
Ideal for	general filtration in laboratories air sampling suspended solid test	
Order number 110V	00829 167311-11	-
Order number 220V	-	00829 167311-22



Model	Rocker 300 C & VF12	
Scope of supply	oil free / chemical-resistant vacuum pump Rocker 300 C filtration set VF12 silicone tube (1m)	
Pump	110V (23 l/min, -680 mm Hg)	220V (20 l/min, -680 mm Hg)
Filtration set	300 ml glass funnel VF12 1000 ml glass receiving flask GL45aluminium clamp connection	
Ideal for	purification of organic, corrosive liquids such as solvent (mobile phase) for HPLC analysis	
Order number 110V	00829 189302-11	-
Order number 220V	-	00829 189302-22



Model	Chemker 300 & VF12	
Scope of supply	chemical-resistant vacuum pump Chemker 300 filtration set VF12 silicone tube (1m)	
Pump	110V (20 l/min, -670 mm Hg)	220V (18 l/min, -670 mm Hg)
Filtration set	300 ml glass funnel VF12 1000 ml glass receiving flask GL45aluminium clamp connection	
Ideal for	purification of organic, corrosive liquids such as solvent (mobile phase) for HPLC analysis	
Order number 110V	00829 169302-11	-
Order number 220V	-	00829 169302-22