

Filtration | stainless steel

Filter Ø47mm, 100ml / 500ml stainless steel filter holder

Ideal for:

- filtration requiring the sterilization such as microbiology test or use for purification of corrosive liquid, solvent, etc.




VF2 - silicone stopper
Ø47mm, 100ml



VF11 - silicone stopper
Ø47mm, 500ml


SPIN-LOCK CONNECTIONS



VF2

VF2 - 100ml funnel						
#	Item	Material	Size	Order number		
	47mm Stainless steel filtration set, complete				00829 167200-02	
	47mm Stainless steel filtration set, filter holder				00829 167103-10	
	Filter holder	Stainless steel	Ø47mm			
1	Funnel	Stainless steel	100 ml	00829 167103-61		
2	Support screen	Stainless steel		00829 167100-31		
3	Support base	Stainless steel		00829 167103-60		
4	Stopper	Silicone	No. 16	00829 167110-17		
5	Receiver flask	Borosilicate glass	1000 ml	00829 167300-08		
6	Fixing sucker	Silicone		00829 167200-45		

VF11 - 500ml funnel						
#	Item	Material	Size	Order number		
	47mm Stainless steel filtration set, complete				00829 167200-11	
	47mm Stainless steel filtration set, filter holder				00829 167103-50	
	Filter holder	Stainless steel	Ø47mm			
1	Funnel	Stainless steel	500 ml	00829 167103-65		
2	Support screen	Stainless steel		00829 167100-31		
3	Support base	Stainless steel		00829 167103-60		
4	Stopper	Silicone	No. 18	00829 167110-18		
5	Receiver flask	Borosilicate glass	2000 ml	00829 167200-13		



VF11