


# Overhead stirrer | Jar tester

Accessories: Stirring impellers and plate stands

## Plate stand with stainless steel rod

- suitable for all overhead stirrer and homogenizer
- plate with silicone rubber pad (ST110, ST120) or heavy duty safety stand (ST200, KA1177)


Type	Rod	Plate	Article	Order number
Silicone rubber pad	Ø16xh800 mm with thread	20 x 35 cm	ST110	DH.WOS501510
Silicone rubber pad	Ø23xh1000 mm with thread	30 x 40 cm	ST120	DH.WOS501515
Heavy duty safety stand	Ø23xh1200 mm	39 x up to 75 cm adjustable	ST200	DH.WOS501520
Heavy duty safety stand	Ø23xh800 mm	34 x 29 cm	KA1177	DH.WOS1177



## Universal jumbo clamp holder for plate stand

- suitable for all overhead stirrer and homogenizer plate stands

Type	Material	Mounting	Article	Order number
Jumbo clamp holder	cast aluminium	for Ø23 mm rod	CL200	DH.WOS502020



## Rods for plate stands

- with and without thread
- made of stainless steel

Type	Rod	Article	Order number
without thread	Ø16 x h800 mm	RD501	DH.WOS502001
without thread	Ø16 x h1000 mm	RD502	DH.WOS502002
without thread	Ø16 x h1200 mm	RD503	DH.WOS502003
without thread	Ø23 x h800 mm	RD531	DH.WOS502031
without thread	Ø23 x 1200 mm	RD533	DH.WOS502033
with thread	Ø16 x h800 mm	RD511	DH.WOS502011
with thread	Ø16 x h1000 mm	RD512	DH.WOS502012
with thread	Ø16 x h1200 mm	RD513	DH.WOS502013
with thread	Ø23 x h800 mm	RD521	DH.WOS502021
with thread	Ø23 x h1000 mm	RD522	DH.WOS502022
with thread	Ø23 x h1200 mm	RD523	DH.WOS502023



## Flexible couplings for safety stirring

Type	for shafts	Article	Order number
basic rubber / metal	Ø6-8 mm	W1-COR-6	1 445 000
basic rubber / metal	Ø8-10 mm	W1-COR-8	1 445 010
PTFE / Polyamid	Ø8mm		1 445 500
PTFE / Polyamid	Ø10mm		1 445 510
PTFE / Polyamid	Ø14mm		1 445 514
PTFE / Polyamid	Ø16mm		1 445 516
heavy-duty stainless steel	Ø6-8 mm	PKK-10	KG.40477
heavy-duty stainless steel	Ø12-14 mm	PKK-16	KG.40478

