


WITOPETTE plus

Stepper Pipette

Your advantage at a glance:

- ideal instrument for serial dosing in microliter range
- up to 48 pipetting steps in intervals of 1 second
- mechanical volume setting, 120 dosing volumes possible
- ergonomic design and low weighted for fatigue-proof dosing
- precise and accurate dosing of volumes from 1 to 5000 µl
- suitable for all established brands of dispensing tips from 0.05 - 50 ml
- robust & stable construction for a long life time
- quality certificate and serial number
- DE-M marked according to the German calibration law



WITOPETTE plus 1-5000 µl

Type	Order number
Stepper Pipette	5 385 205

WITIPS Tips for WITOPETTE plus

Volume (ml)		Packaging unit	Order number
0,05		100	5 385 222
0,5		100	5 385 223
0,5	sterile	100	5 385 252
1,25	sterile	100	5 385 224
1,25		100	5 385 253
2,5		100	5 385 225
2,5	sterile	100	5 385 254
5,0		100	5 385 226
5,0	sterile	100	5 385 255
12,5		100	5 385 227
12,5	sterile	100	5 385 256
25,0		25	5 385 228
25,0	sterile	25	5 385 257
50,0		25	5 385 229
50,0	sterile	25	5 385 258