


Rocker

5 - 50 rpm, tilt angle $\pm 6^\circ$ / $\pm 9^\circ$, platform included



RK-1D



RK-1D with optional stackable second platform

RK 5 - 50 rpm, settable angle

Ideal for:

- gel staining, destaining, blotting, cell culture, hybridization and washing

Features:

- 5 - 50 rpm
- 290x200 mm or 300x300 mm
- tilt angle can be set to $\pm 6^\circ$ and $\pm 9^\circ$
- digital LCD display
- smooth start and change of rocking speed while operation
- timer for 0-99 hours
- stainless steel platform with non-slip rubber pad included
- optional stackable second platform available
- 14W geared DC motor
- CE certified and unique serial number for tracing

Safety mechanism:

- timer 99h 59min with timer-end signal
- operation LEDs
- overload protection
- locking mode

Controller:

- programmable **Jog-Shuttle** controller with touch-button

| Model | RK-1D | RK-2D |
|--------------------------|---|-----------------------|
| Design | rocker | |
| Platform | 290x200 mm (included) | 300x300 mm (included) |
| Speed range | 5-50 rpm | |
| Angle | $\pm 6^\circ$ or $\pm 9^\circ$ ($\pm 9^\circ$ factory default) | |
| Programs | 6 programs, 10 steps | |
| Max. load | 8 kg | 10 kg |
| Timer & alarm | 99hr 59min, error status timer-end | |
| Display | digital LCD with backlight function | |
| Controller | digital feedback controller (Jog-Shuttle) | |
| Safety circuit | overload protection | |
| Dimension | 290x240x150 mm | 300x360x155 mm |
| Packaging size | 450x350xh250 mm, 7 kg | 450x350xh250 mm, 8 kg |
| Power supply | 1 Phase AC 120V, 60 Hz or AC 230V, 50/60 Hz | |
| Order number 230V | DH.WRK01010 | DH.WRK01020 |
| Order number 120V | DH.WRK02010 | DH.WRK02020 |

Accessories:



Second platform with 4 pillars and non-slip rubber pad

| RK5011 for RK-1D | RK5021 for RK-2D |
|------------------|------------------|
| 290x200 mm | 300x300 mm |
| DH.WRK05011 | DH.WRK05021 |