

# Rotator

5 - 60 rpm, 0 - 90° angle, multiple discs available

## RT-10 5 - 60 rpm, 0 - 90° angle

### Ideal for:

- general mixing, diffusing, dialyzing, hybridization and aerating culture
- keeping biological samples in suspension

### Features:

- 5 - 60 rpm
- 0-90° angle adjustable
- digital LCD display
- smooth start and change of stirring speed while operation
- timer for 0-99 hours
- 5 different disks available
- durable construction for long life time
- CE certified and unique serial number for tracing

### Safety mechanism:

- timer 99h 59min with timer-end signal
- operation LEDs
- overload protection
- locking mode

### Controller:

- programmable **Jog-Shuttle** controller with touch-button



RT-10 with disk WRT120

<b>Model</b>	<b>RT-10</b>
<b>Design</b>	turning motion
<b>Speed range</b>	5-60 rpm
<b>Angle range</b>	0-90°
<b>Programs</b>	6 programs, 10 steps
<b>Timer &amp; alarm</b>	99hr 59min, error status timer-end
<b>Display</b>	digital LCD with backlight function
<b>Controller</b>	digital feedback controller (Jog-Shuttle)
<b>Temperature resistance</b>	lower than 65°C (dry air)
<b>Safety circuit</b>	overload protection
<b>Dimension &amp; weight</b>	200x250xh260 mm, 6 kg
<b>Packaging size</b>	370x240xh430 mm, 8.5 kg
<b>Power supply</b>	1 Phase AC 120V, 60 Hz or AC 230V, 50/60 Hz
<b>Order number 230V</b>	<b>DH.WRT00010</b>
<b>Order number 120V</b>	<b>DH.WRT00020</b>

Please note: no disks included - please order separately!

### Scope of supply:


- Rotator RT-10
- 4x knobs, 4x posts WRT300



4x knobs, 4x posts for installing disks  
(included in scope of supply)

<b>Article</b>	<b>Order number</b>
WRT300	DH.WRT00300

### Accessories:




Rotator disc for 8x 50ml conical tubes

<b>Article</b>	<b>Order number</b>
WRT100	DH.WRT00100



Rotator disc for 16x 15ml conical tubes

<b>Article</b>	<b>Order number</b>
WRT110	DH.WRT00110



Rotator disc for 32x 1.5ml conical/centrifuge tubes

<b>Article</b>	<b>Order number</b>
WRT120	DH.WRT00120



Rotator disc for 60x 1.5ml conical/centrifuge tubes

<b>Article</b>	<b>Order number</b>
WRT130	DH.WRT00130



Rotator disc for 4x 250ml flasks

<b>Article</b>	<b>Order number</b>
WRT140	DH.WRT00140