


Consumables

Racks for test tubes, culture tubes and reaction vessels

- made of PP
- different colors available
- autoclavable up to 121°C
- detachable

Holes	Hole Ø mm	Dimensions mm	Order number			
			● white	● blue	● yellow	● red
100	micro tubes 1.5/2.0ml	262 x 108 x 45	7 073 102	7 073 102 BL	7 073 102 G	7 073 102 R
90	13	247 x 104 x 60	7 073 103	7 073 103 BL	7 073 103 G	7 073 103 R
60	16	247 x 104 x 70	7 073 104	7 073 104 BL	7 073 104 G	7 073 104 R
40	20	247 x 104 x 70	7 073 105	7 073 105 BL	7 073 105 G	7 073 105 R
40	25	297 x 124 x 85	7 073 106	7 073 106 BL	7 073 106 G	7 073 106 R
24	30	301 x 111 x 85	7 073 107	7 073 107 BL	7 073 107 G	7 073 107 R



Racks for test tubes and culture tubes

- made of polished wood



Holes	Depth of holes mm	Deck height mm	Decks	Dimensions mm	Order number
1 x 6	18	80	1	210 x 70 x 120	9 000 001
2 x 6				210 x 70 x 120	9 000 002
1 x 6			2	200 x 65 x 120	9 010 001
2 x 6				220 x 75 x 120	9 010 002
2 x 3	30	160	2	120 x 70 x 190	9 020 001
2 x 6				210 x 70 x 190	9 020 002
3 x 8			1	290 x 90 x 190	9 020 004
1 x 6				290 x 65 x 200	9 030 001
2 x 6	18	80	1	290 x 90 x 200	9 030 002
2 x 3				2	120 x 70 x 190
2 x 6			210 x 90 x 190		9 040 002
2 x 12			370 x 90 x 190	9 040 004	
4 x 12				370 x 190 x 150	9 040 006




Tray for reaction vessels

- made of PP
- 20 places for 1.5ml vessels
- numeration from 1-20

Order number

5 416 000



Racks for test tubes

- made of PP

Holes	Hole Ø mm	Dimensions mm	Order number
10	19	250 x 50 x 75	7 068 001
24	21	375 x 65 x 85	7 068 002
12	21	190 x 60 x 80	7 069 002

