

# Consumables

## Reagent bottles - round - stopper - made of soda glass

- conical shoulder
- narrow or wide neck
- with glass or PE stopper
- clear glass or amber stained

| Capacity<br>ml | Socket<br>ST | Neck   | Order number  |               |             |             |
|----------------|--------------|--------|---------------|---------------|-------------|-------------|
|                |              |        | Glass stopper |               | PE stopper  |             |
|                |              |        | clear glass   | amber stained | clear glass | amber glass |
| 50             | 14/15        | narrow | 5 868 050     | 5 869 050     | 5 870 050   | 5 871 050   |
| 100            | 14/23        | narrow | 5 868 100     | 5 869 100     | 5 870 100   | 5 871 100   |
| 250            | 19/26        | narrow | 5 868 250     | 5 869 250     | 5 870 250   | 5 871 250   |
| 500            | 24/29        | narrow | 5 868 500     | 5 869 500     | 5 870 500   | 5 871 500   |
| 1000           | 29/32        | narrow | 5 868 001     | 5 869 001     | 5 870 001   | 5 871 001   |
| 2000           | 29/32        | narrow | 5 868 002     | 5 869 002     | 5 870 002   | 5 871 002   |
| 50             | 24/20        | wide   | 5 873 050     | 5 874 050     | 5 875 050   | 5 876 050   |
| 100            | 29/22        | wide   | 5 873 100     | 5 874 100     | 5 875 100   | 5 876 100   |
| 250            | 34/24        | wide   | 5 873 250     | 5 874 250     | 5 875 250   | 5 876 250   |
| 500            | 45/27        | wide   | 5 873 500     | 5 874 500     | 5 875 500   | 5 876 500   |
| 1000           | 60/64        | wide   | 5 873 001     | 5 874 001     | 5 875 001   | 5 876 001   |



## Bottles - square - screw thread - narrow neck - made of soda glass

- with or without PE-coating **NEW**
- clear glass or amber stained
- suitable tamper-proof caps: 5 897 032, 5 897 045

| Capacity<br>ml | Thread<br>GL | Diameter<br>mm | Height<br>mm | Order number |                            |              |
|----------------|--------------|----------------|--------------|--------------|----------------------------|--------------|
|                |              |                |              | clear glass  | clear glass<br>PE-coating* | amber glass* |
| 100            | 32           | 49             | 119          | 5 896 100    | 5 897 400                  | 5 897 100    |
| 250            | 32           | 64             | 155          | 5 896 250    | 5 897 450                  | 5 897 250    |
| 500            | 32           | 77             | 186          | 5 896 500    | 5 897 600                  | 5 897 500    |
| 1000           | 45           | 97             | 223          | 5 896 001    | 5 897 701                  | 5 897 001    |

\* with pouring ring

## Bottles - square - screw thread - wide neck - made of soda glass

- with or without PE-coating
- clear glass or amber stained
- suitable tamper-proof caps: 5 899 032, 5 899 045, 5 899 054, 5 899 060

| Capacity<br>ml | Thread<br>GL | Diameter<br>mm | Height<br>mm | Order number |                           |             |
|----------------|--------------|----------------|--------------|--------------|---------------------------|-------------|
|                |              |                |              | clear glass  | clear glass<br>PE-coating | amber glass |
| 100            | 32           | 49             | 111          | 5 898 100    | 5 899 110                 | 5 899 100   |
| 250            | 45           | 64             | 146          | 5 898 250    | 5 899 260                 | 5 899 250   |
| 500            | 54           | 76             | 173          | 5 898 500    | 5 899 510                 | 5 899 500   |
| 1000           | 60           | 98             | 213          | 5 898 001    | 5 899 011                 | 5 899 001   |



## Tamper-evident caps for square bottles

- made of PP
- blue

| Thread<br>GL | Suitable for                                     | Order number |           |
|--------------|--|--------------|-----------|
|              |  | narrow neck  | wide neck |
| 32           | narrow neck: 100, 250, 500ml<br>wide neck: 100ml | 5 897 032    | 5 899 032 |
| 45           | narrow neck: 1000ml<br>wide neck: 250ml          | 5 897 045    | 5 899 045 |
| 54           | wide neck: 500ml                                 | -            | 5 899 054 |
| 60           | wide neck: 1000ml                                | -            | 5 899 060 |

