

Compressor (oil-free)

ROCKER
make filtration easy

Rocker 320 / 420 100 psi, 14 l/min, 50dB



Rocker 320



Rocker 420

Ideal for:

- aeration
- pressure filtration

Features:

- no air pollution
- maintenance-free
- quiet operation and low vibration
- direct driven motor
- air outlet filter device for filtering particles, moisture and supply clean air
- CE certified

Safety mechanism:

- thermal protection device ensures shutting down the pump when overheated and resumes working when the temperature cooled down

| Model | Rocker 320 | | Rocker 420 | |
|-----------------------------------|---------------------------------------------|-----------------|-----------------------|-----------------|
| Max. power (W) | 90 | 85 | 147 | 125 |
| Max. current (A) | 0.9 | 0.4 | 1.5 | 0.6 |
| Max. pressure (w/o moisture trap) | 100 psi | | 100 psi | |
| Max. pressure (w. moisture trap) | 80 psi | | 80 psi | |
| Max. flow rate | 14 l/min | 12 l/min | 27 l/min | 25 l/min |
| Moisture trap | yes | | yes | |
| Motor ration | 1700 rpm | 1400 rpm | 1700 rpm | 1400 rpm |
| Horse power | 1/8 HP | | 1/6 HP | |
| Noise level | 50 dB | | 52 dB | |
| Port thread | 5/16 inch (8 mm) | | 5/16 inch (8 mm) | |
| Net weight | 4.1 kg | | 5.2 kg | |
| Dimension (LxWxH) | 27.2 x 13.6 x 16.5 cm | | 29.8 x 14.3 x 17.2 cm | |
| Power supply | 1 Phase AC 120V, 60 Hz or AC 230V, 50/60 Hz | | | |
| Order number 110V | 00829 167320-11 | - | 00829 167420-11 | - |
| Order number 230V | - | 00829 167320-22 | - | 00829 167420-22 |