

Compressor (oil-free)



Rocker 320 / 420 100 psi, 14 l/min, 50dB



Rocker 320



Rocker 420

Ideal for:

- aeration
- pressure filtration

Features:

- no air pollution
- maintenance-free
- quiet operation and low vibration
- direct driven motor
- air outlet filter device for filtering particles, moisture and supply clean air
- CE certified

Safety mechanism:

- thermal protection device ensures shutting down the pump when overheated and resumes working when the temperature cooled down

Model	Rocker 320		Rocker 420	
Max. power (W)	90	85	147	125
Max. current (A)	0.9	0.4	1.5	0.6
Max. pressure (w/o moisture trap)	100 psi		100 psi	
Max. pressure (w. moisture trap)	80 psi		80 psi	
Max. flow rate	14 l/min	12 l/min	27 l/min	25 l/min
Moisture trap	yes		yes	
Motor ration	1700 rpm	1400 rpm	1700 rpm	1400 rpm
Horse power	1/8 HP		1/6 HP	
Noise level	50 dB		52 dB	
Port thread	5/16 inch (8 mm)		5/16 inch (8 mm)	
Net weight	4.1 kg		5.2 kg	
Dimension (LxWxH)	27.2 x 13.6 x 16.5 cm		29.8 x 14.3 x 17.2 cm	
Power supply	1 Phase AC 120V, 60 Hz or AC 230V, 50/60 Hz			
Order number 110V	00829 167320-11	-	00829 167420-11	-
Order number 230V	-	00829 167320-22	-	00829 167420-22