

# Vacuum pump (oil-free)

**ROCKER**  
make filtration easy

## Rocker 801 / 811

up to -735 mmHg, 80 l/min, 68dB



Rocker 801



Rocker 811

### Ideal for:

- filtration for microbiology test
- air sampling
- suspended solid test
- vacuum filtration

### Features:

- lubricant and oil free piston
- no air pollution
- maintenance-free
- Prote™ filter prevents liquid from being sucked into the pump and filters particles as well as moisture
- CE certified

### Safety mechanism:

- thermal protection device ensures shutting down the pump when overheated and resumes working when the temperature cooled down

Model	Rocker 801		Rocker 811	
Max. power (W)	220	200	220	200
Max. current (A)	2.0	0.9	2.0	0.9
Max. vacuum	-670 mmHg		-735 mmHg	
Max. flow rate	80 l/min	68 l/min	48 l/min	40 l/min
Moisture trap / vacuum regulator	Yes		Yes	
Overflow protection	No		No	
Motor ration	1750 rpm	1450 rpm	1750 rpm	1450 rpm
Horse power	1/3 HP		1/3 HP	
Noise level	68 dB		68 dB	
Port thread	3/8 inch (10 mm)		3/8 inch (10 mm)	
Net weight	7.9 kg		7.9 kg	
Dimension (LxWxH)	34.5 x 17.5 x 20 cm		34.5 x 17.5 x 20 cm	
Power supply	1 Phase AC 120V, 60 Hz or AC 230V, 50/60 Hz			
Order number 110V	00829 167801-11	-	00829 167811-11	-
Order number 230V	-	00829 167801-22	-	00829 167811-22