

Verbrauchsartikel

Schmelztiegel - Porzellan

- hergestellt aus Porzellan
- glasiert
- Tiegel und Deckel (muss einzeln bestellt werden)

Niedrige Form

Ø Oben	Höhe	Kapazität	Bestellnummer	
mm	mm	ml	Tiegel	Deckel
30	19	5	8 518 002	8 518 102
35	22	10	8 518 003	8 518 103
40	25	17	8 518 004	8 518 104
45	28	21	8 518 005	8 518 105
50	32	34	8 518 006	8 518 106
70	44	91	8 518 007	8 518 107
83	55	180	8 518 009	8 518 109



Mittlere Form

Ø Oben	Höhe	Kapazität	Bestellnummer	
mm	mm	ml	Tiegel	Deckel
25	20	5	8 519 001	8 519 101
30	25	10	8 519 002	8 519 102
35	28	12	8 519 003	8 519 103
40	32	20	8 519 004	8 519 104
45	36	30	8 519 005	8 519 105
50	40	45	8 519 006	8 519 106
52	43	60	8 519 007	8 519 107
60	48	80	8 519 008	8 519 108
63	55	110	8 519 009	8 519 109
70	56	120	8 519 010	8 519 110
82	65	230	8 519 011	8 519 111
89	75	300	8 519 012	8 519 112



Hohe Form

Ø Oben	Höhe	Kapazität	Bestellnummer	
mm	mm	ml	Tiegel	Deckel
26	34	11	8 520 003	8 520 103
30	38	15	8 520 004	8 520 104
35	44	26	8 520 005	8 520 105
40	50	35	8 520 006	8 520 106
45	56	50	8 520 007	8 520 107
50	62	72	8 520 008	8 520 108
60	75	130	8 520 009	8 520 109



Tiegel - Quarzglas

- hergestellt aus Quarzglas
- komplett (Tiegel und Deckel)



Niedrige Form

Ø Oben	Höhe	Kapazität	Bestellnummer
mm	mm	ml	
30	20	6	8 630 001
40	26	15	8 630 002
50	33	29	8 630 003
60	40	50	8 630 004
70	56	150	8 630 005

Mittlere Form

Ø Oben	Höhe	Kapazität	Bestellnummer
mm	mm	ml	
35	28	15	8 631 001
40	32	20	8 631 002
45	36	33	8 631 003
40	50	40	8 631 004
50	63	75	8 631 005

Hohe Form

Ø Oben	Höhe	Kapazität	Bestellnummer
mm	mm	ml	
40	40	30	8 633 001
45	45	40	8 633 002
50	50	50	8 633 003
60	60	90	8 633 004
65	65	100	8 633 005

