

## Reagent bottles - round - screw cap - made of soda glass

- narrow or wide neck
- with screw cap (PP)
- clear glass or amber glass

Capacity ml	Order number			
	narrow neck		wide neck	
	clear glass	amber glass	clear glass	amber glass
25	-	-	5 894 025	5 895 025
30	5 888 030	5 889 030	5 894 030	5 895 030
50	5 888 050	5 889 050	5 894 050	5 895 050
75	-	-	-	5 895 075
100	5 888 100	5 889 100	5 894 100	5 895 100
125	-	-	5 894 125	5 895 125
150	-	5 889 150	5 894 150	5 895 150
200	-	5 889 200	5 894 200	5 895 200
250	5 888 250	5 889 250	5 894 250	5 895 250
300	-	-	5 894 300	5 895 300
500	5 888 500	5 889 500	5 894 500	5 895 500
1000	5 888 001	5 889 001	5 894 001	5 895 001
2000	-	-	-	5 895 002



## Laboratory bottles - square - screw thread - graduated

- with screw cap and pouring ring (PP) (only clear glass version)
- clear glass or amber stained
- made of Duran® tubing

Capacity ml	Thread GL	Diameter mm	Height mm	Order number	
				clear glass	amber stained*
100	32	50	109	5 523 311	-
100	45	50	109	-	5 523 331
250	45	64	143	5 523 312	5 523 332
500	45	78	181	5 523 313	5 523 333
1000	45	94	222	5 523 314	5 523 334

\* delivered without cap and pouring ring



## Screw caps for bottles with thread

- made of PBTP or PP
- top closed or with bore

Thread GL	Bore	Material	Color	Seal	Order number
25	-	PBTP	red	PTFE liner	0 175 025*
32	-	PBTP	red	PTFE liner	0 175 032*
45	-	PBTP	red	PTFE liner	0 175 045*
50	-	PP	neutral	PTFE liner	0 175 050**
80	-	PP	neutral	PTFE liner	0 175 080**
100	-	PP	neutral	PTFE liner	0 175 099**
25	-	PP	blue	lip seal	0 175 125**
32	-	PP	blue	lip seal	0 175 132**
45	-	PP	blue	lip seal	0 175 145**
45	-	PP	yellow	-	5 526 066**
45	-	PP	green	-	5 526 077**
45	-	PP	red	-	5 526 088**
45	-	PP	neutral	-	5 526 099**
GLS80	-	PP	blue	lip seal	5 527 680**
GLS100	-	PP	blue	lip seal	5 527 690**
25	Ø15mm	PBTP	red	-	0 176 025*
32	Ø20mm	PBTP	red	-	0 176 032*
45	Ø34mm	PBTP	red	-	0 176 045*
50	Ø38mm	PP	neutral	-	0 176 050**
80	Ø68mm	PP	neutral	-	0 176 080**
100	Ø85mm	PP	neutral	-	0 176 100**



## Pouring rings for bottles

- made of PP or ETFE
- different colors
- autoclavable



Thread GL	Color	Material	Order number
32	blue	PP	0 175 201**
45	blue	PP	0 175 202**
45	yellow	PP	5 526 166**
45	green	PP	5 526 177**
45	red	PP	5 526 188**
45	white	PP	5 526 199**
32	red	ETFE	0 175 301*
45	red	ETFE	0 175 302*

\*autoclavable up to 180°C

\*\* autoclavable up to 140°C