

# Shaker

Accessories: Platforms, spring rack and flask holder


## Platforms

- suitable for Shaker SHO and SHR
- available in different types

Type	For device...	Plate	Article	Order number
Universal platform *	SHO-1D, SHR-1D	230x230 mm	SP100	DH.WSO504100
Universal platform **	SHO-2D, SHR-2D	350x350 mm	SP200	DH.WSO504200
Universal platform	OB-2, RP-2	251x251 mm	-	DH.WSO604100
Universal platform	OB-2, RP-2	346x346 mm	-	DH.WSO604200
Rubber mat platform	SHO-1D, SHR-1D	230x230 mm	SP110	DH.WSO504110
Rubber mat platform	SHO-2D, SHR-2D	370x370 mm	SP210	DH.WSO504210
Rubber mat platform	OB-2, RP-2	295x295 mm	-	DH.WSO604110
Rubber mat platform	OB-2, RP-2	370x370 mm	-	DH.WSO604210

\* max. 9x FH100 / FHS100, max. 7x FH250 / FHS250, max. 4x FH500 / FHS500

\*\* max. 20x FH100 / FHS100, max. 16x FH250 / FHS250, max. 9x FH500 / FHS500



## Platforms for flasks and separatory funnels

- suitable for Shaker SHO and SHR
- available in different types

Type	For device...	Dimension	Article	Order number
Platform for flasks	SHO-1D, SHR-1D	350x230x153 mm	SP120	DH.WSO504120
Platform for flasks	SHO-2D, SHR-2D	350x350x153 mm	SP220	DH.WSO504220
Platform for separatory funnels	SHO-1D, SHR-1D	350x235x174 mm	SP130	DH.WSO504130
Platform for separatory funnels	SHO-2D, SHR-2D	351x355x174 mm	SP230	DH.WSO504230
Platform for round flasks	OB-2, RP-2	350 260x153 mm	-	DH.WSO604130
Platform for separatory funnels	OB-2, RP-2	350x260x153 mm	-	DH.WSO604230


suitable separatory funnels can be found on page 221



## Universal spring rack

- suitable for Shaker SHO and SHR

Type	For device...	Dimension	Article	Order number
Universal spring rack	SHO-1D, SHR-1D	270x270x100 mm	SP140	DH.WSO504140
Universal spring rack	SHO-2D, SHR-2D	370x370x100 mm	SP240	DH.WSO504240
Universal spring rack	OB-2, RP-2	295x295x100 mm	-	DH.WSO604140
Universal spring rack	OB-2, RP-2	370x370x100 mm	-	DH.WSO604240



## Flask holder ABS or STS

- suitable for Shaker SHO and SHR
- holder clamps variable for different flask sizes (ABS flask holder with variable size)
- mountable on universal platform

Type	For device...	For flasks...	Article	Order number
ABS flask holder	SHO, SHR	50-100 ml	FH100	DH.WIS504030
ABS flask holder	SHO, SHR	200-300 ml	FH250	DH.WIS504035
ABS flask holder	SHO, SHR	500-1000 ml	FH500	DH.WIS504040
STS wire spring flask holder	SHO, SHR	100 ml	FHS100	DH.WIS600100
STS wire spring flask holder	SHO, SHR	250 ml	FHS250	DH.WIS600250
STS wire spring flask holder	SHO, SHR	500 ml	FHS500	DH.WIS600500
STS wire spring flask holder	SHO, SHR	1000 ml	FHS1000	DH.WIS601000

