

## Staining jar complete

- with glass box, cover, slide frame and wire holder
- for 10 slides 76x26mm

Description	Order number	
	Soda glass	Borosilicate 3.3
Complete	4 640 000	4 640 000 B
Spare tray (slide frame)	4 640 001	4 640 001 B
Spare wire holder	4 640 002	4 640 002 B
Spare glass box and cover	4 640 003	4 640 003 B



## Staining jar acc. to Hellendahl

- with glass box and cover
- standard and with enlargement
- for 16 slides 76x26mm

Description	Order number	
	Soda glass	Borosilicate 3.3
Standard	4 631 001	4 631 001 B
With enlargement	4 630 000	4 630 000 B



## Staining jar acc. to Coplin

- with glass box and cover (glass or PP)
- for 10 slides 76x26mm

Description	Dimension	Order number	
		Soda glass	Borosilicate 3.3
Standard with glass cover	115 x 68 mm	4 650 000	4 650 000 B
Standard with PP screw cap	113 x 68 mm	4 651 001	4 651 001 B
Low form with PP screw cap	98 x 68 mm	4 652 002	4 652 002 B



## Staining jar acc. to Schiefferdecker

- with glass box and cover
- for 10 slides 76x26mm

Description	Order number	
	Soda glass	Borosilicate 3.3
Complete	4 660 000	4 660 000 B



Description	Order number
Screw cap PP 55mm	4 652 100

## Petri dish

- made of glass
- for highest requirements, free from bubbles and imperfections
- bottom and lid (inside / outside) absolutely plane

Height mm	Outer Ø mm	Order number		Order number	
		Pack	Soda glass	Pack	Borosilicate 3.3
12	40		5 830 001		-
12	50	18	5 830 012		-
15	60	18	5 830 002		-
15	80	18	5 830 003	24	5 830 003 B
15	100	18	5 830 004	36	5 830 004 B
20	100	18	5 830 005	36	5 830 005 B
20	120	18	5 830 006		-
25	150	18	5 830 007	36	5 830 007 B
30	180	9	5 830 008		-
30	200	10	5 830 009		-
50	200	10	5 830 010		-

## Petri dish bacteriological

- made of clear PS
- packed in PE bags
- machine sterile or gamma radiated

Outer Ø mm	Machine sterile	Gamma radiation	with vents	Quantity in bag	Order number
55	X			15	5 488 001
55	X		X	15	5 488 002
55		X		15	5 488 003
55		X	X	15	5 488 004
90	X			20	5 488 011
90	X		X	33	5 488 012
90		X		20	5 488 013
90		X	X	33	5 488 014
100	X			11	5 488 031
100		X		11	5 488 032

