Bottle-top dispenser

Choose the right dispenser for your application

		LAB max premium	LAB max eco	LAB max airless	LAB max airless HF	MINISPENSOR	TITREX
Suitability	Aqueous solutions	Х	Х	Х	Х	Х	Х
	Crystalline solutions	X	X			X	X
	Solvents	X *	X *	X	X	X *	X*
	HF				X		
Volume setting	Quick-lock knob	X	X			X	
	Adjustment screw			X	X		
Special features	Air-purging	X	X	Х	X		X
	180° drip-free system	Х	X	Х	X		X
	Zero loss of reagent	Х	Х				Х
	360° safety turn	Х	х	Х	х		Х
	Built-in recalibration	Х					Х
Autoclavable	at 121°C	Х	X	Х	Х	X	X**
Dispensing system	Direct displacement	Х	Х			Х	Х
	Contactless piston			х	Х		
Certifications	Individual certificate	Х	Х	х	Х	Х	Х
	DE-M marked	X	X	X	X	X	X

^{*}For detailed information regarding the chemical resistance of the dispenser, please see page 16.

^{**}Only valve block autoclavable.

