

Miscellaneous

Tripod, for gas burner, made of stainless steel or tempered cast iron

| Outer Ø (mm) | Height (mm) | Order number stainless steel | Order number tempered cast iron |
|--------------|-------------|------------------------------|---------------------------------|
| 80 | 150 | 7 401 300 | - |
| 100 | 180 | 7 401 301 | 7 401 311 |
| 120 | 210 | 7 401 302 | 7 401 312 |
| 140 | 210 | - | 7 401 313 |
| 140 | 220 | 7 401 304 | - |
| 160 | 230 | 7 401 306 | - |




Triangles, with quartz tubes, made of stainless steel wire or nickel-chrome wire

| Length (mm) | Order number stainless steel wire | Order number nickel-chrome wire |
|-------------|-----------------------------------|---------------------------------|
| 40 | 7 401 320 | 7 401 330 |
| 50 | 7 401 321 | 7 401 331 |
| 60 | 7 401 322 | 7 401 332 |
| 70 | 7 401 323 | 7 401 333 |




Wire gauzes, stainless steel, flat

| Dimension (mm) | Order number |
|----------------|--------------|
| 135 x 135 | 7 401 350 |
| 155 x 155 | 7 401 351 |
| 175 x 175 | 7 401 352 |
| 195 x 195 | 7 401 353 |




Wire gauzes, iron wire with ceramic centre

| Dimension (mm) | Order number |
|----------------|--------------|
| 100 x 100 | 7 401 360 |
| 120 x 120 | 7 401 361 |
| 140 x 140 | 7 401 362 |
| 160 x 160 | 7 401 363 |
| 180 x 180 | 7 401 364 |
| 200 x 200 | 7 401 365 |




Lab-jacks

| Material | Plate (mm) | Height min (mm) | Height max (mm) | Load max (kg) | Order number |
|------------------|------------|-----------------|-----------------|---------------|--------------|
| stainless steel | 100x100 | 55 | 120 | 10 | 7 402 000 |
| stainless steel | 160x160 | 60 | 275 | 15 | 7 402 001 |
| stainless steel | 200x200 | 60 | 275 | 20 | 7 402 002 |
| stainless steel | 240x240 | 60 | 275 | 30 | 7 402 003 |
| stainless steel | 300x300 | 130 | 470 | 25 | 7 402 004 |
| aluminium, green | 100x100 | 55 | 120 | 10 | 7 402 010 |
| aluminium, green | 160x160 | 60 | 275 | 15 | 7 402 011 |
| aluminium, green | 200x200 | 60 | 275 | 20 | 7 402 012 |
| aluminium, green | 240x240 | 60 | 275 | 30 | 7 402 013 |
| aluminium, green | 300x300 | 130 | 470 | 25 | 7 402 014 |



stainless steel



aluminium green

Stopwatches, counter and timer

| Model | | Order number |
|--------------------------------------|---|--------------|
| digital stopwatch | 59min 59sec 1/100 start/stop, addition, split, reset, ABS housing, lanyard | 8 484 150 |
| counter/timer with alarm | timer up to 99min 99sec count-down 99min 99sec with magnet and clip | 8 484 160 |
| counter / timer with alarm and clock | timer up to 99min 59sec 1/100 count-down 99min 59sec 1/100 with magnet and clip | 8 484 200 |
| counter / timer with alarm and clock | timer up to 23h 59min 59sec count-down 23h 59min 59sec with magnet and clip | 7 425 100 |
| timer (mechanical) | 0 - 60min | 8 480 112 |



8 484 150



8 484 200



7 425 100



8 484 160



8 480 112