

# Consumables

## Volumetric flask with socket, class A



- made of Duran® tubing
- tolerance according to DIN EN ISO 1042
- DE-M marked according to the German calibration law
- standard and special sizes available
- inclusive batch certificate

| Capacity | Tolerance | Socket | Order number                    |                                |                                              |
|----------|-----------|--------|---------------------------------|--------------------------------|----------------------------------------------|
| ml       | ml        | ST     | brown graduated with PE-stopper | blue graduated with PE-stopper | brown graduated with ST-hollow glass stopper |
| 5        | ±0,025    | 7/16   | 3 664 000                       | 3 664 000 BL                   | 3 668 000                                    |
| 5        | ±0,040    | 10/19  | 3 664 015                       | 3 664 015 BL                   | 3 668 015                                    |
| 10       | ±0,025    | 7/16   | 3 664 001                       | 3 664 001 BL                   | 3 668 001                                    |
| 10       | ±0,040    | 10/19  | 3 664 012                       | 3 664 012 BL                   | 3 668 012                                    |
| 10*      | ±0,060    | 14/23  | 3 664 00114                     | 3 664 00114 BL                 | 3 668 00114                                  |
| 10*      | ±0,060    | 12/21  | -                               | 3 664 01212 BL                 | -                                            |
| 15*      | ±0,040    | 10/19  | 3 664 016                       | 3 664 016 BL                   | 3 668 016                                    |
| 20*      | ±0,040    | 7/16   | 3 664 021                       | 3 664 021 BL                   | 3 668 021                                    |
| 20       | ±0,040    | 10/19  | 3 664 002                       | 3 664 002 BL                   | 3 668 002                                    |
| 20       | ±0,060    | 12/21  | 3 664 013                       | 3 664 013 BL                   | 3 668 013                                    |
| 20*      | ±0,060    | 14/23  | 3 664 00214                     | 3 664 00214 BL                 | 3 668 00214                                  |
| 25       | ±0,040    | 10/19  | 3 664 003                       | 3 664 003 BL                   | 3 668 003                                    |
| 25       | ±0,060    | 12/21  | 3 664 014                       | 3 664 014 BL                   | 3 668 014                                    |
| 25*      | ±0,060    | 14/23  | 3 664 00314                     | 3 664 00314 BL                 | 3 668 00314                                  |
| 50       | ±0,060    | 12/21  | 3 664 005                       | 3 664 005 BL                   | 3 668 005                                    |
| 50       | ±0,100    | 14/23  | 3 664 006                       | 3 664 006 BL                   | 3 668 006                                    |
| 75*      | ±0,060    | 12/21  | 3 664 0075                      | -                              | 3 668 0075                                   |
| 100      | ±0,100    | 12/21  | 3 664 010                       | 3 664 010 BL                   | 3 668 010                                    |
| 100      | ±0,100    | 14/23  | 3 664 011                       | 3 664 011 BL                   | 3 668 011                                    |
| 100*     | ±0,150    | 19/26  | 3 664 135                       | -                              | 3 668 135                                    |
| 150*     | ±0,100    | 14/23  | 3 664 150                       | -                              | 3 668 150                                    |
| 192,2*   | ±0,150    | 14/23  | 3 664 192                       | 3 664 192 BL                   | 3 668 192                                    |
| 200      | ±0,150    | 14/23  | 3 664 020                       | 3 664 020 BL                   | 3 668 020                                    |
| 250      | ±0,150    | 14/23  | 3 664 025                       | 3 664 025 BL                   | 3 668 025                                    |
| 300*     | ±0,150    | 14/23  | 3 664 030                       | -                              | 3 668 030                                    |
| 350*     | ±0,150    | 19/26  | 3 664 035                       | -                              | 3 668 035                                    |
| 368,3*   | ±0,150    | 19/26  | -                               | 3 664 368 BL                   | -                                            |
| 400*     | ±0,250    | 19/26  | 3 664 040                       | -                              | 3 668 040                                    |
| 500*     | ±0,250    | 14/23  | 3 664 049                       | -                              | 3 668 049                                    |
| 500      | ±0,250    | 19/26  | 3 664 050                       | 3 664 050 BL                   | 3 668 050                                    |
| 600*     | ±0,400    | 24/29  | 3 664 060                       | -                              | 3 668 060                                    |
| 800*     | ±0,400    | 24/29  | 3 664 080                       | -                              | 3 668 080                                    |
| 900*     | ±0,400    | 24/29  | 3 664 090                       | -                              | 3 668 090                                    |
| 1000     | ±0,400    | 24/29  | 3 664 100                       | 3 664 100 BL                   | 3 668 100                                    |
| 1000     | ±0,600    | 29/32  | 3 664 101                       | 3 664 101 BL                   | 3 668 101                                    |
| 2000     | ±0,600    | 29/32  | 3 664 200                       | 3 664 200 BL                   | 3 668 200                                    |
| 3000*    | ±0,800    | 29/32  | 3 664 300                       | 3 664 300 BL                   | 3 668 300                                    |
| 5000*    | ±1,200    | 29/32  | 3 664 499                       | -                              | 3 668 499                                    |
| 5000     | ±1,200    | 34/35  | 3 664 500                       | 3 664 500 BL                   | 3 668 500                                    |

\* special size in addition to the ISO 1042 range, with conformity certificate



### NEW

All volumetric flasks with ground joint neck can be ordered with spout.

Just add „A“ at the end of the order number.

**Example:** „3 664 005 A“ or „3 664 100 BL A“

### NEW

All brown graduated volumetric flasks with ST-hollow glass

stoppers can be ordered with blue graduation.

Just add „BL“ at the end of the order number.

**Example:** „3 668 014 BL“ or „3 668 500 BL“

