

Consumables

Volumetric flask with socket, class A



- made of Duran® tubing
- tolerance according to DIN EN ISO 1042
- DE-M marked according to the German calibration law
- standard and special sizes available
- inclusive batch certificate

Capacity	Tolerance	Socket	Order number		
ml	ml	ST	brown graduated with PE-stopper	blue graduated with PE-stopper	brown graduated with ST-hollow glass stopper
5	±0,025	7/16	3 664 000	3 664 000 BL	3 668 000
5	±0,040	10/19	3 664 015	3 664 015 BL	3 668 015
10	±0,025	7/16	3 664 001	3 664 001 BL	3 668 001
10	±0,040	10/19	3 664 012	3 664 012 BL	3 668 012
10*	±0,060	14/23	3 664 00114	3 664 00114 BL	3 668 00114
10*	±0,060	12/21	-	3 664 01212 BL	-
15*	±0,040	10/19	3 664 016	3 664 016 BL	3 668 016
20*	±0,040	7/16	3 664 021	3 664 021 BL	3 668 021
20	±0,040	10/19	3 664 002	3 664 002 BL	3 668 002
20	±0,060	12/21	3 664 013	3 664 013 BL	3 668 013
20*	±0,060	14/23	3 664 00214	3 664 00214 BL	3 668 00214
25	±0,040	10/19	3 664 003	3 664 003 BL	3 668 003
25	±0,060	12/21	3 664 014	3 664 014 BL	3 668 014
25*	±0,060	14/23	3 664 00314	3 664 00314 BL	3 668 00314
50	±0,060	12/21	3 664 005	3 664 005 BL	3 668 005
50	±0,100	14/23	3 664 006	3 664 006 BL	3 668 006
75*	±0,060	12/21	3 664 0075	-	3 668 0075
100	±0,100	12/21	3 664 010	3 664 010 BL	3 668 010
100	±0,100	14/23	3 664 011	3 664 011 BL	3 668 011
100*	±0,150	19/26	3 664 135	-	3 668 135
150*	±0,100	14/23	3 664 150	-	3 668 150
192,2*	±0,150	14/23	3 664 192	3 664 192 BL	3 668 192
200	±0,150	14/23	3 664 020	3 664 020 BL	3 668 020
250	±0,150	14/23	3 664 025	3 664 025 BL	3 668 025
300*	±0,150	14/23	3 664 030	-	3 668 030
350*	±0,150	19/26	3 664 035	-	3 668 035
368,3*	±0,150	19/26	-	3 664 368 BL	-
400*	±0,250	19/26	3 664 040	-	3 668 040
500*	±0,250	14/23	3 664 049	-	3 668 049
500	±0,250	19/26	3 664 050	3 664 050 BL	3 668 050
600*	±0,400	24/29	3 664 060	-	3 668 060
800*	±0,400	24/29	3 664 080	-	3 668 080
900*	±0,400	24/29	3 664 090	-	3 668 090
1000	±0,400	24/29	3 664 100	3 664 100 BL	3 668 100
1000	±0,600	29/32	3 664 101	3 664 101 BL	3 668 101
2000	±0,600	29/32	3 664 200	3 664 200 BL	3 668 200
3000*	±0,800	29/32	3 664 300	3 664 300 BL	3 668 300
5000*	±1,200	29/32	3 664 499	-	3 668 499
5000	±1,200	34/35	3 664 500	3 664 500 BL	3 668 500

* special size in addition to the ISO 1042 range, with conformity certificate



NEW

All volumetric flasks with ground joint neck can be ordered with spout.

Just add „A“ at the end of the order number.

Example: „3 664 005 A“ or „3 664 100 BL A“

NEW

All brown graduated volumetric flasks with ST-hollow glass

stoppers can be ordered with blue graduation.

Just add „BL“ at the end of the order number.

Example: „3 668 014 BL“ or „3 668 500 BL“

