

Volumetric pipettes

Class AS


- color-code to avoid confusion
- made of soda-lime glass
- tolerance according to DIN EN ISO 648
- DE-M marked according to the German calibration law
- calibrated to deliver (EX), 5 sec waiting time
- inclusive batch certificate



Capacity ml	Tolerance ml	Color-Code	Order number	
			brown graduated	blue graduated
0.5	±0.005	black x2	3 535 005	3 535 005 BL
1	±0.008	blue	3 535 010	3 535 010 BL
2	±0.010	orange	3 535 020	3 535 020 BL
5	±0.015	white	3 535 050	3 535 050 BL
10	±0.020	red	3 535 100	3 535 100 BL
20	±0.030	yellow	3 535 200	3 535 200 BL
25	±0.030	blue	3 535 250	3 535 250 BL
50	±0.050	red	3 535 500	3 535 500 BL
100	±0.080	yellow	3 535 900	3 535 900 BL

Class AS - special sizes

- color-code to avoid confusion
- made of soda-lime glass
- tolerance on basis of DIN EN ISO 648
- DE-M marked according to the German calibration law
- calibrated to deliver (EX), 5 sec waiting time
- inclusive batch certificate



Capacity ml	Tolerance ml	Color-Code	Order number	
			brown graduated	blue graduated
3	±0.010	black	3 535 030	3 535 030 BL
3.5	±0.010	-	3 535 035	3 535 035 BL
4	±0.015	red x2	3 535 040	3 535 040 BL
4.5	±0.015	-	3 535 045	3 535 045 BL
5.5	±0.015	-	3 535 055	3 535 055 BL
5.85	±0.015	-	3 535 0585	3 535 0585 BL
6	±0.015	orange x2	3 535 060	3 535 060 BL
7	±0.015	green x2	3 535 070	3 535 070 BL
7.5	±0.015	-	3 535 075	3 535 075 BL
8	±0.020	blue	3 535 080	3 535 080 BL
9	±0.020	black	3 535 090	3 535 090 BL
11	±0.020	-	3 535 110	3 535 110 BL
12	±0.020	-	3 535 120	3 535 120 BL
12.5	±0.020	-	3 535 125	3 535 125 BL
13	±0.020	-	3 535 130	3 535 130 BL
14	±0.020	-	3 535 140	3 535 140 BL
15	±0.020	green	3 535 150	3 535 150 BL
16	±0.030	-	3 535 160	3 535 160 BL
17	±0.030	-	3 535 170	3 535 170 BL
18	±0.030	-	3 535 180	3 535 180 BL
19	±0.030	-	3 535 190	3 535 190 BL
21	±0.030	-	3 535 210	3 535 210 BL
30	±0.030	black	3 535 300	3 535 300 BL
32	±0.030	-	3 535 320	3 535 320 BL
35	±0.030	-	3 535 350	3 535 350 BL
38.5	±0.050	-	3 535 038	3 535 038 BL
39.06	±0.050	-	3 535 039	3 535 039 BL
40	±0.050	white	3 535 400	3 535 400 BL
60	±0.050	-	3 535 600	3 535 600 BL
70	±0.050	-	3 535 700	3 535 700 BL
75	±0.050	green	3 535 750	3 535 750 BL
80	±0.080	-	3 535 800	3 535 800 BL
90	±0.080	-	3 535 901	3 535 901 BL
200	±0.120	blue	3 535 920	3 535 920 BL
250	±0.120	-	3 535 925	3 535 925 BL

further sizes on request