


Balance

high-precision lab balance, 0.001g - 320/640g, percent weighing mode

0.001 g - 320 g
- 640g



WBA-620

WBA -320 /-620 /-620A high-precision

Ideal for:

- high precision balancing

Features:

- percent (%) weighing mode available
- with Ø90 mm/Ø118 mm/Ø128 mm weighing plate and glass draft shield (included)
- weighing response time adjustable
- electromagnetic base
- numeric linearity calibration
- built-in RS232C & USB port for printer
- super-size backlight LCD display, can display present temperature
- time and date setting for recording data
- stabilization time (typical at fast) ≤ 8 seconds
- battery operation possible (6V, 1.2A)
- CE certified and unique serial number for tracing

Safety mechanism:

- four level anti-shock function adjustable
- default and auto lock present working mode

200g (WBA-320) / 500g (WBA-620)
calibration weight included



NEW

Model	WBA-320 Ext-CAL	WBA-620 Ext-CAL	WBA-620A Int-CAL
Weighting range, weighing plate	0.001 - 320 g, Ø 90 mm	0.001 - 620 g, Ø 128 mm	0.001 - 620 g, Ø 118 mm
Read out, reproducibility & linearity	0.001 g		
Permissible temperature	5.0 - 35°C		
Weight modes	% mode / g / ct / GN / dwt / etc.		
Counting mode	available counting & weighing mode		
Calibration	external calibration, 200g weight included	external calibration, 500g weight included	internal calibration, AUTO
Dimension & weight	223 x 345 x 331 mm, 3.6 kg	223 x 345 x 331 mm, 5.7 kg	
Packaging size	290 x 385 x 365 mm, 5.1 kg	390 x 510 x 520 mm, 10.1 kg	
Power supply	DC 12V with AC 230V adapter, 50/60 Hz	DC 12V with AC 230V adapter, 50/60 Hz	
Order number	DH.WBA0320	DH.WBA0620	DH.WBA0620A



Thermal printer, for balance WBA, RS232 connection
for paper rolls (1 paper roll included)

Article	Order number
WBA-100 printer	DH.WBA100P
5x spare paper roll for WBA-100	DH.WBA0084



Thermal printer for balance WBA, RS232 connection
for label paper rolls (1 label paper roll included)

Article	Order number
WBA-420 printer	DH.WBA420P
5x spare label paper roll for WBA-420	DH.WBA0771



Sample plates, 10 pieces

Diameter	Order number
Ø90mm	DH.WBA0005