

WITOPET Mini

Quick pipetting

Your advantage at a glance:

- 6 fix-volumina from 5µl to 100µl
- compact size with an overall length of 130mm
- maximum accuracy
- feeling of a regular pipette
- two-stage plunger to ensure correct pipetting volumes
- smart tip cone design allows perfect fit for 200 µl as well as 20 µl tips
- fully autoclavable
- colour-code for better determination
- with quality certificate per batch (individual certificate on request)

The most competitive way of quick pipetting

WITOPET Mini

Volume µl	Volume tested µl	Systematic errors ± %	Random errors ≤ %	Order number
5	5	4.0	2.5	5 405 005
10	10	3.2	1.4	5 405 010
20	20	3.0	0.5	5 405 020
25	25	1.5	1.1	5 405 025
50	50	1.5	0.4	5 405 050
100	100	1.5	0.3	5 405 090