

Ultrasonic cleaner

analog, with lid, 1.2 l - 22 l, integrated controller, HF frequency 40 kHz

WUC-A analog, with lid

Ideal for:

- cleaning to remove dirt, grease, waxes and oil from all lab instruments by ultrasonic frequency

Features:

- made of stainless steel, finished surface
- convenient preset function of heating and timer
- HF frequency 40 kHz
- bottom-mounted drain valve (WUC-A10H and WUC-A22H)
- maximum 30 minutes cleaning time
- CE certified and unique serial number for tracing

Controller:

- **analog** controller for timer, temperature and output control




WUC-A06H with lid and valve



| Model | WUC-A01H* | WUC-A02H* | WUC-A03H | WUC-A06H | WUC-A10H | WUC-A22H |
|-----------------------------|---|-------------------------|-------------------------|-------------------------|-------------------------|--------------------------|
| Capacity | 1.2 l | 1.8 l | 3.3 l | 6 l | 10 l | 22 l |
| HF-peak-out | 72 W | 75 W | 124 W | 168 W | 265 W | 379 W |
| Heating temperature & power | 65°C fixed, 86 W | | ambient to 85°C, 172 W | | ambient to 85°C, 258 W | ambient to 85°C, 275 W |
| Timer | - | 0 - 30 min. | | | | |
| Drain-valve | - | | | | 1/2" NPT | |
| Inner tank size (WxDxH) | 150 x 135 x 60 mm | 150 x 135 x 100 mm | 240 x 140 x 100 mm | 290 x 150 x 150 mm | 290 x 240 x 150 mm | 500 x 300 x 150 mm |
| Weight | 1.8 kg | 2.5 kg | 4 kg | 5.3 kg | 7 kg | 11 kg |
| Packing size & gross weight | 200 x 200 x 190mm, 3 kg | 200 x 200 x 240mm, 4 kg | 320 x 200 x 290mm, 5 kg | 370 x 220 x 350mm, 6 kg | 434 x 320 x 360mm, 7 kg | 640 x 380 x 370mm, 11 kg |
| Power consumption | 158 W | 161 W | 296 W | 340 W | 523 W | 654 W |
| Power supply | 1 Phase AC 120V, 60 Hz or AC 230V, 50/60 Hz | | | | | |
| Order number 230V | DH.WUC.A01H | DH.WUC.A02H | DH.WUC.A03H | DH.WUC.A06H | DH.WUC.A10H | DH.WUC.A22H |
| Order number 120V | DH.WUC.A101H | DH.WUC.A102H | DH.WUC.A103H | DH.WUC.A106H | DH.WUC.A110H | DH.WUC.A122H |

Note: no basket included, please order separately.

* WUC-A01H and WUC-A02H with flip switch



Stainless steel wire baskets for ultrasonic cleaner WUC-A

| For device | Capacity | Measurements | Article | Order number |
|------------|----------|----------------|---------|--------------|
| WUC-A02H | 1.8 l | 120x105xh65 mm | B02 | DH.WUC.B02 |
| WUC-A03H | 3.3 l | 215x115xh65 mm | B03 | DH.WUC.B03 |
| WUC-A06H | 6 l | 247x123xh96 mm | B06 | DH.WUC.B06 |
| WUC-A10H | 10 l | 262x206xh93 mm | B10 | DH.WUC.B10 |
| WUC-A22H | 22 l | 452x263xh93 mm | B22 | DH.WUC.B22 |