

Spare filter

Description	For device...	Filter replacement period	Order number
14" Super A/C filter	-	approx. 3 - 4 month	DH.WWP.PR14AC
10" A/C filter	-		DH.WWP.PR10AC
10" prefilter	-	approx. 20 - 40 days	DH.WWP.10PF
Pure water Pack I	DH.WWP.RO180	approx. 1 year	DH.WWP.ROM1
Pure water Pack II	DH.WWP.RO280		DH.WWP.ROM2
Pure water Pack III	DH.WWP.RO380		DH.WWP.ROM3
0.2 µm final capsule filter	-	approx. 6 - 8 month	DH.WWP.0.2CF
Ultra pure water Pack, Scholar / Bio	DH.WWP.POWER.S , DH.WWP.POWER.B,	approx. 1 year	DH.WWP.UPS
Ultra pure water Pack, Scholar-UV / Integrate	DH.WWP.POWER.U, DH.WWP.POWER.I		DH.WWP.UPU
MW 5000 UF filter (M)	-	approx. 6 - 8 month	DH.WWP.MWM
MW 5000 UF filter (L)	-	approx. 6 - 8 month	DH.WWP.MWL

The above mentioned filter replacement period can change according to various factors such as feed water condition, producing volume and so on.



14" Super A/C filter



10" Prefilter



Pure water Pack I



Final capsule filter



10" A/C Filter



Ultra Pure water Pack, Scholar / Bio



Ultra pure water Pack, Scholar-UV / Integrate



Pure water Pack II

Spare UV lamps

Description	Using time	Order number
254/185 UV lamp	approx. 2,000 h	DH.WWP.UV

Bottle, HDPE, aspiration

Capacity	Neck I.Ø	Dimensions	Order number white	Order number brown
5 l	Ø70 mm	Ø200 x 360 mm	DH.Bot3102	DH.Bot3106
10 l	Ø70 mm	Ø250 x 460 mm	DH.Bot3103	DH.Bot3107
20 l	Ø95 mm	Ø300 x 555 mm	DH.Bot3104	DH.Bot3108
50 l	Ø95 mm	Ø380 x 700 mm	DH.Bot3105	DH.Bot3109

