

LF 3/3a/5/5a

Filter Holder

LF 30/31/32/33

Vacuum Filtration Set

Important notice

- Please study this instruction manual carefully before using this product.
- Please discard packing material according to related regulations.
- If you have any question in the product, please contact your local distributor or dealer for immediate help
- Filtration accessories are applicable to most of low concentrated acid, base, organic solution, but not all solution.
- Please clean and maintain the filter holder according to the maintenance method after every use.
- For more product details, please refer to our latest catalog and company website.


Unpacking

When unpacking, notice if all supplied parts are in good situation. If not, please contact the distributor or contact us by e-mail: export@rocker.com.tw immediately for assistance


Model	LF 3, PES Filter Holder
1	PP lid kit
2	300 ml PES funnel
3	Support base
4	Silicone stopper (No.16)




Model	LF 5, PES Filter Holder
1	PP lid kit
2	500 ml PES funnel
3	Support base
4	Silicone stopper (No.16)




Model	LF 3a, PES Filter Holder
1	PP lid kit
2	300 ml PES funnel
3	Tissue bottle adaptor




Model	LF 5a, PES Filter Holder
1	PP lid kit
2	500 ml PES funnel
3	Tissue bottle adaptor




Model	LF 30, PES Filtration Set
1	PP lid kit
2	300 ml PES funnel
3	Support base
4	Overflow protection
5	PC waste bottle




Model	LF 31, PES Filtration Set
1	PP lid kit
2	300 ml PES funnel
3	Support base
4	Silicone stopper (No.16)
5	1000 ml receiver flask with drain
6	Silicone fixing sucker



Model	LF 32, S.S. Filtration Set
1	100 ml S.S. filter holder
2	Silicone stopper (No.16)
3	Overflow protection
4	PC waste bottle



Model	LF 33, S.S. Filtration Set
1	300 ml S.S. filter holder
2	Silicone stopper (No.16)
3	Overflow protection
4	PC waste bottle



Application

Complete vacuum filtration system is ideal for suspended solid test and general lab filtration.



Ordering Information

- 197000-01 LF 3, 300 ml PES filter holder
- 197000-02 LF 3a, 300 ml PES filter holder
- 197000-03 LF 5, 500 ml PES filter holder
- 197000-04 LF 5a, 500 ml PES filter holder
- 197010-30 LF 30, PES filtration set
- 197010-31 LF 31, PES filtration set
- 197010-32 LF 32, Stainless steel filtration set
- 197010-33 LF 33, Stainless steel filtration set
- 167103-10 100 ml Stainless steel filter holder
- 167103-30 300 ml Stainless steel filter holder
- 197000-11-PC PC Waste Bottle (1200 ml)
- 167100-15 1000ml Receiver Glass Flask with draining and silicone fixing sucker
- 197000-36 PP lid kit for LF 3/3a
- 197000-39 Overflow protection
- 167110-17 Silicone stopper (No.16)
- 167200-45 Silicone fixing sucker