

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 010 053

Burettes, 5 ml, straight stopcock, AS, conf. cert.

Batch: 10-20




Serial:

Measuring Nr. 72525

Batch size:

10

Volume:	5,0000 ml	3,0000 ml	1,0000 ml
Tolerance:	$\pm 0,0100$ ml	$\pm 0,0100$ ml	$\pm 0,0100$ ml
Average value:	5,0001 ml	2,9959 ml	0,9984 ml
Deviation:	0,0021 ml	0,0019 ml	0,0031 ml
Precision (CV):	0,0425 %	0,0627 %	0,3111 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

03.04.2025

Balance:

Sartorius -
LE244S-OCE

DE-M