

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 010 503 B

Burettes, 50 ml, str. stopc., amber grad., AS, conf. cert.




Batch: 17-25

Serial:

Measuring Nr. 73366

Batch size: 199

Volume:	50,0000 ml	30,0000 ml	10,0000 ml
Tolerance:	$\pm 0,0500$ ml	$\pm 0,0500$ ml	$\pm 0,0500$ ml
Average value:	50,0153 ml	30,0000 ml	10,0066 ml
Deviation:	0,0068 ml	0,0085 ml	0,0053 ml
Precision (CV):	0,0136 %	0,0282 %	0,0533 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

23.06.2025

Balance:

Sartorius -
LE244S-OCE

DE-M