

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 012 052

**Burettes, 5 ml, straight stopc., glass, AS, conf. cert.**

**Batch:** 46-24




**Serial:**

**Measuring Nr.** 73663

**Batch size:**

5

|                        |                  |                  |                  |
|------------------------|------------------|------------------|------------------|
| <b>Volume:</b>         | <b>5,0000 ml</b> | <b>3,0000 ml</b> | <b>1,0000 ml</b> |
| <b>Tolerance:</b>      | $\pm 0,0100$ ml  | $\pm 0,0100$ ml  | $\pm 0,0100$ ml  |
| <b>Average value:</b>  | 5,0025 ml        | 3,0033 ml        | 1,0021 ml        |
| <b>Deviation:</b>      | 0,0011 ml        | 0,0009 ml        | 0,0003 ml        |
| <b>Precision (CV):</b> | 0,0222 %         | 0,0288 %         | 0,0323 %         |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:**

7

**Date:**

21.07.2025

**Balance:**

Sartorius -  
LE244S-OCE

DE-M